

Product datasheet for **SC300155**

Pancreatic Lipase (PNLIP) (NM_000936) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Pancreatic Lipase (PNLIP) (NM_000936) Human Untagged Clone
Tag:	Tag Free
Symbol:	PNLIP
Synonyms:	PL; PNLIPD; PTL
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC300155 representing NM_000936.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCTGCCACTTTGGACTCTTTCACTGCTGGGAGCAGTAGCAGGAAAAGAAGTTTGCTACGAAAGA
CTCGGCTGCTTCAGTGATGACTCCCATGGTCAGGAATTACGAAAGACCCTCCATATATTGCCTTGG
TCTCCAAAAGATGTCAACACCCGCTTCCCTCTATACTAATGAGAACCACCAACTTTCAAGAAGTT
GCCGCAGATTCATCAAGCATCAGTGGCTCCAATTTCAAAACAAATAGAAAAACTCGCTTTATTATTCAT
GGATTCATAGACAAGGGAGAAGAAAAGTGGCTGGCCAATGTGTGCAAGAACTGTTCAGGTGGAAGT
GTGAACGTGATCTGTGTGGACTGGAAAGGTGGCTCCCGAACTGGATACACACAAGCCTCGCAGAACATC
AGGATCGTGGGAGCAGAAGTGGCATATTTTGTGAATTTCTTCAGTCGGCGTTTCGGTTACTCACCTTCC
AACGTGCATGTCATTGGCCACAGCCTGGGTGCCACGCTGCTGGGAGGCTGGAAGGAGAACCAATGGG
ACCATTGGACGCATCACAGGGTTGGACCCAGCAGAACCTTGCTTTCAGGGCACACCTGAATTAGTCCGA
TTGGACCCAGCGATGCCAAATTTGTGGATGTAATTCACACGGATGGTGCCTCCATAGTCCCAATTTG
GGGTTTGGAAATGAGCCAAGTCGTGGGCCACCTAGATTTCTTTCCAAATGGAGGAGTGAAATGCCTGGA
TGTAAGAAAGAACATTCTCTCAGATTGTGGACATAGACGGAATCTGGGAGGGACTCGAGACTTTGCG
GCCTGTAATCACTTAAGAAGCTACAAAATATTACACTGATAGCATCGTCAACCCTGATGGCTTTGCTGGA
TTCCCCTGTGCCTCTTACAACGCTTCACTGCAAAACAAGTGTTCCTTGTCCAAGTGGAGGCTGCCCA
CAGATGGGTCACTATGCTGATAGATATCCTGGGAAAACAATGATGTGGCCAGAAATTTTATCTAGAC
ACTGGTGTGCCAGTAATTTTGCACGTTGGAGGTAAAGGTATCTGTCACTGTCTGGAAAAAGGTT
ACAGGACACATACTAGTTTCTTTGTTTCGAAATAAAGGAACTCTAAGCAGTATGAAATTTCAAGGGC
ACTCTCAAACCAGATAGTACTCATTCCAATGAATTTGACTCAGATGTGGATGTTGGGGACTGCAGATG
GTTAAATTTATTGGTATAACAATGTGATCAACCAACTTTACCTAGAGTGGGAGCATCCAAGATTATA
GTGGAGACAAATGTTGGAAAACAGTTCAACTTCTGTAGTCCAGAAACCGTCAGGGAGGAAGTTCTGCTC
ACCCTCACACCGTGTTAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites:

Sgfl-MluI

Plasmid Map:



ACCN: NM_000936

Insert Size: 1398 bp

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_000936.3
RefSeq Size:	1510 bp
RefSeq ORF:	1398 bp
Locus ID:	5406
UniProt ID:	P16233
Cytogenetics:	10q25.3
Protein Families:	Druggable Genome, Secreted Protein
Protein Pathways:	Glycerolipid metabolism, Metabolic pathways
MW:	51.2 kDa
Gene Summary:	<p>This gene encodes a member of the lipase family of proteins. The encoded enzyme is secreted by the pancreas and hydrolyzes triglycerides in the small intestine, and is essential for the efficient digestion of dietary fats. Inhibition of the encoded enzyme may prevent high-fat diet-induced obesity in mice and result in weight loss in human patients with obesity. Mutations in this gene cause congenital pancreatic lipase deficiency, a rare disorder characterized by steatorrhea. [provided by RefSeq, Jul 2016]</p>