

Product datasheet for **SC300197**

ZDHHC13 (NM_001001483) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZDHHC13 (NM_001001483) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZDHHC13
Synonyms:	HIP3RP; HIP14L
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC300197 representing NM_001001483.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGTCATATTACTCCAGCATGGTGCAGACCCCACTCTTATTGATGGAGAGGGATTCCAGCAGCATC
CACCTGGCAGTATTGTTTCAACACATGCCTATTATAGCATATCTCATCTCAAAGGGACAGAGTGTGAAT
ATGACAGATGAAATGGGCAGACACCTCTCATGTTATCAGCTCACAAAGTAATTGGGCCAGAACCAACT
GGATTTCTTTAAAGTTTAACTCTCTCAATGTGGTTGATAAAATACACCAAAACACTCCAATTAC
TGGGCAGTTGCAGCAGGAAATGTTAATGCAGTTGATAAGCTTTTGGAAAGCTGGTTCTAGCCTGGATATC
CAGAATGTTAAGGGGAAAAACACCTCTTGATATGGCTCTACAAAACAAAAATCAGCTCATTATTCATATG
CTAAAAACAGAAGCCAAATGAGAGCCAACCAAAAGTTCCAGACTTTGGAGGTGGCTGCAGAAATGCGAG
CTCTTCTGCTGCTGATGCTTTCTGTGATTACCATGTGGGCTATTGGATACATATTGGACTTCAATTCA
GATTCTTGGCTTTTAAAGGATGCTTCTAGTAACACTGTTTTTCTGACATCTTTGTTTCCAAGGTTTCT
TTGGTTGGGTATAAGAACCTGTATACTTACCAACAGCCTTTCTGCTAAGTTCTGTTTTTGGATATTT
ATGACTTGGTTTCACTTATTTTTTCTGATTTAGCAGGAGCCCTTTCTATTTTCAGTTTCATTTTTACG
ATAGTAGCCTTTCTATACTTTTTCTATAAGACTTGGGCAACTGATCCAGGCTTCACTAAGCTTCTGAA
GAAGAAAAGAAAGTGAATATCATCACCTTGCAGAAACTGGCTCTCTGGACTTCAGAACATTTTGTACA
TCATGTCTTATAAGGAAGCCATTAAGGTCCTCCACTGCCATGTATGCAACTGCTGTGTGGCTCGATAT
GATCAACACTGCCTGTGGACTGGACGGTGCATAGGTTTTGGCAACCATCACTATTACATATTTCTTTG
TTTTTCTTTCCATGGTATGTGGCTGGATTATATGGATCTTTCATCTATTTGTCCAGTCATTGTGCC
ACAACATTCAAAGAAGATGGATTATGGACTTACCTCAATCAGATTGTGGCCTGTTCCCTTGGGTTTTA
TATATCTTGATGCTAGCAACTTTCCATTTCTCATGGTCAACATTTTATTATTAATCAACTTTTCAG
ATTGCCTTTCTGGCCTGACCTCCATGAGAGAATCAGCCTGCAGAAGCAGAGCAAGCATATGAACAG
ACGTTGTCCCTCAGGAAGACACCATAAATCTTGATTATGCAGAACCTGGCAGATTTCTTTCAGTGT
GGCTGCTTTGGCTTGGTGAAGCCCTGTGTGGTAGATTGGACATCACAGTACACCATGGTCTTTCACCCA
GCCAGGGAGAAGGTTCTTCGCTCAGTATGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites:

SgfI-MluI

Plasmid Map:



ACCN:

NM_001001483

Insert Size:	1479 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001001483.2
RefSeq Size:	2302 bp
RefSeq ORF:	1479 bp
Locus ID:	54503
UniProt ID:	Q8IUH4
Cytogenetics:	11p15.1
Protein Families:	Transmembrane
MW:	56.5 kDa
Gene Summary:	Palmitoyltransferase for HTT and GAD2. May play a role in Mg(2+) transport. [UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (2) lacks an alternate exon in the 5' region which results in translation initiation at a downstream start codon, compared to variant 1. The resulting isoform (2) has a shorter N-terminus, compared to isoform 1.