

Product datasheet for **SC300217**

GGA1 (NM_001001560) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GGA1 (NM_001001560) Human Untagged Clone
Tag:	Tag Free
Symbol:	GGA1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC300217 representing NM_001001560.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGAGCCCGCGATGGAGCCGGAGACTCTGGAGGCGCGAATCAATAGAGCCACGAACCCCTGAACAAG
GAGCTCGACTGGCCAGCATCAACGGCTTCTGCGAGCAGCTCAACGAGGACTTTAGGGGCTCCACTC
GCCACCCGGTCTGGCCACAAGATCCAGTCCCACAGGAGTGGGAGGCGATCCAGGCCTTGACGGTG
CTGAAAACATGCATGAAGAGCTGCGGCAAGCGGTTCCACGACGAAGTGGGCAAGTTCCGCTTTCTCAAC
GAGCTCATCAAGGTCGTGTCTCCCAAGTATCTGGGCTCTCGGACATCGGAGAAGGTGAAGAACAAGATC
TTGGAGCTCCTCTACAGCTGGACAGTGGGCTGCCCGAGGAGGTGAAAATCGCAGAGGCCTACCAGATG
CTAAAGAAGCAGGGGATTGTAAGTCCGACCCCAAGCTTCCAGATGACTACTCTTCCCCTTCTCCT
CCACGGCCGAAGAATGTGATCTTTGAAGATGAGGAGAAATCCAAGATGCTGGCCCGCTGCTGAAGAGC
TCCCATCCGAAGACCTCCGCGCAGCCAATAAGCTCATCAAAGAGATGGTGCAGGAGGACCAGAAGCGG
ATGGAGAAGATCTCGAAGAGGGTGAATGCCATCGAGGAGGTGAACAACAATGTGAAACTGCTCACGGAG
ATGGTGATGAGCCACAGCCAGGGCGGCGCAGCAGCTGGCAGCAGGAGGACCTCATGAAGGAACTGTAC
CAGCGCTGTGAGCGGATGCGGCCACGCTCTTCCGACTGGCGAGTGACACAGAGGACAATGATGAGGCC
TTAGGCCTCAGTGACCCACACCCCTTTCAGGCCAAGCCTGGATGGTACCGGATGGAACAGCTTCCAG
TCGTGGATGCCACTGAGCCCCAGCCCTGCTCTGGCCAGGCCCCAGTATGGAAGCCGACCCCA
GCGCAGACATCCTGCCAGCAAGCAGCGGTCTGGACGACCTAGACCTCCTGGGAAGACCCCTCTGCAG
CAGTCGCTGCCCCGGAATCCCAGCAAGTGGTGGGAGAAGCAGCAGCCAACCCCGGCTCACACTC
CGGGACCTGCAGAATAAGAGCAGCAGCTGCAGCTCCCCAGCTCCAGGCCACCAGCCTTCTCCACACC
GTGTCCCAGAGCCCCCAGCCCTCCGACGACCCCGTACCAACCGAGCTCTACTGGCCAGCATCACT
GTGCCCTGGAGTCCATCAAACCCAGCAACATCCTGCCGTGACTGTGTATGACCAGCAGCGCTTCCGC
ATCCTTCTCCATTTTGGCCGGGACCCACTGCCAGGGCGCTCCGACGTGCTGGTGGTGGTGGTTCCATG
CTGAGCACCGCCCCAGCCATCCGCAACATCGTGTCCAGTCAGTGTCCCAAGTTATGAAGGTG
AAGCTGCAGCCACCCTCGGGACGGAGCTGCCAGCTTTTAACCCATCGTCCACCCTCAGCAATCACC
CAGGTCTGTGCTTGCCAACCCCAAGGAGAAGGTTGCCTCCGCTACAAGCTCACCTTACCATG
GGTGACCAGACCTACAACGAGATGGGGATGTGGACCAGTCCCCCACCTGAAACCTGGGGTAGCCTC
TAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites:

SgfI-MluI

Plasmid Map:



ACCN:	NM_001001560
Insert Size:	1659 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001001560.2
RefSeq Size:	2924 bp
RefSeq ORF:	1659 bp
Locus ID:	26088
UniProt ID:	Q9UJY5
Cytogenetics:	22q13.1
Protein Families:	Druggable Genome
Protein Pathways:	Lysosome
MW:	61.4 kDa
Gene Summary:	<p>This gene encodes a member of the Golgi-localized, gamma adaptin ear-containing, ARF-binding (GGA) protein family. Members of this family are ubiquitous coat proteins that regulate the trafficking of proteins between the trans-Golgi network and the lysosome. These proteins share an amino-terminal VHS domain which mediates sorting of the mannose 6-phosphate receptors at the trans-Golgi network. They also contain a carboxy-terminal region with homology to the ear domain of gamma-adaptins. Multiple alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (2) lacks 2 in-frame exons in the coding region, as compared to variant 1. The encoded isoform 2 thus lacks an internal segment, as compared to isoform 1.</p>