

Product datasheet for **SC300218**

GGA1 (NM_001001561) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GGA1 (NM_001001561) Human Untagged Clone
Tag:	Tag Free
Symbol:	GGA1
Synonyms:	ADP-ribosylation factor binding protein 1; gamma-adaptin related protein 1; golgi associated, gamma adaptin ear containing, ARF binding protein 1; OTTHUMP00000028975; OTTHUMP00000042200
Vector:	<u>pCMV6 series</u>
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001001561, the custom clone sequence may differ by one or more nucleotides ATGGAGCCCGGATGGAGCCGGAGACTCTGGAGGCGGAATCAATAGAGCCACGAACCCC CTGAACAAGGAGCTCGACTGGGCCAGCATCAACGGCTTCTGCGAGCAGCTCAACGAGGAC TTTGAGGGGCTCCACTCGCCACCCGGCTGCTGGCCACAAGATCCAGTCCCCACAGGAG TGGGAGGCGATCCAGGCCTTGACGGTGAGAAGGGGAGAGGCCACCATCCGTCCCCGCCA TGTGACGACACCAAGGGAGGCCAAGACTGA
Restriction Sites:	Please inquire
ACCN:	NM_001001561
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001001561.1, NP_001001561.1</u>
RefSeq Size:	1723 bp
RefSeq ORF:	270 bp
Locus ID:	26088
Cytogenetics:	22q13.1
Protein Families:	Druggable Genome
Protein Pathways:	Lysosome
Gene Summary:	<p>This gene encodes a member of the Golgi-localized, gamma adaptin ear-containing, ARF-binding (GGA) protein family. Members of this family are ubiquitous coat proteins that regulate the trafficking of proteins between the trans-Golgi network and the lysosome. These proteins share an amino-terminal VHS domain which mediates sorting of the mannose 6-phosphate receptors at the trans-Golgi network. They also contain a carboxy-terminal region with homology to the ear domain of gamma-adaptins. Multiple alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (3) lacks multiple 3' exons and has an alternate 3' end, as compared to variant 1. The encoded isoform 3 has a shorter and distinct C-terminus, as compared to isoform 1.</p>