

## Product datasheet for **SC300222**

### **PDE9A (NM\_001001569) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	PDE9A (NM_001001569) Human Untagged Clone
Tag:	Tag Free
Symbol:	PDE9A
Synonyms:	HSPDE9A2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


**Fully Sequenced ORF:** >SC300222 representing NM\_001001569.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGGATCCGGCTCCTCCAGCTACCGGCCAAGGCCATCTACCTGGACATCGATGGACGCATTTCAGAAG
CACTCCGTACAAAGTGAGACCTGTGGCCATCAAGCAACTCTCCGAGCATTCAAAATCAATGAACTGAAA
GCTGAAGTTGCAAATCACTTGGCTGTCTAGAGAAACGGTGGAATTGGAAGGACTAAAAGTGGTGGAG
ATTGAGAAATGCAAGAGTGACATTAAGAAGATGAGGGAGGAGCTGGCGGCCAGAAGCAGCAGGACCAAC
TGCCCTGTAAGTACAGTTTTTTGGATAACCACAAGAAGTTGACTCCTCGACGCGATGTTCCCACTTAC
CCCAAGTACCTGCTCTCTCCAGAGACCATCGAGGCCCTGCGGAAGCCGACCTTTGACGTCTGGCTTTGG
GAGCCCAATGAGATGCTGAGCTGCCTGGAGCAGTGTACCAGCAGCTCGGGCTGGTCAGGGACTTCAGC
ATCAACCTGTCAACCTCAGGAGGTGGCTGTTCTGCGTCCAGCAACTACAGAAACAACCCCTCCAC
AACTCCGGCACTGCTTCTGCGTGGCCAGATGATGTACAGCATGGTCTGGCTCTGCAGTCTCCAGGAG
AAGTTCTCACAACGGATATCCTGATCCTAATGACAGCGGCCATCTGCCAGCATCTGGACCATCCCGGC
TACAACAACACGTACCAGATCAATGCCCGCACAGAGCTGGCGGTCCGCTACAATGACATCTCACCGCTG
GAGAACCACCACTGCGCGTGGCCTTCCAGATCCTCGCCGAGCCTGAGTGCAACATCTTCTCCAACATC
CCACCTGATGGGTTCAAGCAGATCCGACAGGGAATGATCACATTAATCTTGGCCACTGACATGGCAAGA
CATGCAGAAATATGGATTCTTTCAAAGAGAAAATGGAGAATTTTGACTACAGCAACGAGGAGCAGATG
ACCCTGCTGAAGATGATTTTGATAAAATGCTGTGATATCTCTAACGAGGTCGGTCCAATGGAAAGTCGCA
GAGCCTTGGTGGACTGTTTATTAGAGGAATATTTTATGCAGAGCGACCGTGAGAAGTCAGAAGGCCTT
CCTGTGGCACCCTCATGGACCGAGACAAAGTGACCAAGGCCACAGCCAGATTGGGTTCCATCAAGTTT
GTCCTGATCCCAATGTTTGAACAGTGACCAAGCTCTTCCCATGTTTGGAGAGATCATGTCCAGCCA
CTTTGGGAATCCCGAGATCGCTACGAGGAGCTGAAGCGGATAGATGACGCCATGAAAGAGTTACAGAAG
AAGACTGACAGCTTGACGTCTGGGGCCACCGAGAAGTCCAGAGAGAGAAGCAGAGATGTGAAAAACAGT
GAAGGAGACTGTGCCTGA
ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATTACAAG
GATGACGACGATAAGGTTTAAACGGCCGGCCGGCT
  
```

**Restriction Sites:** SgfI-NotI

**Plasmid Map:**



**ACCN:** NM\_001001569

**Insert Size:** 1398 bp

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_001001569.1</a></u>
<b>RefSeq Size:</b>	1719 bp
<b>RefSeq ORF:</b>	1398 bp
<b>Locus ID:</b>	5152
<b>UniProt ID:</b>	<u><a href="#">O76083</a></u>
<b>Cytogenetics:</b>	21q22.3
<b>Protein Families:</b>	Druggable Genome
<b>Protein Pathways:</b>	Progesterone-mediated oocyte maturation, Purine metabolism
<b>MW:</b>	54.1 kDa
<b>Gene Summary:</b>	<p>The protein encoded by this gene catalyzes the hydrolysis of cAMP and cGMP to their corresponding monophosphates. The encoded protein plays a role in signal transduction by regulating the intracellular concentration of these cyclic nucleotides. Multiple transcript variants encoding several different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (4) lacks four alternate exons compared to variant 1. The resulting predicted isoform (d) is shorter and has a distinct N-terminus compared to isoform a.</p>