

Product datasheet for **SC300227**

PDE9A (NM_001001574) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PDE9A (NM_001001574) Human Untagged Clone
Tag:	Tag Free
Symbol:	PDE9A
Synonyms:	HSPDE9A2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC300227 representing NM_001001574.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGGATCCGGCTCCTCCAGCTACCGGCCAAGGCCATCTACCTGGACATCGATGGACGCATTCAAG
GTAATCTTCAGCAAGTACTGCAACTCCAGCGACATCATGGACCTGTTCTGCATCGCCACCGCCTGCCT
CGCACTCCGTACAAAGTGAGACCTGTGGCCATCAAGCAACTCTCCGAGAGAGAAGAATTAATCCAGAGC
GTGCTGGCGCAGGTTGCAGAGCAGTTCTCAAGAGCATTCAAATCAATGAACTGAAAGCTGAAGTTGCA
AATCACTTGGCTGTCCTAGAGAAACGCGTGAATTGGAAGGACTAAAAGTGGTGGAGATTGAGAAATGC
AAGAGTGACATTAAGAAGATGAGGGAGGAGCTGGCGCCAGAAGCAGCAGGACCAACTGCCCTGTAAG
TACAGTTTTTTGGATAACCACAAGAAGTTGACTCCTCGACGCGATGTTCCCACTTACCCCAAGTACCTG
CTCTCTCCAGAGACCATCGAGGCCCTGCGGAAGCCGACCTTTGACGTCTGGCTTTGGGAGCCCAATGAG
ATGCTGAGCTGCCTGGAGCAGATGTACCACGACCTCGGGTGGTCAGGGACTTCAGCATCAACCCTGTC
ACCCTCAGGAGGTGGCTGTTCTGCGTCCACGACAACTACAGAAACAACCCCTTCCACAACCTCCGGCAC
TGCTTCTGCGTGGCCAGATGATGTACAGCATGGTCTGGCTCTGCAGTCTCCAGGAGAAGTTCTCACAA
ACGGATATCCTGATCCTAATGACAGCGGCCATCTGCCACGATCTGGACCATCCCGCTACAACAACAG
TACCAGATCAATGCCGCACAGAGCTGGCGGTCCGCTACAATGACATCTCACCGCTGGAGAACCACCAC
TGCGCCGTGGCCTTCCAGATCCTCGCCGAGCCTGAGTGCAACATCTTCTCCAACATCCCACCTGATGGG
TTCAAGCAGATCCGACAGGGAATGATCACATTAATCTTGGCCACTGACATGGCAAGACATGCAGAAAT
ATGATTTTGATAAAATGCTGTGATATCTAACGAGGTCCTCCAATGGAAGTCGAGAGCCTTGGGTG
GACTGTTTATTAGAGGAATATTTTATGCAGAGCGACCCGTGAGAAGTCAGAAGGCCCTTCTGTGGCACC
TTATGAGCCGAGACAAAGTGACCAAGGCCACAGCCAGATTGGGTTTCAAGTTTGTCTCTGATCCCA
ATGTTTGAACAGTGACCAAGCTTCTCCCATGGTTGAGGAGATCATGCTGCAGCCACTTTGGGAATCC
CGAGATCGCTACGAGGAGCTGAAGCGGATAGATGACGCCATGAAAGAGTTACAGAAGAAGACTGACAGC
TTGACGTCTGGGGCCACCAGAGAAGTCCAGAGAGAGAAGCAGAGATGTGAAAAACAGTGAAGGAGACTGT
GCCTGA
ACGCGGCCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATTACAAG
GATGACGACGATAAGGTTTAAACGGCCGGCCGGGT
  
```

Restriction Sites: SgfI-NotI

Plasmid Map:



ACCN: NM_001001574

Insert Size:	1524 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001001574.1</u>
RefSeq Size:	1845 bp
RefSeq ORF:	1524 bp
Locus ID:	5152
UniProt ID:	<u>O76083</u>
Cytogenetics:	21q22.3
Protein Families:	Druggable Genome
Protein Pathways:	Progesterone-mediated oocyte maturation, Purine metabolism
MW:	58.9 kDa
Gene Summary:	<p>The protein encoded by this gene catalyzes the hydrolysis of cAMP and cGMP to their corresponding monophosphates. The encoded protein plays a role in signal transduction by regulating the intracellular concentration of these cyclic nucleotides. Multiple transcript variants encoding several different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (9) lacks two alternate in-frame exons compared to variant 1. The resulting isoform (h) has the same N- and C-termini but lacks two internal segments compared to isoform a.</p>