

Product datasheet for **SC300231**

PDE9A (NM_001001578) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PDE9A (NM_001001578) Human Untagged Clone
Tag:	Tag Free
Symbol:	PDE9A
Synonyms:	HSPDE9A2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC300231 representing NM_001001578.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGGATCCGGCTCCTCCAGCTACCGGCCAAGGCCATCTACCTGGACATCGATGGACGCATTCCAGAAG
GAACACGACCATCTCCCTGTGACCACCGACGACGCCATGGTCTCCATCGACCCACCATGCCCGCGAA
TTCAGAACGCACTCCGTACAAAGTGAGACCTGTGGCCATCAAGCAACTCTCCGAGCATTCAAAATCAAT
GAACTGAAAGCTGAAGTTGCAAATCACTTGGCTGTCTAGAGAAAACGCGTGAATTTGGAAGGACTAAAA
GTGGTGGAGATTGAGAAATGCAAGAGTGACATTAAGAAGATGAGGGAGGAGCTGGCGCCAGAAGCAGC
AGGACCAACTGCCCTGTAAGTACAGTTTTTTGGATAACCACAAGAAGTTGACTCCTCGACGCGATGTT
CCCCTTACCCCAAGTACCTGCTCTCTCCAGAGACCATCGAGGCCCTGCGGAAGCCGACCTTTGACGTC
TGGCTTTGGGAGCCCAATGAGATGCTGAGCTGCCTGGAGCACATGTACCACGACCTCGGGCTGGTCAGG
GACTTCAGCATCAACCTGTACCCTCAGGAGGTGGCTGTTCTGCGTCCACGACAACTACAGAAACAAC
CCCTTCCACAACCTCCGGCACTGCTTCTGCGTGGCCAGATGATGTACAGCATGGTCTGGCTCTGCAGT
CTCCAGGAGAAGTTCTCACAAACGGATATCCTGATCCTAATGACAGCGGCCATCTGCCACGATCTGGAC
CATCCCGCTACAACAACACGTACCAGATCAATGCCCGCACAGAGCTGGCGGTCCGCTACATGACATC
TACCCGCTGGAGAACCACCACTGCGCCGTGGCCTTCCAGATCCTCGCCGAGCCTGAGTGCAACATCTTC
TCCAACATCCCACCTGATGGGTTCAAGCAGATCCGACAGGGAATGATCACATTAATCTTGGCCACTGAC
ATGGCAAGACATGCAGAAATTATGGATTCTTTCAAAGAGAAAAATGGAGAATTTTGACTACAGCAACGAG
GAGCACATGACCTGTGAAGATGATTTTGATAAAATGCTGTGATATCTTAACGAGGTCGTCCTCAATG
GAAGTCGACAGCCTTGGGTGGACTGTTTATTAGAGGAATATTTTATGCAGAGCGACCCGATGAGAAGTCA
GAAGGCCCTTCTGTGGCACCCTTTCATGGACCGAGACAAAGTGACCAAGGCCACAGCCAGATTGGGTTT
ATCAAGTTTGTCTGATCCCAATGTTTGAACAGTGACCAAGCTCTTCCCATGGTTGAGGAGATCATG
CTGCAGCCACTTTGGGAATCCCAGATCGCTACGAGGAGCTGAAGCGGATAGATGACGCCATGAAAGAG
TTACAGAAGAAGACTGACAGCTTACGCTCTGGGCCACCGAGAAGTCCAGAGAGAGAAGCAGAGATGTG
AAAAACAGTGAAGGAGACTGTGCCTGA
ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATTACAAG
GATGACGACGATAAGGTTTAAACGGCCGGCCGCGGT
  
```

Restriction Sites: SgfI-NotI

Plasmid Map:



ACCN: NM_001001578

Insert Size:	1476 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001001578.1</u>
RefSeq Size:	1797 bp
RefSeq ORF:	1476 bp
Locus ID:	5152
UniProt ID:	<u>O76083</u>
Cytogenetics:	21q22.3
Protein Families:	Druggable Genome
Protein Pathways:	Progesterone-mediated oocyte maturation, Purine metabolism
MW:	57.3 kDa
Gene Summary:	<p>The protein encoded by this gene catalyzes the hydrolysis of cAMP and cGMP to their corresponding monophosphates. The encoded protein plays a role in signal transduction by regulating the intracellular concentration of these cyclic nucleotides. Multiple transcript variants encoding several different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (13) lacks three alternate exons compared to variant 1. The resulting isoform (l) has the same N- and C-termini but lacks two internal segments compared to isoform a.</p>