

Product datasheet for **SC300233**

PDE9A (NM_001001580) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: PDE9A (NM_001001580) Human Untagged Clone
Tag: Tag Free
Symbol: PDE9A
Synonyms: HSPDE9A2
Vector: pCMV6 series

Fully Sequenced ORF: >NCBI ORF sequence for NM_001001580, the custom clone sequence may differ by one or more nucleotides

```
ATGACCAACTGCCCTGTAAGTACAGTTTTTTGGATAACCCACAAGAAGTTGACTCCTCGA
CGCGATGTTCCCACTTACCCCAAGTACCTGCTCTCTCCAGAGACCATCGAGGCCCTGCGG
AAGCCGACCTTTGACGTCTGGCTTTGGGAGCCCAATGAGATGCTGAGCTGCCTGGAGCAC
ATGTACCACGACCTCGGGCTGGTCAGGGACTTCAGCATCAACCCTGTACCCTCAGGAGG
TGGCTGTTCTGCGTCCACGACAACTACAGAAACAACCCCTTCCACAACCTCCGGCACTGC
TTCTGCGTGGCCAGATGATGTACAGCATGGTCTGGCTCTGCAGTCTCCAGGAGAAGTTC
TCACAAACGGATATCCTGATCCTAATGACAGCGGCCATCTGCCAGATCTGGACCATCCC
GGCTACAACAACAGTACCAGATCAATGCCCGCACAGAGCTGGCGTCCGCTACAATGAC
ATCTCACCGCTGGAGAACCACCACTGCGCCGTGGCCTTCCAGATCCTCGCCGAGCCTGAG
TGCAACATCTTCTCCAACATCCCACCTGATGGGTTCAAGCAGATCCGACAGGGAATGATC
ACATTAATCTTGGCCACTGACATGGCAAGACATGCAGAAATTAATGGATTCTTTCAAAGAG
AAAATGGAGAATTTTACTACAGCAACGAGGAGCACATGACCCTGCTGAAGATGATTTTG
ATAAAATGCTGTGATATCTTAACGAGGTCCTGCAATGGAAGTCGCAGAGCCTTGGGTG
GACTGTTTATTAGAGGAATTTTTATGCAGAGCGACCGTGAGAAGTCAGAAGGCCTTCTC
GTGGCACCGTTCATGGACCGAGACAAAGTGACCAAGGCCACAGCCAGATTGGGTTTCATC
AAGTTTGCCTGATCCCAATGTTTGAACAGTGACCAAGCTCTTCCCATGGTTGAGGAG
ATCATGCTGCAGCCACTTTGGGAATCCCAGATCGCTACGAGGAGCTGAAGCGGATAGAT
GACGCCATGAAAGAGTTACAGAAGAAGACTGACAGCTTGACGTCTGGGGCCACCGAGAAG
TCCAGAGAGAGAAGCAGAGATGTGAAAAACAGTGAAGGAGACTGTGCCTGA
```

Restriction Sites: Please inquire

ACCN: NM_001001580

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001001580.1, NP_001001580.1</u>
RefSeq Size:	1590 bp
RefSeq ORF:	1131 bp
Locus ID:	5152
UniProt ID:	<u>O76083</u>
Cytogenetics:	21q22.3
Protein Families:	Druggable Genome
Protein Pathways:	Progesterone-mediated oocyte maturation, Purine metabolism
Gene Summary:	<p>The protein encoded by this gene catalyzes the hydrolysis of cAMP and cGMP to their corresponding monophosphates. The encoded protein plays a role in signal transduction by regulating the intracellular concentration of these cyclic nucleotides. Multiple transcript variants encoding several different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (15) lacks six alternate exons compared to variant 1. The resulting predicted isoform (j) is shorter and has a distinct N-terminus compared to isoform a. Variants 11 and 15 both encode isoform j.</p>