

Product datasheet for **SC300236**

PDE9A (NM_001001583) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PDE9A (NM_001001583) Human Untagged Clone
Tag:	Tag Free
Symbol:	PDE9A
Synonyms:	HSPDE9A2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC300236 representing NM_001001583.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGGATCCGGCTCCTCCAGCTACCGGCCAAGGCCATCTACCTGGACATCGATGGACGCATTGAGAA
GTAACTTTCAGCAAGTACTGCACTCCAGCGACATCATGGACCTGTTCTGCATCGCCACCGCCTGCCT
CGCACTCCGTACAAAGTGAACCTGTGGCCATCAAGCACTCTCCGCTGGTGTGAGGACAAGAGAACC
ACAAGCCGTGGCCAGTCTGCTGAGAGACCACTGAGGGACAGACGGTGTGGCCTGGAGCAGCCCCGG
AGGGAAGGAGCATTGAAAGTGGACAGGTAGAGCCAGGCCAGAGAGCCCCAGGGCTGCTACCAGGAA
GGCCAGCGCATCCCTCCAGAGAGAGAAGAATTAATCCAGAGCGTGTGGCGCAGTTGCAGAGCAGTTC
TCAAGAGCATTCAAAATCAATGAACTGAAAGCTGAAGTTGCAATCACTTGGCTGTCTAGAGAAACGC
GTGGAATTGGAAGGACTAAAAGTGGTGGAGATTGAGAAATGCAAGAGTGACATTAAGAAGATGAGGGAG
GAGCTGGCGGCCAGAAGCAGCAGGACCACTGCCCTGTAAGTACAGTTTTTTGGATAACCACAAGAAG
TTGACTCCTCGACGCGATGTTCCCACTTACCCCAAGTACCTGCTCTCCAGAGACCATCGAGGCCCTG
CGGAAGCCGACCTTTGAGTCTGGCTTTGGGAGCCCAATGAGATGCTGAGCTGCCTGGAGCACATGTAC
CAGGACCTCGGGCTGGTCAAGGACTTCAAGCATCAACCCTGTACCCCTCAGGAGGTGGCTGTTCTGCGTC
CAGCAACTACAGAAACAACCCCTTCCACAACCTCCGCACTGTTCTGCGTGGCCAGATGATGTAC
AGCATGGTCTGGCTCTGCAGTCTCCAGGAGAAGTTCTCACAACGGATATCCTGATCCTAATGACAGCG
GCCATCTGCCAGCATCTGGACCATCCCGGCTACAACAACACGTACCAGATCAATGCCCGCACAGAGCTG
GCGGTCCGCTACAATGACATCTCACCGCTGGAGAACCACCACTGCGCCGTGGCCTTCCAGATCCTCGCC
GAGCCTGAGTGCAACATCTTCTCAACATCCCACCTGATGGGTTCAAGCAGATCCGACAGGGAATGATC
ACATTAATCTTGGCCACTGACATGGCAAGACATGCAGAAATTATGGATTCTTTCAAAGAGAAAATGGAG
AATTTTGACTACAGCAACGAGGAGCACATGACCCCTGCTGAAGATGATTTTGATAAAATGCTGTGATC
TCTAACGAGTCCGTCATGGAAGTGCAGAGCCTTGGGTGGACTGTTTATTAGAGGAATATTTTATG
CAGAGCGACCGTGAGAAGTCAAGGACCTTCTGTGGCACCGTTCATGGACCGAGACAAAGTGACCAAG
GCCACAGCCAGATTGGGTTTCAAGTTTGTCTGATCCCAATGTTTGAACAGTGACCAAGCTCTTC
CCCATGGTTGAGGAGATCATGCTGCAGCCACTTTGGGAATCCCGAGATCGCTACGAGGAGCTGAGCGG
ATAGATGACGCCATGAAAGAGTTACAGAAGAAGACTGACAGCTTGACGTCTGGGGCCACCGAGAAGTCC
AGAGAGAGAAGCAGAGATGTGAAAAACAGTGAAGGAGACTGTGCTGA
ACGGCGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATTACAAG
GATGACGACGATAAGGTTTAAACGGCCGGCCGCGGT
```

Restriction Sites: SgfI-NotI

Plasmid Map:



ACCN: NM_001001583

Insert Size: 1704 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001001583.1](#)

RefSeq Size: 2025 bp

RefSeq ORF: 1704 bp

Locus ID: 5152

UniProt ID: [O76083](#)

Cytogenetics: 21q22.3

Protein Families:	Druggable Genome
Protein Pathways:	Progesterone-mediated oocyte maturation, Purine metabolism
MW:	65.7 kDa
Gene Summary:	<p>The protein encoded by this gene catalyzes the hydrolysis of cAMP and cGMP to their corresponding monophosphates. The encoded protein plays a role in signal transduction by regulating the intracellular concentration of these cyclic nucleotides. Multiple transcript variants encoding several different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (18) lacks an alternate in-frame exon compared to variant 1. The resulting isoform (o) has the same N- and C-termini but lacks an internal segment compared to isoform a.</p>