

Product datasheet for **SC300254**

ARHGEF37 (NM_001001669) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ARHGEF37 (NM_001001669) Human Untagged Clone
Tag:	Tag Free
Symbol:	ARHGEF37
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC300254 representing NM_001001669.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGCCAAGCATGGAGCCGACGAGCCATCTCCAGGTCAGGGAGTCCGGACAGGGAAGTAGGGCCTCT
GAGGACAGATCGCTGCTTCATCAGAGGCTGGCTGTCGGGAGCTCATCGACACTGAGGTCTCCTACTTG
CACATGCTCCAGCTCTGTGCCTCTGACATCAGGAGCCGCCTCCAGCAGTTGCCGCAGGGAGATCTGGAT
GTCCTGTTCTCAAACATTGATGATATCATCAAAGTGAACAGCAGATTCTCCATGATCTGCAGGAGACA
GCCTCCAAGGAAGAGGAACAAGTGCAGCTAGTTGGTAACATATTTCTGGAATTCGAAGAGGAGTTGGAG
CAAGTCTATAAGGTCTACTGTGCCAGCTACGACCAGGCCTTGCTACTGGTGGACACGTACCGGAAGGAG
CCGGAGCTGCAGCGGCACATCCAGGGCATCGTTGAGGCGGTGGTCCGCAAGCTGGATCTTCAGGCCCTC
AGTTTCTGTGGTAATTCCTCTGCAGAGGATCACCAGGTACCCACTGCTGCTGCAGAAAATCCTGGAG
AACACAGTCCCTGATGCCAGTGCCTATCCTGTCCTTCAGAGGGCTGTCTGCCCCTCAGGACGTGAAC
ACCAATATCAATGAGTACAAGATGCGCAAGGAAGTGGCCTCCAAGTACACCAAGGTAGAGCAGCTGACC
CTCCGGGAGCGGCTGGCCCGCATCAACACACACACCCTCTCCAAGAAGACCACCCGGCTGAGCCAGCTG
CTGAAGCAGGAGGCGGGGCTGATCCCCAGGACAGAAGACAAGGAATTTGATGATTTAGAAGAGAGTTTC
CAGTGGGTGTCTGTGTGTGACTGAGCTGAAGAAACAGTGGCTGCTTACCTGGACAATCTGCAGGCT
TTCCTCTACTTACGGCCGACGAATAACAATCTGGACATCCCCGAGGGGCTGCAGTGCAGTATTGAAT
TTGGCAAGAGACCTTACCTTGAGGCCCTTCTGAAATTTAAGCAACGGCTAGAAGGCCTGGTGTGGCAG
CCACTGTGCAGCTGGCCAAAGCCCTGCTTGGCCCTCAGAACCTGATCAAGAAGCGTCTGGACAAGCTA
CTGGACTTTGAGCGGGTGAAGAGAAGCTGCTGGAGGTGGCAGTGTGACCTACCAGGAGGAGCCCGCC
CGGCACACATACCAGGCACTCAACTCGTGTAGTGGCTGAGCTCCACAGTTTAAACCAGTGGTCATG
CAGTGGCTGGGCCAGATCATGTGCACATTCGTGACCCTCCAGAGGGACCTTGCAAAGCAAGTGTGCAG
AGGGCAGAGGGAAGCATGGCCCAGCTGCCCCACCACCAGTCCCAGAGCCTGCCTTCAGGAAGCTGGTG
GAGGACGCACTGGGCCGACGAGTAACCAGCTTCGCTCCTTTCAAGAGACCTTTGAGAAAGTGCAGCCA
CCTCCCACCACACAACCCTCCTCCAGGGTCTGAACGCCAGGTGCAGGCTCTCCTGAGCAGGTATGGC
CCTGGGAAGCTGTACCAGGTGACAAGCAACATCAGTGGGACTGGGACTCTGGACCTGACTCTGCCTCGG
GGCCAAATCGTGGCCATCCTTCAAACAAGGACACCAAAGGCAACAGCGGCCGCTGGCTGGTGGACACC
GGGGGACATCGTGGGTATGTGCCGGCTGGGAAACTACAGCTGTACCATGTGGTCCCCAGTGCAGAGGAG
CTCAGAAGGCAGGCGGGGCTGAACAAAGACCCCGATGTCTAACACCGGAGCCCAGCCAGCTCTAGTG
CCCTCTATTTCCACCATGAACCAGGTATAGCCGCGTACCCTTTTGTGGCCAGAAGCAGCCATGAAGTG
AGCCTGCAGGCAGGCCAGCCTGTGACCATCCTGGAGGCCAGGACAAGAAGGGGAACCCTGAGTGGAGC
CTGGTGGAAAGTGAATGGACAGAGGGTTATGTGCCTTCTGGCTTCTTGGCCAGGGCTCGGAGCCAGTT
CTGTGGGGCTGGAGTCTGCCCTTAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001001669

Insert Size: 2028 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001001669.2
RefSeq Size:	4865 bp
RefSeq ORF:	2028 bp
Locus ID:	389337
UniProt ID:	A1IGU5
Cytogenetics:	5q32
MW:	76.3 kDa
Gene Summary:	May act as a guanine nucleotide exchange factor (GEF).[UniProtKB/Swiss-Prot Function]