

Product datasheet for **SC300349**

MTUS1 (NM_001001931) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MTUS1 (NM_001001931) Human Untagged Clone
Tag:	Tag Free
Symbol:	MTUS1
Synonyms:	ATBP; ATIP; ATIP3; ICIS; MP44; MTSG1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC300349 representing NM_001001931.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGACTGTACCAGGAGGATTTCCGAGTTGTACAGAACTGATATTTCTTCAAAAATCTTTATCAATTCT
ACACTCACCCACCGGCTGGTTCAGAGAGGCACTATGATGCTACCTTATTGACACTGCTGGTCGTGGGA
TCGTACAGCCTTTGTATAATTCTTTGTTAGCCAGTTTTACTGGGAAAAAACTGGTAATGCCGCTGTC
ATCAAATATGAGGAGAAAACCTCCAAAACAGCATTTCAGAAATGGTTCCTCAGGATCCTTTTATTTGAAG
CCTTTGGTATCCAGGGCTCATGTTCACTTGATGAAAACCTCTCCAAAAGGTCCTTCGAGAAAAAATTTA
TTTACAGCTCTTAATGCAGTTGAAAAGAGCAGGCAAAAGAATCCTCGAAGCTTATGTATCCAGCCACAG
ACAGCTCCCGATGCGCTGCCCTGAGAAAACACTTGAATTGACGCAATATAAAACAAAATGTGAAAAC
CAAAGTGGATTTATCCTGCAGCTCAAGCAGCTTCTTGCCTGTGGTAATACCAAGTTTGAGGCATTGACA
GTTGTGATTGAGCACCTGCTGTCTGAGCGGGAGGAAGCACTGAAACAACACAAAACCTATCTCAAGAA
CTTGTTAACTCCGGGGAGAGCTAGTCACTGCTTCAACCCTGTGAGAAATTAGAAAAAGCCAGGAAT
GAGTTACAACAGTGTATGAAGCATTCTGCCAGCAGCACCAGGCTGAAAAACAGAACGAGAGAATCGG
CTTAAAGAGTTTTACACCAGGGAGTATGAAAAGCTTCGGGACACTTACATTGAAGAAGCAGAGAAGTAC
AAAATGCAATTGCAAGAGCAGTTTGACAACCTAAATGCTGCGCATGAAACCTCTAAGTTGGAAATGAA
GCTAGCCACTCAGAGAACTTGAATTGCTAAAGAAGGCCATGAAGCCTCCCTTTTCAGAAATTAAGAAA
GGCCATGAAATAGAAAAGAAATCGCTTGAAGATTTACTTTCTGAGAAGCAGGAATCGCTAGAGAAGCAA
ATCAATGATCTGAAGAGTGAAAATGATGCTTTAAATGAAAATGAAATCAGAAGAACAAAAAGAAGA
GCAAGAGAAAAAGCAAATTTGAAAAATCCTCAGATCATGTATCTAGAACAGGAGTTAGAAAGCCTGAAA
GCTGTGTTAGAGATCAAGAATGAGAACTGCATCAACAGGACATCAAGTTAATGAAAATGGAGAACTG
GTGGACAACAACACAGCATTGGTTGACAAATTGAAGCGTTTTCCAGCAGGAGAATGAAGAATTGAAAGCT
CGGATGGACAAGCACATGGCAATCTCAAGGCAGCTTTCACGGAGCAGGCTGTTCTGCAAGAGTCGCTG
GAGAAGGAGTCGAAAGTCAACAAGCGACTCTATGGAAAACGAGGAGCTTCTGTGGAACTGCACAAT
GGGACCTGTGTAGCCCCAAGAGATCCCCCACATCCTCCGCCATCCCTTTGCAGTCACCAAGGAATTCG
GGCTCCTCCCTAGCCCCAGCATTTCACCCAGATGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001001931

Insert Size:	1554 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001001931.2
RefSeq Size:	4022 bp
RefSeq ORF:	1554 bp
Locus ID:	57509
UniProt ID:	Q9ULD2
Cytogenetics:	8p22
Protein Families:	Druggable Genome, Transcription Factors, Transmembrane
MW:	59 kDa
Gene Summary:	<p>This gene encodes a protein which contains a C-terminal domain able to interact with the angiotension II (AT2) receptor and a large coiled-coil region allowing dimerization. Multiple alternatively spliced transcript variants encoding different isoforms have been found for this gene. One of the transcript variants has been shown to encode a mitochondrial protein that acts as a tumor suppressor and participates in AT2 signaling pathways. Other variants may encode nuclear or transmembrane proteins but it has not been determined whether they also participate in AT2 signaling pathways. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (4), also known as ATIP4, lacks several 5' exons but has an alternate 5' exon, compared to variant 1. It encodes isoform 4 which has a much shorter and distinct N-terminus, compared to isoform 1.</p>