

Product datasheet for **SC300372**

FAM13C1 (FAM13C) (NM_001001971) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FAM13C1 (FAM13C) (NM_001001971) Human Untagged Clone
Tag:	Tag Free
Symbol:	FAM13C1
Synonyms:	FAM13C1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC300372 representing NM_001001971.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTTTTCTTGTCTGTTTCAGCCTTCAGGATAATTCTTCAGCAGCACCCTGTAACAGAGTGTGAC
GAAGATCCAGTCTCTACATGAAGACCAGACTGATTGCTCCAGTCTCAGAGATGAAACAATAAGAG
AACTACCCCGACGCAGGGGCTCTGGTAGAAGAGCACGCGCCGCCCTCTTGGGAGCCGCAGCAGCAGAAT
GTAGAGGCGACCGTGCTGGTGGACAGCGTATTGCGACCCAGCATGGGCAACTTCAAGTCCAGGAAGCCC
AAGTCCATCTTCAAAGCGGAGAGCGGGAGGAGCCACGGAGAAAGTCAGGAGACAGAGCATGTGGTATCC
AGCCAGTCAGAGTGTGAGGTGAGAGCAGGAACACCAGCTCATGAGAGTCCACAAAACAATGCCTTCAAG
TGCCAAGAAACAGTGCAGCTTCAACCAAGAATAGACCAGAGGACTGCCATTTGCGCAAAGGATGCTTTT
GAAACTCGGCAGGACTTAAATGAGGAAGAAGCTGCTCAGGTGCATGGAGTCAAGGACCCGGCGCCAGCA
TCAACCCAGAGCGTGCTTGCCGATGGGACAGATTCTGCAGACCCCTCACCAGTCCACAAAGATGGGCAG
AATGAGGCCGACAGTGCACCAGAAGACCTCCACTCTGTGGGACCAGCAGGCTGCTCTATCACATCACT
GATGGTGATAACCCACTGCTGTGCCACGATGCTCCATCTTCAGCCAAAGCCAGAGATTTCACTTAGAC
CCCGAGTCAGCCCCATCTCCACCCAGCACTCAGCAGTTTATGATGCCGCGGAGTTCTTACGCTGCAGC
TGTGGAGATGGCAAGGAGCCACAGACCATCACCCAGCTCACCAAGCACATCCAGAGCCTCAAGCGGAAA
ATTCGGAATTTGAAGAAAAATTTGAACAAGAAAAAGAAATACCGGTAACTAAGCAAGACAAGAACCTC
ATAAAGCCGCTTTATGACCGATACAGAATTATCAAGCAAACTTTGTCAACACCTTCCCTTATTCCAACA
ATTCAGGAGGAAGAGGACTCTGATGAAGACCGTCCACAGGAAGCAACAACCTTCTTTGGCAGATCCA
GCATCTCACCTTCTGTTGGTGACCACCTCACCTACTCTAATGAGACTGAGCCTGTTAGGGCCCTTTTA
CCAGATGAAAAGAAAGAGTAAAACCACAGCTCTCCATGTCTAATTTACATGAGGCTACCATGCCT
GTACTTCTTGACCATCTCCGAGAACTAGGGCTGACAAGAAGAGACTGCGGAAAGCCTTAAGAGAATTT
GAAGAACAGTTTTTTAAACAAACAGGAAGAAGTCCACAAAAGGAAGATAGGATACCAATGCAGATGAG
TATTATGAATATAAGCACATAAAAGCCAAACTGAGACTATTAGAGGTCTCATCAGCAAGCAAGATGTG
GCCAAACTATTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: Sgfl-MluI

ACCN: NM_001001971

Insert Size: 1464 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001001971.2](#)

RefSeq Size: 3101 bp

RefSeq ORF: 1464 bp

Locus ID: 220965

UniProt ID: [Q8NE31](#)

Cytogenetics: 10q21.1

MW: 55 kDa