

Product datasheet for **SC300376**

ATE1 (NM_001001976) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ATE1 (NM_001001976) Human Untagged Clone
Tag:	Tag Free
Symbol:	ATE1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC300376 representing NM_001001976.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCTTCTGGCGGGGGTTCGCCAGCGTCGTGGACTATTTCCCTAGCGAGGACTTCTACCGCTGC
GGCTACTGCAAGAACGAGTCGGGCAGCCGCTCCAATGGCATGTGGGCACATTCCATGACAGTACAGGAT
TATCAGGATCTCATAGACCGAGGATGGCGAAGAAGTGGAAAAATGTGTACAAACCTGTCATGAATCAA
ACATGTTGTCCTCAGTACACAATAAGGTGCCGACCTTTACAATTTAGCCTTCAAATCTCACAAGAAG
GTTTTGAAAAAATGTTGAAATTTCTAGCTAAAGGGGAGGTTCCCAAAGGAAGTTGTGAGGATGAGCCC
ATGGATCCACAATGGATGATGCTGTTGCGGGTGACTTTGCATTGATAAAATAACTGGATATACAGTGT
GATCTTAAACACTCAGTATGACATCAAAGAGAGTTTAGAGAGTGAAGGAAAAAATCAAAGAAAGAA
GAACCTCAGGAATTACTCAGTACAAGATTTCTGAGGAGAGAAGTTGGGCTCTGGTGAACCGTCACAT
TCAGTTAAAGTTCACACAGTTCTAAGCCAGGCAAAGGGGCTGATTTGAGTAAGCCTCCATGTCGAAAA
GCAAAGGAAATCCGGAAGAAAGGAAAGGTTAAACTAATGCAGCAGAACCCAGCTGGAGAAGTTGAG
GGTTTCCAGGCTCAAGGTCACCCACCATCTTTGTTTCCACCAAAGGCTAAATCCAACCCAGCCAAAATCA
CTCGAAGATTTAATTTTGTGACTTTACAGAGAATGCATCACACAAGTTAGAGGTGAGGGTGGTGAGA
TCATCTCCACCAAGTTCGCAGTTCAAAGCCACACTTCTGGAGTCTTACCAGGTCTATAAACGTTACCAG
ATGGTTATTCACAAGAACCACCTGATACGCCAACCGAAAGCCAGTTCACAAGATTCTTTGCGAGTTCA
CCCTTGGAGGCAGAGACTCCCCCTAATGGGCCAGATTGTGGCTATGGCTCCTTTACCAGCAGTACTGG
CTTGACGAAAGATCATTGCTGTGGGGTGATTGACATCCTCCAAACTGTGTATCATCTGTGATTTG
TACTACGATCCTGATTATTCGTTTTTGTCTTTGGGCGTACTCTGCACTACGAGAAATGCTTTTACT
AGGCAGCTTCATGAGAAAACTTCAACTCAGCTATTATTATATGGGTTTCTACATTCATTCATGTCCC
AAGATGAAATATAAGGGTCAGTATAGACCTTCTGATTTGCTGTGCCTGAGACATATGTTGGGTACCC
ATTGAGCAATGCCTGCCTTCACTTGAAAACCTCAAGTACTGCCGTTTCAACCAGGACCCAGAAGCAGTG
GATGAGGATCGCAGTACGGAACCTGACCGATTGCAGGTGTTTACAAGAGAGCCATCATGCCTTACGGT
GTTTATAAGAAACAGCAGAAAGACCCAAGTGAAGGAGGCTGCTGTTCTGCAGTACGCCAGCCTGGTGGG
CAGAAGTGCTCCGAGCGGATGCTGCTGTTTCAAGAACTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites:

SgfI-MluI

Plasmid Map:



ACCN:

NM_001001976

Insert Size:	1557 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001001976.2
RefSeq Size:	4930 bp
RefSeq ORF:	1557 bp
Locus ID:	11101
UniProt ID:	O95260
Cytogenetics:	10q26.13
MW:	59.1 kDa
Gene Summary:	<p>This gene encodes an arginyltransferase, an enzyme that is involved in posttranslational conjugation of arginine to N-terminal aspartate or glutamate residues. Conjugation of arginine to the N-terminal aspartate or glutamate targets proteins for ubiquitin-dependent degradation. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2013]</p> <p>Transcript Variant: This variant (1) encodes isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>