

Product datasheet for **SC300379**

USP16 (NM_001001992) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	USP16 (NM_001001992) Human Untagged Clone
Tag:	Tag Free
Symbol:	USP16
Synonyms:	UBP-M; UBPM
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC300379 representing NM_001001992.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGGGAAAAGAAACGGACAAAGGGAAAAACTGTTCCAATCGATGATTCTCTGAAACTTTAGAACCTGTG
TGCAGACACATTAGAAAAGGATTGGAACAAGGTAATTTGAAAAGGCTTTAGTGAATGTGGAATGGAAT
ATCTGCCAAGACTGTAAGACTGACAATAAAGTGAAGATAAAGCTGAAGAAGAAACAGAAGAAAAGCCT
TCAGTTTGGCTGTGTCTTAAATGTGCCATCAGGGCTGTGGCAGAAATTCTCAGGAGCAGCATGCCTTG
AAGCACTACTGACGCCAAGATCTGAACCTCACTGTCTGGTTCTTAGTTTGGACAACCTGGAGTGTATGG
TGTTACGTATGTGATAATGAGGTCCAGTATTGTAGTTCAAACAGTTGGGTCAAGTGGTTGATTATGTC
AGAAAACAAGCCAGCATTACAACCCAAAGCCAGAGAAAGATAATGAAAATATTGAACTTAAAAATAAA
AAATTAGAAAAGAGAGTAAGAATGAACAAGAGAGAGAAAAGAAGAAAACATGGCTAAAGAGAATCCT
CCCATGAATTCTCTTGGCAAATAACCGTGAAAGGACTCAGTAATTTGGGAAACACATGTTTCTTCAAT
GCAGTTATGCAGAACTTGTCAACAACACCAAGTCTTAGAGAACTACTAAAAGAAGTGA AAAATGTCTGGA
ACAATTGTAAAAATTGAACCACCTGATTTGGCATTAAACAGAACCTTAGAAAATAAACCTTGAGCCTCCA
GGCCCTCTTACTTTAGCCATGAGCCAGTTTCTTAAAGAGATGCAAGAGACCAAAAAGGGGTTGTGACA
CCGAAAGAACTCTTTTCTCAGGTCTGTAAAAAGCAGTGCCGTTTAAAGGCTATCAGCAGCAAGACAGC
CAGGAGCTGCTTCGCTACTTATTGGATGGGATGAGAGCAGAAGAACCAAAAGAGTGAGTAAAGGAATA
CTTAAAGCATTGGTAATTTCTACTGAAAAGTTGGATGAAGAACTAAAAATAAAGTTAAAGATTATGAG
AAGAAAAATCAATGCCAAGTTTGTGGACCGCATCTTTGGTGGTGAACCTAACTAGTATGATCATGTGT
GATCAATGCAGAACTGTCTCCTTGGTTCATGAATCTTTCCTTGATTTGTCCCTCCAGTTTTAGATGAT
CAGAGTGGTAAGAAAAGTGA AATGATAAAAACTGAAAAGACAGTGGAGGATGAAGATCAAGATAGT
GAGGAAGAAAAAGATAACGCAGATTACATAAAAGAGAGAAGTGATATTCCTTCTGGAACAAGTAAGCAC
TTACAGAAAAAGCAAAGAAACAAGCCAAAAAGCAAGCCAAAGAACCAACGAAGACAACAAAAAATCAA
GGAAAAGTTCTTCATTTAAATGATATTTGACTATTGACCATCCTGAAGACAGTGAATATGAAGCTGAA
ATGTCACTTCAAGGAGAAGTAAATATTAATCCAACCATATTTCAACAAGAGGGTGTATGCATAAAGAA
TATTGTGTCAACCAGAAAGATTTGAATGGCCAAGCAAAAATGATCGAAAGTGAACCTGACAATCAAAAA
TCCACAGAGGAAGTAGATATGAAAATATCAACATGGATAATGATCTGGAGGTTTAAACATCTTCTCCC
ACTAGGAATTTAAATGGTGCCTACCTAACGGAAGGGAGCAATGGAGAAGTGGACATTTCCAATGGTTTC
AAAAACCTAAATTTGAATGCTGCTCTTCATCCTGATGAAATAAATATAGAGATTCTGAATGATAGTCAT
ACTCCTGGAACAAGGTGTATGAGGTTGTAATGAAGATCCAGAACTGCTTTCTGTA CTCTTGCAAC
AGGGAAGTTTTCAATACTGATGAGTGTTCAATCCAACATTGTTTATATCAGTTCACCCGTAATGAGAAA
CTTCGAGATGCGAATAAACTGCTTTGTGAAGTATGCACACGGAGACAGTGAATGGACCAAAAGGCAAAAT
ATAAAAGGTGAAAGGAAGCATGTTTACACCAATGCCAAAAAGCAGATGCTAATTTCTCTTGCTCCTCCT
GTTCTTACTCTTCATTTAAAGAGATTTACAGCAGGCTGGTTTTAACCTACGCAAAGTTAACAAACACATA
AAGTTCCGGAAATCTTAGATTTGGCTCCTTTTGCACCCTTAAATGTAAGAATGTTGCAGAGAAAAT
ACAAGGGTACTCTATTCTTATATGGAGTTGTTGAACACAGTGGTACTATGAGGTGCGGGCATTACACT
GCCTATGCCAAGGCAAGAACCGCAAATAGTCATCTCTAATCTTGTTCACGGTGATATCCACAAA
GATTTTGAATGGAATCAAAAGGGCAGTGGTTTCACATCAGCGACACACATGTGCAAGCTGTGCCTACA
ACTAAAGTACTAACTCAAGCGTACCTCCTATTTTATGAGAGAATACTGTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: Sgfl-Mlul

Plasmid Map:


ACCN: NM_001001992

Insert Size: 2469 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001001992.1](#)

RefSeq Size: 2986 bp

RefSeq ORF: 2469 bp

Locus ID: 10600

UniProt ID: [Q9Y5T5](#)

Cytogenetics: 21q21.3

Protein Families: Protease

MW: 93.5 kDa

Gene Summary: This gene encodes a deubiquitinating enzyme that is phosphorylated at the onset of mitosis and then dephosphorylated at the metaphase/anaphase transition. It can deubiquitinate H2A, one of two major ubiquitinated proteins of chromatin, in vitro and a mutant form of the protein was shown to block cell division. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. [provided by RefSeq, Jul 2008]
Transcript Variant: This variant (2) uses an alternate in-frame splice site compared to variant 1. The resulting isoform (b) has the same N- and C-termini but is shorter compared to isoform a.