

Product datasheet for **SC300383**

GMPR2 (NM_001002000) Human Untagged Clone

Product data:


| | |
|---------------------------|--|
| Product Type: | Expression Plasmids |
| Product Name: | GMPR2 (NM_001002000) Human Untagged Clone |
| Tag: | Tag Free |
| Symbol: | GMPR2 |
| Synonyms: | GMPR 2 |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |
| Fully Sequenced ORF: | >SC300383 representing NM_001002000. Blue=Insert sequence Red=Cloning site Green=Tag(s) |

```
GCTCGTTT TAGTGAACCGTCAGAATTTTGT AATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCC GCGATCGCC
ATGCCTCATATTGACAACGATGTGAACTGGACTTCAAGGATGTCCTTTTGAGGCCAAACGCAGTACC
CTTAAGTCTCGAAGTGAGGTGGATCTCACAAGATCCTTTTCATTCGGAACCTAAAGCAGACATACTCT
GGGGTTCCCATCATTGCTGCCAATATGGATACTGTGGCCACCTTTGAGATGGCCAAGTTCTCTGTAAG
TTCTCTCTTCACTGCTGCCATAAGCACTATAGCCTCGTTCAGTGGCAAGAGTTTGTGGCCAGAAAT
CCTGACTGTCTTGAGCATCTGGCTGCCAGCTCAGGCACAGGCTCTTCTGACTTTGAGCAGCTGGAACAG
ATCCTGGAAGCTATTCAGGTGAAGTATATATGCCTGGATGTGGCAAATGGCTACTCTGAACACTTT
GTTGAATTTGAAAAGATGTACGGAAGCGCTTCCCCAGCACACCATCATGGCAGGGAATGTGTTAACA
GGAGAGATGGTAGAAGAGCTCATCTTTCTGGGGCTGACATCATCAAAGTGGGAATTTGGGCCAGGCTCT
GTGTGTAATACTCGAAGAAAAGTGGAGTGGGGTATCCACAGCTCAGCGCAGTGTGGAGTGTGCAGAT
GCTGCTCATGGCCTCAAAGGCCACATCATTTTCAGATGGAGGTTGCAGCTGTCTGGGGATGTGGCCAAG
GCTTTTGGGGCAGGAGCTGACTTCGTGATGCTGGTGGCATGCTGGTGGGCACAGTGAATCAGGTGGT
GAGCTCATCGAGAGGGATGGCAAGAAGTACAAGCTCTTCTATGGAATGAGTTCTGAAATGGCCATGAAG
AAGTATGCTGGGGCGTGGCTGAGTACAGAGCCTCAGAGGAAAGACAGTGGAAAGTTCCTTTAAAGGA
GATGTGGAACATACCATCCGAGACATCCTAGGAGGATCCGCTCTACGTGTACCTATGTGGGAGCAGCT
AAGCTCAAAGAGTTGAGCAGGAGAACTACCTTCATCCGAGTCAACCCAGCAGGTGAATCCAATCTTCAGT
GAGGCGTGCTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_001002000

Insert Size: 1047 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001002000.2](#)

RefSeq Size: 1989 bp

RefSeq ORF: 1047 bp

Locus ID: 51292

UniProt ID: [Q9P2T1](#)

Cytogenetics: 14q12

Protein Families: Druggable Genome

Protein Pathways: Purine metabolism

MW: 37.9 kDa

Gene Summary: This gene encodes an enzyme that catalyzes the irreversible and NADPH-dependent reductive deamination of guanosine monophosphate (GMP) to inosine monophosphate (IMP). The protein also functions in the re-utilization of free intracellular bases and purine nucleosides. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Apr 2017]

Transcript Variant: This variant (2) differs in the 5' UTR, lacks a portion of the 5' coding region, and initiates translation at a downstream in-frame start codon, compared to variant 1. The encoded isoform (2) has a shorter N-terminus than isoform 1. Variants 2, 3, and 4 encode the same isoform (2).