

Product datasheet for SC300427

ZFYVE27 (NM_001002262) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	ZFYVE27 (NM_001002262) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZFYVE27
Synonyms:	PROTRUDIN; SPG33
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC300427 representing NM_001002262. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCAGACATCAGAACGTGAGGGGAGTGGGCGGAGCTGAGCCCCAGCGTGATGCCCGAGGCTCCCTG
GAGTCTCCACCTTTCTACCAAGTCCCGAGCGTTTGACCTTTTCAACTTGGTTCTCTCTACAAGAGG
CTGGAGATCTACCTGGAACCCCTGAAGGATGCAGGTGATGGTGTTCGATACTTGCTCAGGTGCCAGATG
CCTTTGTGTTCTTGCTGACCTGCCTGGGCTCAACGTCTTGTTCCTCACTTTGAATGAGGGTGCATGG
TACTCAGTAGGTGCCCTGATGATTTAGTGCCCGCCTGCTGGGCTACCTTCAGGAGTTTGCCGGGCA
CGGCTGCCTGATCCGAGCTGATGCGGAGGAAGTATCATAGCGTGAGGCAGGAGGACCTGCAGAGAGGT
CGCCTGTCTCGTCCCGAGGCCGTGGCTGAGGTGAAGAGCTTCTTGATCCAGCTGGAGGCCTTCTGAGC
CGCCTGTGCTGCACATGTGAAGCCGCTACCGCGTGTGCACTGGGAGAACCCCGTGTGCTCCTACAG
TTCTATGGGGCTCTTCTGGGCACAGTCTGCATGCTGATTTGCTGCCACTCTGCTGGGTTCTCACCCCT
TTAAACAGCAGCCTCTTTCTGGGGAATGTGGAGTTCTTCCGAGTTGTGTCTGAGTACAGGGCATCTCTG
CAGCAGAGGATGAACCCAAAGCAGGAAGAGCATGCCTTTGAGAGTCTCCACCACAGATGTTGGGGGG
AAGGATGGTCTGATGGACAGCACGCTGCCCTCACACCCACGGAGGACCTCACACCGGCAGCGTGGAG
GAGGCTGAGGAGGCTGAGCCAGATGAAGAGTTTAAAGATGCGATTGAGGAGGATGATGAGGGCGCCCG
TGCCAGCAGAGGATGAGCTGGCCCTGCAGGACAACGGTTCTGAGCAAGAATGAGGTGCTGCGCAGC
AAGGTGTCTCGGCTCACGGAGCGGCTCCGCAAGCGCTACCCACCAACAACCTTCGGGAACCTGCACGGC
TGCTCGGCCACCTTCTCAGTGCTGAAGAAGAGGGCGGAGCTGCAGTAATTGTGAAACAGCTTCTGTCT
CGATGCTGCTCCTTCAAGGTGCCAAGTCTCCATGGGGCCACAGCCCTGAAGCCAGAGGGAGACT
GTGTTTGTGTGCCTCGTGAACCCAGACCTTGAGCAAGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-MluI



[View online >](#)

Plasmid Map:


ACCN: NM_001002262

Insert Size: 1215 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001002262.3](#)

RefSeq Size: 3038 bp

RefSeq ORF: 1215 bp

Locus ID: 118813

UniProt ID: [Q5T4F4](#)

Cytogenetics: 10q24.2

Protein Families: Transmembrane

MW: 45.1 kDa

Gene Summary: This gene encodes a protein with several transmembrane domains, a Rab11-binding domain and a lipid-binding FYVE finger domain. The encoded protein appears to promote neurite formation. A mutation in this gene has been reported to be associated with hereditary spastic paraplegia, however the pathogenicity of the mutation, which may simply represent a polymorphism, is unclear. [provided by RefSeq, Mar 2010]
Transcript Variant: This variant (3) differs in the 5' UTR, uses an alternate in-frame splice site in the middle portion of the coding region, and lacks an in-frame alternate exon in the 3' coding region, compared to variant 1. This results in a shorter protein (isoform c), compared to isoform a.