

Product datasheet for **SC300455**

PIB5PA (INPP5J) (NM_001002837) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PIB5PA (INPP5J) (NM_001002837) Human Untagged Clone
Tag:	Tag Free
Symbol:	INPP5J
Synonyms:	INPP5; PIB5PA; PIPP
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC300455 representing NM_001002837.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGAGGGCCAGAGCAGCAGGGGCAGCAGGAGGCCAGGGACCCGGGCTGGCCTGGGTTCCCTGCCCATG
CCCCAGGGTGTGGCCAAACTGGGGCACCCCTCCAAGGTGGACTCAAGTTTTTCAGCTCCAGCAAAGAAG
AACGCAGCCCTAGGACCCTCGGAACCAAGGATCACTGTGGTCACATGGAACGTGGGCACTGCCATGCC
CCAGACGATGTCACATCCCTCCTCCACCTGGGCGGTGGTGACGACAGCGACGGCGAGACATGATCGCC
ATAGGGTTGCAGGAAGTGAACCTCATGCTCAACAAGCGACTCAAGGACGCCCTCTTCACGGACAGTGG
AGTGAGCTGTTTATGGATGCGCTAGGGCCCTTCAACTTCGTGCTGGTGAGTTCGGTGAGGATGCAGGGT
GTCATCCTGCTGCTGTTCCGCAAGTACTACCACCTGCCCTTCTGCGAGACGTGCAGACCGACTGCACG
CGCACTGGCCTGGGCGGCTACTGGGTAACAAGGGTGGCGTGAGCGTGCCTGGCGCCCTTCGGGCAC
ATGCTCTGCTTCTGAAGTGCCTTGCCTGCGCATATGGACAAGGCGGAGCAGCGCAAAGACAACCTTC
CAGACCATCCTCAGCCTCCAGCAGTTCGAAGGGCCGGGCGCACAGGGCATCCTGGATCATGACCTCGTG
TTCTGGTTCGGGGACCTGAACCTCCGCATTGAGAGCTATGACCTGCACCTTTGTCAAGTTTCCCATCGAC
AGTGACCAGCTCCATCAGCTCTGGGAGAAGGACCAGCTCAACATGGCCAAGAACCTGGCCCATCTTG
AAGGGCTTTCAGGAGGGGCCCTCAACTTCGCTCCACCTTCAAGTTTGTATGTGGGTACCAACAAATAC
GATACCAAGTGCAGAAACCGGAAGCCAGCTTGGACAGACCGTATCCTATGGAAGGTCAAGGCTCCAGGT
GGGGTCCCAGCCCTCAGGACGGAAGACCCAGCTCCAGGTGACGCAGCACAGCTACCGCAGCCAC
ATGGAATACACAGTCAGCGACCACAAGCCTGTGGCTGCCAGTTCCTCCTGCAGTTTGCCTTCAGGGAC
GACATGCCACTGGTGGGCTGGAGGTGGCAGATGAGTGGTGGCGCCGAGCAGGCGGTGGTGAGGTAC
CGCATGGAAACAGTGTTCGCCCGCAGCTCCTGGGACTGGATCGGCTTATACCGGGTGGGTTCCGCCAT
TGCAAGGACTATGTGGCTTATGTCTGGGCCAAACATGAAGATGTGGATGGGAATACCTACCAGGTAACA
TTCAGTGAGGAATCACTGCCCAAGGGCCATGGAGACTTCACTCCTGGGCTACTATAGTCACAACCACAGC
ATCCTCATCGGCATCACTGAACCCTTCCAGATCTCGCTGCCTTCTCGGAGTTGGCCAGCAGCAGCACA
GACAGCTCAGGCACCAGCTCAGAGGGAGAGGATGACAGCACACTGGAGCTCCTTGACCCAAAGTCCCGC
AGCCCCAGTCTGGCAAGTCCAAGCGACACCCGACCCGAGCCCGGGACTGGCCAGGTTCCCTGGGCTT
GCCCTACGGCCCTCATCCCGTGAACGCCGTGGTCCAGCCGTAGCCCTCACCCAGAGCCGCCCTG
TCCCGAGTGGCTCCTGACAGGAGCAGTAATGGCAGCAGCCGGGCGAGTAGTGAAGAGGGGCCCTCTGGG
TTGCCTGGCCCTGGGCCTTCCCACCAGCTGTGCCTCGAAGCCTGGGCCTGTTGCCCGCCTTGCCTA
GAGACTGTAGACCCTGGTGGTGGTGGCTCCTGGGACCTGATCGGGAGGCCCTGGCGCCCAACAGCCTG
TCTCCTAGTCCCCAGGGCCATCGGGGGCTGGAGGAAGGGGGCCTGGGGCCCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001002837

Insert Size: 1917 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001002837.2](#)

RefSeq Size: 2303 bp

RefSeq ORF: 1917 bp

Locus ID: 27124

UniProt ID: [Q15735](#)

Cytogenetics:	22q12.2
Protein Families:	Druggable Genome
Protein Pathways:	Inositol phosphate metabolism, Metabolic pathways, Phosphatidylinositol signaling system
MW:	70.2 kDa
Gene Summary:	<p>Inositol 5-phosphatase, which converts inositol 1,4,5-trisphosphate to inositol 1,4-bisphosphate. Also converts phosphatidylinositol 4,5-bisphosphate to phosphatidylinositol 4-phosphate and inositol 1,3,4,5-tetrakisphosphate to inositol 1,3,4-trisphosphate in vitro. May be involved in modulation of the function of inositol and phosphatidylinositol polyphosphate-binding proteins that are present at membranes ruffles (By similarity).[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (3) lacks an alternate in-frame segment in the 5' coding region, compared to variant 1, resulting in an isoform (c) that is shorter than isoform a.</p>