

Product datasheet for SC300488

POTEA (NM_001002920) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	POTEA (NM_001002920) Human Untagged Clone
Tag:	Tag Free
Symbol:	POTEA
Synonyms:	A26A1; CT104.3; POTE-8; POTE8
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC300488 representing NM_001002920. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGTGGCTGAGGTTTCTCAAAGCTGGCTGCCTCCCCTATGAAGAAGCCCTTTGGCTTCAGGGGCAAG
ATGGGCAAGTGGTGTCTGCTGCTGCTTCCCCTGCTGCAGGGGAAGCGGCAAGAACAACATGGGTGCTTGG
AGAGACCACGACGACAGCGCCTTACGGAGCCAAGGTATCACGTCCGTCGAGAAGATCTGGGCAAGCTC
CACAGAGCTGCCTGGTGGGTGAAGTCCCAGAGCGGATCTCATCGTCATGCTCAGGGGCCCTGGCATT
AACAAAAGGGACAAGAAGAAGAGGACTGCTCTACACTTGGCCTGTCCAATGGCAATTCAGAAGTAGTA
AGTCTTCTGCTGGACAGACAATGTCAATTCATGTCTTTGACAGCAAAAAGAGGACAGCTCTGATAAAG
GCCGTACAATGCCAGGAGGATGAATGTGCGTTAATGTTGCTACAACATGGCACTGATCCAATCTTCCA
GATATGTATGGGAATACTGCTCTACACTATGCAGTCTACAATGAAGATAAATTAATGGCCAAAACACTG
CTTTTATACGGTGTGACATTGAATCAAAAAACAAGGGTGGCCTCACACCACTTTTGCTTGTGATCAT
GGCAAAAACAGCGAATGGTGAATTTTAATCAAGAAAAAGCTAATTTAAATGCACCTTGATAGGTTT
GGAAGAATTTGCCAATTAATTTCTGACTACAAAGAAAACAGATGCCAAACAACCTTTCTGGAAATAGC
AATCCAGAACAAGACTTAAAGCTGACATCAGAGGAAGGCCACAAAGGCTTAAAGGAAGTGAAAAATG
CAGCATGAGAAAGTACTCAAGAACCAGACATAAATAAGGATTGCGATAGAGAGGTTGAAGAAGAAATG
CAGAAGCATGGAAGTAATATGTGGGATTATCAGAAAACCTGACTGATGGTGTGCTGCTGGCAATGGT
GATGGTGGATTAGTTCCACAAAGAAAGAGCAGAAAGCATGAAAAATCAGCAATTTCCCTAACACAGAGATT
GAAGAGTATCACAGGCCGTAAGAAAGAAATCTAATGAAAAGAAACAAGGTCAAAGCCAAATACATTCAGTG
GATAACCTTGATGACATAACATGGCCATCAGAAATAGCGTCAGAGGATTATGATTTGCTTTTCTCTAAT
TATGAGACTTTTACATTGCTCATTGAACAACCTCAAAATGGATTTTAAACGATTCTGCTAGCCTATCAAAA
ATCCAGGATGCAGTTATTTCAGAAGAACAACCTACTAGAACTCAAAAATAGCCACTATGAGCAACTTACA
GTAGAAGTTGAACAAATGAAAAATATGGTTCATGTACTACAAAAATAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```



View online »

Restriction Sites: Sgfl-MluI

Plasmid Map:



ACCN: NM_001002920

Insert Size: 1359 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001002920.1](#)

RefSeq Size: 1838 bp

RefSeq ORF: 1359 bp

Locus ID: 340441

UniProt ID: [Q6S8I7](#)

Cytogenetics: 8p11.1

MW: 51.3 kDa