

Product datasheet for **SC300526**

ARHGEF35 (NM_001003702) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ARHGEF35 (NM_001003702) Human Untagged Clone
Tag:	Tag Free
Symbol:	ARHGEF35
Synonyms:	ARHGEF5L
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC300526 representing NM_001003702.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGAGGCTGAGGAGGCCAGCATGGAGCCTCTCTCCCATCTGCCATAGAGGAATTCAGCATTATC
CCTGAGGCTCCCATGAGGAGCAGCCAGGTCTCTGCCTTGGGGCTTGAAGCTCAAGAAGATGAGGACCA
TCCTATAAGTGGAGAGAGGAACACAGACTCTCAGCAACTCAGCAGAGTGAGTTAAGGGATGTGTGTGAC
TATGCGATTGAGACGATGCCCTCTTTTCCCAAGGAAGTTCTGCAGATGTGGAGCCCAATCAGGAAAGC
CTTGTGGCTGAGGCTGTGACACTCCGGAACACTGGGAGGCAGTACCCAGAGCCTAGCAGGCCGACAA
GCAAGGACTCTAGCTCCCCAGAGCTCTGGGCTGCCCATTCAGAGTGAGCATCTAGACATGGCCCCA
TTTTCCAGTGACCTGGGAAGCGAAGAAGAGGAGGTGGAATTTGGCCAGGACTTACTTCTTTGACATTG
GGATCTGGACAGGCAGAAGAAGAAGAGGAAACCTTTCAGATAACTCTGGTCAGACCAGATATTATCT
CCCTGCGAAGAGCATCTGCAGAGACCAACCAGAATGAAGGCGCTGAAAGTGGGACTATCAGGCAGGGG
GAAGAGCTGCCATCTGAGGAGCTGCAGGAAAGTCAAGGGCTCTTGCATCCCAGGAGGTCCAAGTTCTG
GAGGAGCAGGGACAGCAGGAAGCAGGATTTCCGGGGGAAGGAAGTCTGAGGGAGGATGTTTGTGCCGAT
GGGCTATTAGGGGAGGAACAGATGATAGAGCAGTTAATGATGAAAAGGGAGAACAGAAGCAAAACAG
GAACAGGTACAAGATGTGATGCTTGGGAGACAAGGAGAAAGATGGGGCTCACTGGGAGCCAGAGGGT
CTGAATGACGGTGAGTGGGAGCAGGAGGATATGGAGAGGAAGGCTCAGGGTCAGGGAGGTCCAGAACAG
GGAGAAGAGAGGAAGAGGAGCTGCAGGTGCCAGAAGAGAACAGGGCGGACTCTCAGGACGAAAAGAGT
CAAACCTTTTGGGAAAATCAGAGGAAGTAAGTGAAGCAAGAAGATCATGGTATAAAGGAGAAAGGG
GTCCCAGTCAGCGGCAGGAGGCGAAAGAGCCAGAGAGTTGGGATGGGGGCAGGCTGGGGCAGTGGGA
AGAGCGAGGAGCAGGGAAGAGGAGAATGAGCATCATGGGCTTCAATGCCGCTCTGATAGCCCCTGAG
GACTCTCCTCACTGTGACCTGTTTCCAGGTGCCTCATATCTCGTGACTCAGATTCCCGGACTCAGACA
GAGTCCAGGGCTGAGGAAGTGTCCCCGCACTCTGTCTCCCTTGCTAGAGCCCATCAGATGCTCTCAC
CAGCCCATTCTCTACTGGGCTCCTTTTACTGAGGAGTCACCTGACAAGGAAAACTTCTATCAGTA
CTTTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001003702

Insert Size: 1455 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001003702.2](#)

RefSeq Size: 2474 bp

RefSeq ORF: 1455 bp

Locus ID: 445328

UniProt ID: [A5YM69](#)

Cytogenetics: 7q35

MW: 53.3 kDa