

Product datasheet for **SC300538**

KLHL31 (NM_001003760) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	KLHL31 (NM_001003760) Human Untagged Clone
Tag:	Tag Free
Symbol:	KLHL31
Synonyms:	bA345L23.2; BKLHD6; KBTBD1; KLHL
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC300538 representing NM_001003760.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGCACCCAAAAAGAAGATTGTCAAAAAGAACAAGGAGATATCAATGAGATGACTATAATCGTAGAA
GATAGCCCCCTAAACAACTGAATGCTTTGAATGGGCTCCTAGAGGGAGCAATGGCCTTAGCTGCATT
TCTTCTGAACAAACAGATGCTTCTTATGGCCCCAACCTCTTGAAGGTTTAAAGTAAATGCGGCAGGAG
AACTTCTATGTGACTTAGTCATTGGTACCAAAACCAAATCCTTTGATGTTTCATAAGTCAGTCATGGCT
TCATGCAGTGAGTATTTTACAACATCCTAAAAAAGACCCGTC AATTCAGAGGGTGGATCTCAATGAT
ATCTCACCCTAGGCTGGCCACTGTCATTGCATATGCCACTGAAAAGCTCACTCTCTCCTTGTAT
ACAATAGGAAGCATTATTTCTGCTGCTGTTTATCTTCAGATCCATACTCTGTAAAGATGTGCAGTGAT
TTTCTGATACGGGAGATGAGTGTGAGAATTGCATGTATGTTGTTAATATTGCTGAAACATACTCCCTA
AAAAATGCAAAGCAGCAGCCAGAAATTTATTCGGGATAACTTCCTTGAATTTGCAGAATCGGATCAG
TTTATGAACTTACATTTGAACAAATTAATGAACTTCTTATAGATGATGACTTACAGTTGCCTTCTGAG
ATAGTAGCATTCCAGATTGCAATGAAATGGTTAGAATTTGACCAAAAGAGAGTAAATACGCTGCAGAT
CTTTTGAGCAATATTCGCTTTGGTACCATCTCTGCACAAGACCTGGTCAATTATGTTCAATCCGTACCA
AGAATGATGCAAGATGCTGATTGTCACAGACTTCTCGTAGATGCTATGAACTACCACTTGCCTCCATAT
CATCAAAACACATTGCAATCTAGGCGAACAAGAATCCGAGGTGGCTGCCGAGTCCCTGCTCACTGTTGGG
GGACGCCCAGGCCTTACTGAGAAGTCCCTTAGCAGAGACATCTTGATAGAGACCCTGAAAATGGATGG
AGCAAGCTTACGAAAATGCCAGCCAAAAGTTTTAATCAGTGTGTGGCTGTGATGGATGGATTTCTTTAT
GTAGCCGGTGGTGAAGACCAGAATGATGCAAGAAATCAAGCCAAGCATGCAGTCAGCAATTTCTGCAGA
TAGCATCCCCGCTTCAACACCTGGATACACCTGGCCAGCATGAACCAGAAGCGCACGCACTTCAGCCTG
AGCGTGTTCAACGGGCTCGTGTACGCCGCGGGCGCCGCAACGCAGAAGGAAGCCTGGCCTCGCTGGAG
TGCTACGTGCCCTCCACCAATCAGTGGCAGCCGAGACGCCCTGGAGGTGGCGCCTGCTGCCACGCT
AGCGCGGTGCGCGACGGCCGCTGCTGGTGACCGGAGGCTACATAGCCAACGCTACTCGCGCTCTGTG
TGCGCCTACGACCCGGCCAGCGACTCGTGGCAGGAGCTGCCGAACCTCAGCACACCCGGGGCTGGCAC
TGCGCGGTACGCTGAGCGACAGAGTGTACGTGATGGCGGCAGCCAGCTGGGGCCGCGCGGGGAGCGC
GTGGACGTGCTACCGTGGAGTGTACAGCCCCGCGACCCGCCAGTGGAGCTACGCGGCCGCTGCAG
GTGGGAGTGAGCACTGCGGGCTCTCGGCCTGCATGGCCGCGCCTACCTGGTGGGGGGTGAACGAG
GGCGAGAAGAAGTACAAGAAGTGCATCCAGTGCTTCAGCCCCGAGCTCAACGAGTGGACGGAGGACGAC
GAGCTACCCGAGGCCACTGTCCGGCTGCTGCTGCACCTCTCGATGCCAACAACGTGACTCGGGAA
TCCCGGGCCAGTTCGGTATCTTCTGTGCCAGTCAGTATCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: Sgfl-MluI

ACCN: NM_001003760

Insert Size: 1905 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001003760.4](#)

RefSeq Size: 5743 bp

RefSeq ORF: 1905 bp

Locus ID: 401265

UniProt ID: [Q9H511](#)

Cytogenetics: 6p12.1

Protein Families: Transmembrane

MW: 70.2 kDa

Gene Summary: Transcriptional repressor in MAPK/JNK signaling pathway to regulate cellular functions. Overexpression inhibits the transcriptional activities of both the TPA-response element (TRE) and serum response element (SRE).[UniProtKB/Swiss-Prot Function]