

Product datasheet for SC300541

LYK5 (STRADA) (NM_001003787) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LYK5 (STRADA) (NM_001003787) Human Untagged Clone
Tag:	Tag Free
Symbol:	STRADA
Synonyms:	LYK5; NY-BR-96; PMSE; Stlk; STRAD; STRAD alpha
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC300541 representing NM_001003787. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTCATTTCTTGTAAAGTAAACCAGAGCGAATCAGGCGGTGGGTCTCGGAAAAGTTCATTGTTGAGGGC
TTAAGAGATTTGGAATATTTGGAGAGCAGCCTCCGGGTGACACTCGGAGAAAAACCAATGATGCGAGC
TCAGAGTCAATAGCATCCTTCTCTAAACAGGAGGTGATGAGTAGCTTTCTGCCAGAGGGAGGGTGTAC
GAGCTGCTCACTGTGATAGGCAAAGGATTTGAGGACCTGATGACTGTGAATCTAGCAAGGTACAAACCA
ACAGGAGAGTACGTGACTGTACGGAGGATTAACCTAGAAGCTTGTCCAATGAGATGGTAACATTCTTG
CAGGGCGAGCTGCATGTCTCAAACCTCTTCAACCATCCCAATATCGTGCCATATCGAGCCACTTTTATT
GCAGACAATGAGCTGTGGTTGTACATCATTCATGGCATAACGGTTCTGCAAAAGATCTCATCTGTACA
CACTTCATGGATGGCATGAATGAGCTGGCGATTGCTTACATCCTGCAGGGGGTGTGAAGGCCCTCGAC
TACATCCACCACATGGGATATGTACACAGGAGTGTCAAAGCCAGCCACATCCTGATCTCTGTGGATGGG
AAGGTCTACCTGTCTGGTTTGGCAGCAACCTCAGCATGATAAGCCATGGGCAGCGGCAGCGAGTGGTC
CACGATTTTCCCAAGTACAGTGTCAAGTCTGCCGTGGCTCAGCCCCGAGGTCTCCAGCAGAATCTC
CAGGGTTATGATGCCAAGTCTGACATCTACAGTGTGGGAATCACAGCCTGTGAATGGCCAAACGCCAT
GTCCCCTTAAGGATATGCCTGCCACCCAGATGCTGCTAGAGAACTGAACGGCAGTGCCTGCCCTGCCTG
TTGGATACCAGCACCATCCCCGCTGAGGAGCTGACCATGAGCCCTTCGCGCTCAGTGGCCAACTCTGGC
CTGAGTGACAGCCTGACCACAGCACCCCGGCCCTCCAACGGTGACTCGCCCTCCACCCCTACCAC
CGAACCTTCTCCCCCACTTCCACCCTTTGTGGAGCAGTGCCTTCAGCGCAACCCGGATGCCAGGCC
AGTGCCAGCACCTCCTGAACCACTCTTTCTTCAAGCAGATCAAGCGACGTGCCTCAGAGGCTTTGCC
GAATTGCTTCGTCTGTACCCCATCACCAATTTTGGGGCAGCCAGTCTCAGGACCACAGTGGAAATC
TTTGGCCTGGTAACAACTGGAAGAGCTGGAGGTGGACGATTGGGAGTCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online »](#)

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001003787

Insert Size: 1296 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001003787.2](#)

RefSeq Size: 2367 bp

RefSeq ORF: 1296 bp

Locus ID: 92335

UniProt ID:	Q7RTN6
Cytogenetics:	17q23.3
Protein Families:	Druggable Genome, Protein Kinase
Protein Pathways:	mTOR signaling pathway
MW:	48.4 kDa
Gene Summary:	<p>The protein encoded by this gene contains a STE20-like kinase domain, but lacks several residues that are critical for catalytic activity, so it is termed a 'pseudokinase'. The protein forms a heterotrimeric complex with serine/threonine kinase 11 (STK11, also known as LKB1) and the scaffolding protein calcium binding protein 39 (CAB39, also known as MO25). The protein activates STK11 leading to the phosphorylation of both proteins and excluding STK11 from the nucleus. The protein is necessary for STK11-induced G1 cell cycle arrest. A mutation in this gene has been shown to result in polyhydramnios, megalencephaly, and symptomatic epilepsy (PMSE) syndrome. Multiple transcript variants encoding different isoforms have been found for this gene. Additional transcript variants have been described but their full-length nature is not known. [provided by RefSeq, Sep 2009]</p> <p>Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (1).</p>