

Product datasheet for **SC300560**

DLGAP1 (NM_001003809) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	DLGAP1 (NM_001003809) Human Untagged Clone
Tag:	Tag Free
Symbol:	DLGAP1
Synonyms:	DAP-1; DAP-1-ALPHA; DAP-1-BETA; DAP1; DLGAP1A; DLGAP1B; GKAP; SAPAP1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC300560 representing NM_001003809.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAACTTAATTTCCATAAAGACATTCTGTTTGGCATTCCAGCTAATAAGGTTCCACAAGATGAATGG
ACAGGGTACACCCACGAGGTAAGATGATGAAATTCATGCCGAAGAATGCGGAGTGCCAGTTATATC
AAGGCCATGGGGATGAAGACAGTGGAGACTCAGACACGAGTCTAAGCCTTCTCCAAAAGTTGCTGCG
CGGAGAGAAAGCTATCTCAAGGCTACTCAGCCATCCCTTACAGAACTCACCACACTCAAAAATCTCCAAT
GAACACTACCCAACTCCAGATCCGGAGTCATAGTTACCTGAGGGCAGTGAGTGAAGTCTCCATCAAC
CGGAGCCTGGACAGCCTGGACCTGCAGGCTTGCTCACATCACCAAAGTCCGCTCCAGGAATGAGAGC
TACATGCGAGCCATGAGCACCATCAGCCAGGTGAGCGAGATGGAAGTGAACGGGCAGTTCGAGTCCGTG
TGCGAGTCCGTGTTGAGCGAGTGGAGTGCAGGCCGTTGAAGCGCTGGACCTGCCATGCCCGGCTGC
TTCCGCATGCGGAGCCACAGCTATGTGCGGGCCATTGAGAAAGGCTGCTCCAGGACGACGAGTCCGTG
TCCCTGAGGTCGTCTCGCCGCCGCGCACCACCACCCTTAGGACCATCCAGAGCAGCAGGTTGCA
TCTTGATTACAACATATAAGAAGACACCCTCCAGTCCACCCAGAATACCAGAAAACCTTTTCAAT
TCTATCACAGCCAGAGTAGCACAGAGTCAAGCCAGGATGCCTACATGGACGGACAGGGCCAGCGAGGA
GATATTATCAGCCAGTCTGGACTCAGCAACTCCACCGAGAGCCTGGACAGTATGAAGGCTCTGACAGCC
GCCATCGAAGCTGCAAACGCCAGATCCATGGCCCTGCCAGTCAACACATGGGCAATAAAGTCCACC
GTCACCACCAGACTACCATAGCCACCGTACCACGGAGGACAGGAAGAAGGACCACTTTAAGAAAAAT
CGATGCCTGTCTATCGGGATACAGGTGGATGATGCTGAAGAACCTGACAAAACAGGGGAGAATAAAGCA
CCCAGTAAGTTCAGTCCGTGGGAGTGAAGTGAAGAAGAGAAGTCTCCGCAGGTTCACTCGATCC
AACAGTGTGACGACAGCAGTACAGGCCGACCTGGACTTCCATGATAATCTGGAATAATCTGGAATCT
ATAGAGGACAATTCGTGTCTGGCCCATGGCCAGACAGTCTCCCGCATGCCAGCACCTCCACAGTC
AGCATTACAGGGCTCAGGAAACCATTACCATGCCTGTGCCGCCGATGATGACTTTGACACGGATTTGAC
CCCTCTATTCTGCCTCTCCGGACCCCTGGATTGACTCTATCACTGAAGACCCTCTGGAGGCCGTGCAA
AGGTGAGTGTGCCACCGGATGGCCACTGGTTCCTGAAGCTTCTCCAGGCAGAGCGAGACCCGATGGAG
GGGTGGTGTCAACAGATGGAGCGGGAAGAACGGGAAAACAACCTGCCGAAGACATTCTAGGAAAAATC
CGAACCAGTGGGCAGTGCCCACTTCTCATGGCCAGAAATCTACCAGTTCAGAGAAGTGTGTGAA
GAAAACCTGAATCCTAATGCTCATCAAGACCCACCTCCAGGATTTGGCGGGTTTTGGGACATGCTG
CAGTTGTCCATAGAAAATATTAGTATGAAATTTGATGAACTTCATCAGTTAAAGGCAATAATTGAAA
CAGATGGATCCTCTTGACAAGAAGGAGAGAAGGGCCCTCCTCCAGTGCCAAAGAAGCCGCGAAGGGC
CCCGCGCCGCTGATCCGGGAGCGCTCGCTGGAGAGCTCGCAGCGCCAGGAGGCCGCAAGCGCTGATG
GCCGCAAGCGCGCCGCTCCGTCCGCCAGAACTCGGCCACCGAGAGCGCCGAGAGCATCGAGATCTAC
ATCCCCGAGGCGCAGACCCGGCTCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001003809

Insert Size: 2028 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001003809.2](#)

RefSeq Size: 5271 bp

RefSeq ORF: 2028 bp

Locus ID: 9229

UniProt ID: [O14490](#)

Cytogenetics: 18p11.31

MW: 75.3 kDa

Gene Summary: Part of the postsynaptic scaffold in neuronal cells.[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (2) differs in the 5' UTR, lacks a portion of the 5' coding region, and initiates translation at an alternate start codon, compared to variant 1. The encoded isoform (2) has a distinct N-terminus and is shorter than isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.