

Product datasheet for **SC300566**

RNF21 (TRIM34) (NM_001003827) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RNF21 (TRIM34) (NM_001003827) Human Untagged Clone
Tag:	Tag Free
Symbol:	TRIM34
Synonyms:	IFP1; RNF21
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC300566 representing NM_001003827.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCTTCAAAAATCTTGCTTAACGTACAAGAGGAGGTGACCTGTCCCATCTGCCTGGAGCTGTTGACA
GAACCTTGAGTCTAGACTGTGGCCACAGCCTCTGCCAGCCTGCATCACTGTGAGCAACAAGGAGGCA
GTGACCAGCATGGGAGGAAAAAGCAGCTGTCTGTGTGGTATCAGTTACTCATTGAAACATCTACAG
GCTAATCAGCATCTGGCCAACATAGTGGAGAGACTCAAGGAGGTCAAGTTGAGCCAGACAAATGGGAAG
AAGAGAGATCTCTGTGATCATCATGGAGAGAACTCCTACTCTTCTGTAAGGAGGATAGGAAAGTCATT
TGCTGGCTTTGTGAGCGGTCTCAGGAGCACCGTGGTCACCACACAGTCTCACGGAGGAAGTATTCAAG
GAATGTCAGGAGAACTCCAGGCAGTCTCAAGAGGCTGAAGAAGGAAGAGGAGGAAGCTGAGAAGCTG
GAAGCTGACATCAGAGAAGAGAAAACCTCTGGAAGTATCAGGTACAACTGAGAGACAAAGGATACAA
ACAGAATTTGATCAGCTTAGAAGCATCCTAAATAATGAGGAGCAGAGAGAGCTGCAAAGATTGGAAGAA
GAAGAAAAGAAGACGCTGGATAAGTTTGCAGAGGCTGAGGATGAGCTAGTTCAGCAGAAGCAGTTGGTG
AGAGAGCTCATCTCAGATGTGGAGTGTCCGAGTCAAGTGGTCAACAATGGAGCTGCTGCAGGACATGAGT
GGAATCATGAAATGGAGTGAGATCTGGAGGCTGAAAAAGCAAAAATGGTTTCCAAGAACTGAAGACT
GTATTCATGCTCCAGATCTGAGTAGGATGCTGCAAAATGTTTAGAGAAGTACAGCTGTCCGGTGCTAC
TGGGTGGATGTCACACTGAATTCAGTCAACCTAAATTTGAATCTTGTCTTTGAGAAGATCAGAGACAA
GTGATATCTGTGCCAATTTGGCCTTTTCAAGTGTATAATTATGGTGTCTTGGGATCCCAATTTTCTCC
TCTGGGAAACATTACTGGGAAGTGGACGTGTCCAAGAAAACCTGCCTGGATCCTGGGGTATACTGTAGA
ACATATCCCCCCATATGAAGTATGTTGTTAGAAGATGTGCAAAATCGTCAAAATCTTTACACCAAATAC
AGACCTCTATTTGGCTACTGGGTTATAGGTTACAGAATAAATGTAAGTATGGTGTCTTTGAAGAGTCT
TTGTCCTCTGATCCCGAGGTTTTGACTCTCTCCATGGCTGTGCCTCCCTGCCGTGTGGGTTTTCCCTC
GACTATGAAGCAGGCATTGTCTCATTTTTCAATGTACAAAGCCATGGCTCCCTCATTACAAAGTTCTCT
AAATGTTGCTTTTCTCAGCCTGTTATCCATATTTCAATCCTTGGAACTGTCCAGCTCCCATGACTCTA
TGCCACCAAGCTCTTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001003827

Insert Size:	1467 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001003827.1
RefSeq Size:	2326 bp
RefSeq ORF:	1467 bp
Locus ID:	53840
UniProt ID:	Q9BYJ4
Cytogenetics:	11p15.4
Protein Families:	Druggable Genome
MW:	56.9 kDa
Gene Summary:	<p>The protein encoded by this gene is a member of the tripartite motif (TRIM) family. The TRIM motif includes three zinc-binding domains, a RING, B-box type 1 and B-box type 2 domain, and a coiled-coil region. Expression of this gene is up-regulated by interferon. This gene is mapped to chromosome 11p15, where it resides within a TRIM gene cluster. Alternative splicing results in multiple transcript variants. A read-through transcript from the upstream TRIM6 gene has also been observed, which results in a fusion product from these neighboring family members. [provided by RefSeq, Oct 2010]</p> <p>Transcript Variant: This variant (4) differs in the 5' UTR, compared to variant 1. Both variants 1 and 4 encode the same isoform (1, also known as the middle form).</p>