

## Product datasheet for **SC300587**

### GPRASP2 (NM\_001004051) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	GPRASP2 (NM_001004051) Human Untagged Clone
Tag:	Tag Free
Symbol:	GPRASP2
Synonyms:	DFNX7; GASP2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

**Fully Sequenced ORF:** >SC300587 representing NM\_001004051.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGATGGGGCAGAGATTGAGCCTAGTGCCCAAGCCAGCCTGAAAAGAAGGCTGGGGAAGAGGTTATC
GCTGGGCTGAGAGAGAGAATGATGTCCCTCTGGTGGTCAGACCCAAGGTTAGGACCCAGGCAACTACT
GGGGCAAGGCCAAAACCTGAGACCAAGTCTGTGCCTGCGGCAAGGCCAAAACCTGAGGCCAAGCAATG
TCTGGGGCAAGGCCAAAACCTGAGGTCCAAGTAATGGGTGGTGAAGACCCAAAACCGAGGCTCAAGGA
ATCACAGGGGCCAGGCCAAAACCGATGCCAGGGCAGTAGGTGGCGCTCGTTCTAAAACCTGATGCCAAG
GCAATCCCTGGAGCAAGGCCAAGGATGAGGCCAGGCATGGGCCAGAGTGAATTTGGGACTGAAGCA
GTGTCACAGGCAGAAGGAGTGTCCAGACTAATGCCGTTGCTTGGCCACTGGCCACTGCTGAGTCTGGA
TCAGTTACTAAATCTAAGGGCTGTCTATGGATAGAGAACTAGTCAATGTGGATGCTGAAACCTTCTCT
GGCACCCAGGGTCAGAAAGGAATCCAGCCCTGGTTTGGACCAGGGGAGGAGACTAATATGGGGTCTTGG
TGCTATCCAGGCCAGGGCCAGAGAGGAGGCCTAATGAGTCTGGGTTCTGGTCAGCAGATGAGACC
TCTACAGCGTCTCTTTCTGGACTGGAGAAGAGACAAGTGTGAGATCATGGCCCAGGGAAGAGTCCAAT
ACCAGGTCAGGCACAGGGCTAAACATCAGACTAATCCCAGGTCCAGGCCAGATCCAAGCAAGAAGCC
TATGTTGATTCCTGGTCTGGATCTGAGGATGAGGCCAGCAACCCATTCTCCTTCTGGGTTGGAGAAAAT
ACCAATAACTTGTTCAGGCCAGAGTCAGGGAGGAGGCAATAATCAGGTCCAAGCTCAGGACAAATAGA
GAAGATTGTTTTGAATCTGAGTCTGAAGATGAGTTCTATAAGCAGTCTGGGTTTTGCCTGGAGAAGAG
GCCAATAGTAGATTACAGGCACAGAGACAAAAGAAGTCTAATACTGCCTTGAACTCAGGGCCAGAAA
GATGTTGACAGTGATAGGGTCAAACAAGAACCAGGTTTGAGGAGGAAGTCAATATTGGGCTCTGGTTC
TGGGCAGAAAAGAGGCCAGTTTGGAGGGTGGAGCTTCAGCAATCTGTGAATCTGAGCCAGGAAGTGA
GAGGGGGCCATTGGCGGATCCCGCTACTGGGCTGAGGAAAAGTCCAGTTTGGGGGCTGTGGCCAGAGAA
GAGGCCAAGCCGAGTCTGAAGAAGAGGCCATATTTGGGTCTGGTCTGGGACAGAGATGAGGCCCTGC
TTTGACCTAAATCCCTGTCTGTGTACAAGTCAAGTCAAGTTCAGAGATGCAGCTGAGGAGCTTAAT
GCATCCTCCAGGCCAAAACCTGGGACGAGGTCCTGTTGAATTCAAAACCTGGTCTTTTTTCATGGGGTT
GGCTCCGATCCACAAGCCCTTTGGAATTCGGAAGAGGCTTCTGAAATGCTTGAGGCAAGGCCAAG
AACCTGGAATTAGCCAGAAAGGAGAAGAGCAGGAATCTTTGCTCAGCCTGATCAGCCTAGTCTGAG
TTCACATTTAGTATGATCCTTCTACCGGTCAAGTCCGGGAAATTCGAGAGCATCTTAGGGCCAGGGAG
AGTGCAGAGTCTGAGAGTTGGTCTGCAGCTGCATACAATGTGAGCTGAAAATTGGTTCTGAAGAGTTT
GAAGAATTCCTTTTATAATGGACAAAATTCGGGATCCTTTTATTCATGAAATATCTAAAATTGCAATG
GGTATGAGAAGTCTTCAATTTACCCGAGATTTCAATTCGAGATTCAGGTGTTGTCTCACTTATTGAA
ACCTTGCTTAATATCCATCCTCTAGAGTTAGGACAAGTTTTTTGGAAAATATGATTCACATGGCTCCA
CCTATCCAAATCTAAACATGATTGAGACATTCATATGTCAAGTGTGTGAGGAAACCCCTTGACATAGT
GTGGATTCCTTGAGCAGCTGACTGGAATAAGGATGCTTAGACACCTCACTATGACTATTGACTATCAC
ACACTGATTGCCAACTATATGTCGGGTTTCTCTCCTTATTAACACAGCCAATGCGAGAACGAAGTTT
CACGTTCTGAAAATGCTATTGAATTTGTCTGAAAATCTGCTGTGGCAAAAAAATATTCAGTGCCAAA
GCTCTTTCAATATTTGTGGTCTCTTTAACATAGAAGAGACAAAATGATAATATTCAAATTTGTTATAAA
ATGTTTCAGAATATCAGTAACATTATAAAAAGTGGAAAGATGTCCTTAATTGATGATGATTCAGTCTT
GAGCCGCTTATTTCTGCATTTCTGAATTTGAGGAGTTAGCTAAGCAACTACAAGCCAAAATAGACAAC
CAAAATGATCCTGAGGTGGGACAACAAAGTTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**Plasmid Map:**


**ACCN:** NM\_001004051

**Insert Size:** 2517 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001004051.3](#)

**RefSeq Size:** 3818 bp

**RefSeq ORF:** 2517 bp

**Locus ID:** 114928

**UniProt ID:** [Q96D09](#)

**Cytogenetics:** Xq22.1

**Protein Families:** Druggable Genome

**MW:** 93.8 kDa

**Gene Summary:** The protein encoded by this gene is a member of a family that regulates the activity of G protein-coupled receptors (GPCRs). The encoded protein has been shown to be capable of interacting with several GPCRs, including the M1 muscarinic acetylcholine receptor and the calcitonin receptor. Several transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, May 2010]  
Transcript Variant: This variant (1) differs in the 5' UTR compared to variant 3. All five variants encode the same protein.