

## Product datasheet for **SC300599**

### **PALM2AKAP2 (NM\_001004065) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	PALM2AKAP2 (NM_001004065) Human Untagged Clone
Tag:	Tag Free
Symbol:	PALM2AKAP2
Synonyms:	AKAP-2; AKAP-KL; AKAP2; AKAPKL; MISP2; PALM2; PALM2-AKAP2; PRKA2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC300599 representing NM\_001004065.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGCGCTGGCCCCAGCCCGGGGCTGCCGCTCGCCTTCCCCGGAGTCTCCTGGACCCCCGGAGTCTCCT
GGACCCCCGGAGCGGGAGGCAGCGGCCCGCGCGCTGGACTGGAGCAGAACCCAGACTGCGCCCCC
GGGAGCGGGCGCCAGAGAAGCCTCCCCAGCTTTCTGAGGATGATATCTGGCTAAAAAGCGAGGGAGAC
AACTATAGTGCCACCCTCCTGGAGCCTGCTGCCAGCTCTCTTTCCCCAGATCACAAAAACATGAAATT
GAGGTGTCTGTTGCAGAAATGAAAAGTGTCTGGAATCACCTCTACCCACATCCCATGGACCATCCC
TCCGCTTTCTATTACCCCCGCATAATGGCCTCCTTACTGATCACACGAATCCCTGGATAATGATGTT
GCCAGAGAGATCCGCTATCTAGATGAGGTGCTAGAGGCCAACTGCTGTGATTCTGCTGTGGATGGAACG
TACAATGGAACATCCTCCCAGAGCCTGGTGCAGTGGTTCTGGTGGGCGGCCTAAGCCCCCTGTCCAC
GAGGCGACCCAGCCAGAACCCTGAAAGAACAGCTAGCCGGCAGGCACCTCTCACATCGAGCTCAGT
AATAGCAGCCCTGACCCATGGCAGAGGCAGAAAGAACAATGGCCATCCCCAGCCAGCTAGAGAT
GGCTGGGGGACAGCCTGCAGGTGCCTGTAGCCCCAGCTCCACAACCAGCTCACGGTGTCTTCCCGA
GATGGAGAGTTCACTCTACCACACTGAAAAAGGAGGCCAAGTTTGAGCTGCGTGCCCTCCATGAGGAC
AAGAAGCCCTCCAAGCTCTTTGAGGATGACGAGCATGAGAAAGAACAATACTGCATTAGAAAAGTGAGG
CCTTCAGAGGAGATGCTGGAGCTGAAAAGGAGAGGAGAGAGCTCATCCGAGCCAGGCCGTCAAGAAG
AATCCTGGCATTGCAGCAAAATGGTGAATCCCCCGCAGGAAAAAACCATCGAGGAGCAGCTGGACGAG
GAACATCTGGAGTCGCACAAAAGTACAAGGAGCGCAAAGAGAGAAGGGCACAGCAGGAACAGTTGCTG
CTGCAGAAGCAGTTACAGCAGCAGCAGCAGCAGCCCCATCGCAGCTCTGCACAGCCCCCTGCCTCTCT
CATGAACGCGCAAGCATGATTGACAAAGCAAAGGAGGACATTGTACAGAGCAGATAGATTTCTCTGCT
GCTCGCAAAACAATTTTCTGCTGATGGAGAATTTCCAGGCAAGCGGTGGCCAAGGGCCAGAGTACACCCAGG
CTGTTCTCCATCAAGCCTTTCTACAGGCCCTCTGGGGTCACTCAACTCAGACAAGCCACTGACTAATCCG
AGACCACCTTCTGTGGGGGACCTCCAGAAGACAGTGGTGCCTCAGCCGCAAGGGACAGAAATCCCCC
GGTGCCCTGGAGACCCCATCGGCAGCAGGAAGCCAGGGCAACACAGCCTCTCAGGGGAAGGAAGGGCC
TACAGCGAGCCTTCTAACGTGGGCCCTTATCTAAACTGTGGGCTGAGGATGGAGAATTTACGAGCGCC
CGGGCTGTCTCACTGTGGTCAAGGATGATGACCATGGGATTTTGGATCAGTTCTCAAGATCTGTCAAT
GTCTCCTTGACCAAGAGGAGCTTGACTCTGGTCTGGACGAATTGTCGGTGAAGTCTCAGGATACCACA
GTCCTGGAGACCTATCCAATGATTTTCTGATGGACAACATCAGTGACAGCGGGCATCCAATGAGACA
ACCAATGCCCTCCAGGAAAATTTCACTGGCTGATTTTTCTGCCCCAGACACCACAAACTGACAACCC
TCAGAGGGCCGAGGAGAAGGCGTCTCCAAGTCATTTAGTGATCATGGTTTCTATTTCCCTTCTCCACG
CTGGGGGACTCTCCGTTGGTTGATGACCCCTTGGAGTATCAGGCTGGCCTCCTGGTGCAGAAATGCCATT
CAACAAGCCATAGCCGAGCAGGTGGATAAAGCTGTGTCCAAAACCAGCAGGGATGGAGCAGAGCAACAG
GGACCTGAAGCGACTGTAGAGGAAGCTGAAGCTGCGGCTTTCCGCTCAGAAAAGCCTCAGAGCATGTTT
GAGCCACCTCAGGTGTCTTCTCCTGTTCAGAGAAAAGGGATGTATTACCAAAGATTCTGCCTGTGAA
GACAGGGCGCTCAGGGAAGGGGGCCCCCAGCCACTGCCAGCTGTGCAGCCCAGTGGCCCGATTAAAC
ATGGAGGAGACCAGGCCGAAGGAAGCTATTTTCAAGCAAGTACTCGGAGGCAGCTGAGCTGAGAAGCACA
GCCTCCCTCCTGGCCACTCAAGAATCTGACGTGATGGTTGGGCTTTCAAGCTGAGGTCCAGGAAACAG
CGGACTTTGTCCATGATTGAGGAAGAGATCCGAGCAGCTCAGGAAAGGGAAGAGGAGCTGAAGAGGAG
AGACAAGTCTTGACAGTACGCAGAGCCCCAGGACAAAGAATGCCCATCACTGCCCTCCAGAACATGC
TACAAAAGTCTCCAGGAAAAATAGAGAAAGTCAAACCTCCTCCATCCCCACCCTGAAGGCCCCAGC
TTGCAGCCTGACTTAGCCCTGAAGAGGCTGCCGGAACCCAGCGGCCAAGAATCTGATGCAGACCCCTC
ATGGAAGACTATGAGACACAAATCTAAAAGGCGGAGAGAATGGATGATAGTAGTGTCTCGAGGCC
ACACGGGTTAATCGAAGAAAGAGCGCACTGGCTTTGCGCTGGGAAGCAGGGATCTATGCCAACAGGAG
GAAGAAGACAACGAATAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

<b>ACCN:</b>	NM_001004065
<b>Insert Size:</b>	2847 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001004065.4</a>
<b>RefSeq Size:</b>	6977 bp
<b>RefSeq ORF:</b>	2847 bp
<b>Locus ID:</b>	445815
<b>UniProt ID:</b>	<a href="#">Q9Y2D5</a>
<b>Cytogenetics:</b>	9q31.3
<b>MW:</b>	104.1 kDa
<b>Gene Summary:</b>	<p>This gene belongs to the paralemmin downstream gene (PDG) family defined in PMID:22855693. Paralemmin downstream genes may have evolved contiguously with the paralemmin genes and are associated with other paralemmin paralogs in humans and several other taxa. The gene encodes three distinct protein isoforms, the PALM2 isoform, the AKAP2 isoform and the PALM2-AKAP2 isoform. The biological significance of the PALM2-AKAP2 isoforms is yet unknown. Earlier, PALM2 and AKAP2 were annotated as separate genes and PALM2-AKAP2 was annotated as a readthrough gene. [provided by RefSeq, May 2019]</p> <p>Transcript Variant: This variant (5) encodes isoform 5 (AKAP2). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>