

Product datasheet for **SC300622**

FYB2 (NM_001004303) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FYB2 (NM_001004303) Human Untagged Clone
Tag:	Tag Free
Symbol:	FYB2
Synonyms:	ARAP; C1orf168
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC300622 representing NM_001004303.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGAAGGGGAAGGGTAAGAACTTCAAGGAACCTCGAGCCAAATTTCAAATCTTGATGCTCCACCT
CTTCCAGGACCTATTAATTTCCAGCAGGTGTTCTCCAAAGGGTGACATTGGAGGCACACAGTCAACT
CAAATTTTGGCCAATGGGAAACCCCTCTCATCCAACCACAAGCAGCGCACACCATACTGTTCCAGTAGT
GAGTCCCAGCCTCTTCAACCTCAGAAAATAAAGTTGGCTCAGAAGAGTGAATTCAAAAATGTTCTAAC
TCCCCAGGCGCTCTGGGAAAGTCTACTGTATGTTCTGCAACAAGTTCACAGAAGGCTTCTCTGCTGTTA
GAGGTGACTCAATCAAATGTTGAGATAATCACTAAGGAAAAAGTAAATGGTGGCCAATAGCTTCAGAAAC
AACTCTGGAAGTGGGAGAAGGTTTCATCTCAGAAAAGTAAATGTTCTCAGCCCTTCTCCTTGCCAAC
TATGGAAGTAAGGCCATCCATCTGGAAGGGCAAAAAGGCATGGGGCTTACTCCAGAGGAACCCAGGAAA
AAGCTGGAACAAAAGGAGGCCAGACTTCTCTCCAGAACGACGTGGTGGCCCCCAAAATATTACAT
AACGTCTCTGAAGATCCCTCTTTTGTAAATTTCTCAACATATCAGAAAAAGCTGGGAAAACCCACCTCT
GAGAGGAGCCCGCAAGCAGCCCTGCCAGCCATCTATGAGTGTGAGCTTGCCAGTCAGGCCCCAGAA
AAACAGCCAGATGTCAGGCATCACACCTTCCCAAAACAAAGCCATTGCCCTCCATCGACTCCCTGGGT
CCTCTCCCCAAAGCCTTCAAGACCTCCCATCGTGAACCTCCAGGCCTTTCAGAGGCAGCCAGCTGCT
GTTCCCAAGACTCAGGGGAAGTGACTGTGGAAGAGGGCTCCCTGTCTCCAGAGAGGCTTTTCAATGCA
GAATTTGAAGAACCACATAATTACGAGGCAACAATTTCAATCTGAGACACTCTGGCACTCCATTAAC
CTGTGCACTGCAAAAGAAATGCTGATCCAACCTTATGAAGTTGGAATTGAAGAACTCCAAAGCCTGGG
AAGAATTTTCCCTATCCAGAACCTAGTGCTAAACATGAAGATAAAAAATGAAGGAAAAACAACCATGT
GAATTGAAACCTAAAAACACAGAAAAGGAACCATATTCAAACCATGTTTTCAAGTAGATGCCTGTGAA
GGGACACCTGAAAAAATTCAGATGACCAACGTCCACACAGGTAGAAGGAACATGTTGGCTGGAAAGCAA
GAGGCCATGATTGACATCATCCAGACAAATCCCTGCCCTGAGGGCCCAAGCTGGCCAGGCACTCCCAA
GGCCACTGTGGGCATCTGGAGGTTTTGGAGTCAACTAAAGAACTCCAGACCTAGGGGTCTCTAAGACA
AGTTCCATCTCGGAGGAGATATATGATGATGTCGAGTACTCCAGGAAAGAGGTACCGAAGCTGAACCTAC
TCTAGCTCACTTGCCCTCAAGTAGTGAAGAAAATAGAGAACTGTATGAAGATGCTACAAAACAAAGAAC
AACTACCCAAAGATAGATTTAGATGGAAGAAAGCACTCAAAAGACTGCAGCAATTTCTCAAGAAAGAA
AAGGATAGATTTAAATAAAGAAAACCAAGTCGAAAGAAAACCTAAGTGCATTTTCCATTTTGTGCTCT
GATTTAGAACTTAAGTCTCAGGAAGTTATTATTTATGATGATGTAGACCTGAGTAAAAAGAGTCAAAA
GATGAAGATAAACTGAAATGTGGAAGCCCAAGTTTCTGACACCAAGGAAAAAAGAGAAAAACGGT
GCTGAAGAATCAGAAAGTTTCTCTCCTAGAAATTTCTTCAAAACCAAGAAGCAAACTTAGAAAAAGAAC
AGAATGAAAGAGAGAAAAAATTTAGAGAAAGGTTTAAAGTACGACAAAGAGATTATTGTCATCAAT
ACAGCAGTGGCCTGTTCCAATAATTCAAGAAATGGAATATTTGATTTGCCAATAAGTCTGGAGAAGAA
TTGGAAGTCATTGATACCAACCAACAAATCTAGTGATATGTCGTAATTTCTAAAGGCAAAATATGGATAT
GTGCTCATTGAACATCTAGATTTCAAGCATCAAAGTTGGTCACCTTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: SgfI-MluI

ACCN: NM_001004303

Insert Size: 2187 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001004303.4</u>
RefSeq Size:	3466 bp
RefSeq ORF:	2187 bp
Locus ID:	199920
UniProt ID:	<u>Q5VWT5</u>
Cytogenetics:	1p32.2
MW:	82.1 kDa
Gene Summary:	Adapter protein that plays a role in T-cell receptor (TCR)-mediated activation of signaling pathways. Required for T-cell activation and integrin-mediated T-cell adhesion in response to TCR stimulation (PubMed:27335501).[UniProtKB/Swiss-Prot Function]