

## Product datasheet for **SC300648**

### PH domain containing family G member 7 (PLEKHG7) (NM\_001004330) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	PH domain containing family G member 7 (PLEKHG7) (NM_001004330) Human Untagged Clone
Tag:	Tag Free
Symbol:	PLEKHG7
Synonyms:	FLJ46688
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC300648 representing NM_001004330. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGATCTTTATGAATACTAAGATATCTGCAAACCTCATGAATATCTCCTAGATGTGGATTTATGGAGA
CTTTTTGCAAACCTGGAGGAGTTAACTCAGACAAGCCTTGGTTTTGTGAACAGTCTCTTTGGCATCATC
AAGGACTATGTAGACGCTTCTGAGATTTCTCATCACTGGATTTTATTTCCGTGCTCACAAAGTATTTT
CGAGGGAGTCTGTGCAGAGCCACCAGACCTACTGCCTGAACTATTCAGCTGCTATCTTTATCTTGAG
AGCCTGAGGCAGAGAGATGACTTTGGAATTTATTTAAATGGTGTGAGCAGAATGAACAATGCAGACGG
CTCCACGTGCCAGAGCTGCTAGTGGCCCCACTACAGAGGCTCACTCGATATCCGTTGTTGCTGAAGAA
ATCTGGAAAAGGAGCATGGACTCTGCTGAGAAAATCATGATCTACTCCATCAAGGAAAAGGTGGAAAAG
TCCATCCGGGACCTTGAAGGAAAAGTGAAGTGGCTGGACAATTTCCAAAAATTTAGATATCTACAGGAG
ATTATAGTGTGGCCACCGCTTTGGGATAGAGATAAAAGTTTTTCATTCCAGAGTGTGTAACACATT
TTAAAGAACACATGGCAGAAAACATCTTGTCCAAACCAGCAGACACCTTCTCTATGAAGGAAAATTA
ACTCTTGCAAGAACACGAGATTCTAGATGTTTATCTGTTTCTTCAATGATTTCTCTTAGTTACG
AAAATAAGTGAACAAAAGAACTTGGAGGCTCAGACCCTGGTTAATGTGTCCTTCTTACTCCT
GAGTTGCAAGCAGTAATAAAAGAGGGTGGTTCGTGTACAGTACTCGATCAGCCTATTCCACTAGATAGA
TTGGTAGTCAAAAGTATTGAACCACTCCATGTGTCAGTCTTTGGGCTGAGAAATGCTTTTCTTATACAA
CACGAAAACAGATATCGACAGTGTATAGCAGCATTCTTATTACAAGCCAAACGGAAAACATCAAAAA
ACATGGATGGCACAATAACAACCTGCAATTTCTTGCTTTACCAAGAGTCAGGAAACCAAGAAAATATCT
TTATTCACATTGCCCGCAGAATCCTCTGAAATTTAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online »](#)

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM\_001004330

Insert Size: 1140 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_001004330.2](#)

RefSeq Size: 3267 bp

RefSeq ORF: 1140 bp

Locus ID: 440107

Cytogenetics: 12q22

MW: 44.3 kDa