

Product datasheet for **SC300659**

TRIM67 (NM_001004342) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TRIM67 (NM_001004342) Human Untagged Clone
Tag:	Tag Free
Symbol:	TRIM67
Synonyms:	TNL
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

Fully Sequenced ORF: >OriGene ORF within SC300659 sequence for NM_001004342 edited (data generated by NextGen Sequencing)

```
ATGGAGGAAGAGCTGAAGTGTCCCGTGTGCGGCTCTGTGTTTCGGGAGCCTATCATCCTG
CCCTGTTCCCAATGTCTGCCTGCCTTGCCTCGCACCATCGCGGTGCAGACCCCGGAC
GGTGAGCAGCACCTGCCCCAGCCGCTCTGCTTCCCGGGGATCGGGGCTGCAGGCGGGC
GCCGCCCGCGCTGCCTCTCTGGAGCAGCAGCTGCGGCTGGCCCGGCTGCGGCGGTGCA
GGCGGAGTGCAGCTGGCGGCTCGGCGGCGGTGCGGGAGGTGGCGGAGACCACGCGGAC
AAGCTCAGCTTGTACAGCGAGACAGACGCGGCTACGGGTCTACACCCGAGCCTCAAG
TCCCCAACGGGGTTCGCGTGCTGCCCATGGTGGCCGACCAACCGGCTCCTCGGCTGCG
GCGGCTCGGGGTGCCGCTGCTCCTCGCTGCTCCTCGTCTTCGAGCTCCATCACGTGCCCC
CAGTGCCACCGCAGCGCATCCCTGGACCACCGCGGCTGCGCGGCTTCAGCGCAACCGG
CTGCTCGAGGCCATCGTGACGGTACCAGCAGGGCCGCGGGGCGGTGCCGGGACGTCT
GCAGCCGCGGCGGTGGCCATCTGCCAGCTGTGCGACCGCACCCCGCCAGAGCCAGAGCC
ACGCTCTGCGAGCAGTGCAGCTCCTCTACTGCTCTGCCTGCCAGCTCAAGTGCCATCCA
TCCCGGGGACCTTCGCAAGCATCGCCTGGTGCAGCCGCGCGCCGCCGCCGCCGCCGCC
GCCGAGGAGCCTCCGGGCCACTGGCACCGCCAGGGCGCCCCAGCGAGGCGGCGGCG
TGCAAGAGCCCGGAGGCGCGGGGCGGGGCGACTGGGGGAGCAGCGCCCGCAAGTTC
CCCAGTGTTCCCGAGCATGAAATGGAGAACTACAGCATGTAAGTGCCTGAGCTGTGCAACC
CCGGTGTGTTATCTGTGCTGGAGGAGGCGCGGACGCCAAGCAGAGGTGAAGCCGCTG
GGGGCCATGTGGAAGCAGCACAAGGCACAATATCTCAGGCCTTAAATGGAGTTTCAGAT
AAGGCAAAGGAAGCAAAGGAGTTTCTGGTTCAGCTAAAGAACATATTGCAGCAGATCCAG
GAAAACGGACTGGACTACGAAGCCTGCCTCGTTGCTCAGTGTGATGCCCTTGTGGATGCT
TTAAGTGTGAGAAAGCAAGCTGCTCAGCAAGGTGACTAAAGAGAGGGAACACAAGTTG
AAGATGGTTTGGGACAGATCAATCACTGCACATTGAAGCTGCGTCAGTCCACCGGACTG
ATGGAGTACTGCCTGGAGGTGATCAAGGAGAACGACCCCTCCGGGTTCTTACAGATCTCA
GATGCTCTGATCAAGCGCGTCCAGGTGTCTCAGGAGCAGTGGGTCAAAGGCGCCCTGGAG
CCGAAAGTGTCTGCGGAGTTTGATCTGACTTTGGACAGCGAGCCGCTGCTGCAGGCCATC
CACCAGCTGGACTTCATTAGATGAAATGTAGGGTGCCACCCGTCCTCCCTACTGCAGCTG
GAGAAATGCTGCACCCGTAACAACAGCGTCACGCTGGCCTGGAGGATGCCACCTTCACC
CACAGCCCGTGGACGGTACATCCTGGAGCTGGACGACGGTGCCGGGGGACAGTTCCGG
GAAGTGTACGTGCGTAAGGAGACTTTGTGTACCATCGACGGTCTTCACTTCAACAGCACC
TACAACGCCCGAGTCAAAGCTTTCAACTCTTCTGGTGTGCGGGCCTTACAGTAAAAGTGTG
GTCCTGCAGACATCCGATGTGGCCTGGTTCACATTTGACCCCAACTCTGGGCATCGGGAC
ATCATTTTATCCAATGACAACAGACAGCCACCTGCAGCAGCTATGACGACCGGGTGGTG
CTGGGCACAGCTGCGTTCTCAAGGGCGTGCACTACTGGAGCTGCACGTGGACCGGTAC
GACAACCAACCCAGACCCCGCTTCGGGGTGGCCAGGGCCAGCGTGGTCAAGGACATGATG
CTGGGCAAGGATGACAAGGCTGGGCCATGTATGTGGACAACAACCGCAGCTGGTTCATG
CACTGCAACTCCACACCAACAGGACGGAAGGTGGCGTGTGCAAGGGGGCCACCGTGGGC
GTGCTGCTGGACCTGAATAAGCACACTCTACCTTCTTCATCAACGGGCAGCAGCAGGGC
CCACAGCCTTCAGCCAGTGGACGGGTCTTCATGCCAGCCCTCAGCCTCAACCGCAAC
GTGCAGGTCAACCTGCACACAGGATTGGAAGTGCCGACTAACCTGGGGCGGCCAAAGCTG
TCAGGCAATTAG
```

Clone variation with respect to NM_001004342.3

Restriction Sites: Please inquire

ACCN: NM_001004342

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	The open reading frame of this TrueClone was fully sequenced and found to be a perfect match to the protein associated to this reference.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001004342.2, NP_001004342.2</u>
RefSeq Size:	3929 bp
RefSeq ORF:	2271 bp
Locus ID:	440730
UniProt ID:	<u>Q6ZTA4</u>
Cytogenetics:	1q42.2
Protein Families:	Druggable Genome