

Product datasheet for **SC300667**

C9orf173 (STPG3) (NM_001004353) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	C9orf173 (STPG3) (NM_001004353) Human Untagged Clone
Tag:	Tag Free
Symbol:	STPG3
Synonyms:	C9orf173
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC300667 representing NM_001004353. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTT TAGTGAACCGTCAGAATTTTGT AATACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCC GCGATCGCC
ATGATGAATTCTGACCAGAAGGCAGTGAATTCCTGGCAAATTTTACATCAATGGAGGCAAACTGG
ACCCATGGTCACCTGAGGCAGACACAACCAGAGCCACCCAGCCCAAGGCATCTGTGCTGTTGTTGGGC
CCGGAGCCGGGGATGGCCTGGGATGAGACACAGCCCCAAAGATGAAGGAGATCCCAGTTGGCCTCAGG
TTACAGACGGGCACCCCCAGGAGTCCCTGCCACCTACACCCAGACCCTGAGGGAATTTGCTGGAG
CAACGCCCTGATCACAGCTGACCTGGAAGTCCCAGCCCCACCAGGTACCAGGTGCCAGCCCCGTC
GTACGAGAGTCCCTCCACACCCCCACTACAGCATCGGCTGCAAGCACCAGGGCCGAGAGGGCGGTGCC
CGCAGGGCATGGCAGACTTTGTGGTCCAGAGCGAAAGCCCTTCACGCAGAAAGCTGACTTCGACCAA
GAGCAAAAGTGGCCCTCGCCAGCCCACTACCAGCTGCTCAGCCGGCCCGCCTTCCCGCCTTCAGCTTC
AGAGGCTGCCACTCCGCTTCCAAGACACCTGAAGGCCACACCCACCTAGGGCTGCCTGGGGCTAGGGGG
CTGGGCCCTCAGGGTGCAGCCCCAGTCCCTGCTTCAGGCCCTTTTGCAGGCACCTGGCAAGAGATGCCCT
GGCCCCAACACCTACAATATCCTTCCCTGGGAGCCGCTGCAGAGCCCCGCTCGCCGGCCTTCTCGATG
AGCCGCTCCCCTGCGTTCACCTCCTGGCTCAGCACCTCCCGAACCCCTGGCCAGCCGCCTACCACGTG
GAGGACTGCAACTCAGCTTCCCTTCGGCGCCTGGCGTGGTATCCAGGGTGTACGCAGACCCAAGCGC
CACGACACAGGCCCTTCTGCACGCTC TAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_001004353

Insert Size: 927 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001004353.3](https://www.ncbi.nlm.nih.gov/RefSeq/record/NM_001004353.3)

RefSeq Size: 1205 bp

RefSeq ORF: 927 bp

Locus ID: 441476
UniProt ID: [Q8N7X2](#)
Cytogenetics: 9q34.3
MW: 34.1 kDa