

Product datasheet for **SC300680**

LINGO4 (NM_001004432) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LINGO4 (NM_001004432) Human Untagged Clone
Tag:	Tag Free
Symbol:	LINGO4
Synonyms:	DAAT9248; LRRN6D; PRO34002
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC300680 representing NM_001004432.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGATGCAGCCACAGCTCCAAAGCAAGCCTGGCCCCCATGGCCCCGCTCCTTTTCTCCTCCTCCTA
CCTGGAGGGAGCGGTGGCAGCTGCCCTGCTGTGTGTGACTGCACCTCCCAGCCCCAGGCTGTGCTCTGT
GGCCACAGGCAACTGGAGGCTGTACCTGGAGGACTCCCACTGGACACTGAGCTCCTGGACCTGAGTGGG
AACCCGCTGTGGGGGCTCCAGCAGGGAATGCTCTCCCGCTGAGCCTGCTCCAGGAATTGGACCTCAGC
TACAACCAGCTCTCAACCCTTGAGCCTGGGGCCTTCCATGGCCTACAAAGCCTACTCACCTGAGGCTG
CAGGGCAATCGGCTCAGAATCATGGGGCTGGGGTCTTCTCAGGCCTCTCTGCTCTGACCCTGCTGGAC
CTCCGCTCAACCAGATTGTTCTTCTTAGTGGAGCTTTTGGGAGCTAGGCAGCCTCCAGAAGCTG
GAGGTTGGGGACAACCACCTGGTATTTGTGGCTCCGGGGCCTTTCAGGGCTAGCCAAGTTGAGCACC
CTCACCTGGAGCGCTGCAACCTCAGCACAGTGCCTGGCCTAGCCCTGCCGTCTCCCGCACTAGTG
GCCCTAAGGCTTAGAGAAGTGGATATTGGGAGGCTGCCAGCTGGGGCCCTGCGGGGGCTGGGGCAGCTC
AAGGAGCTGGAGATCCACCTCTGGCCATCTCTGGAGGCTCTGGACCTGGGAGCCTGGTTGGGCTCAAT
CTCAGCAGCCTGGCCATCACTCGCTGCAATCTGAGCTCGGTGCCCTTCCAAGCACTGTACCACCTCAGC
TTCTCAGGGTCTGGATCTGTCCAGAAATCCCATCTCAGCCATCCCAGCCCCAAGGCTCAGCCCCCTG
GTGCGGCTCCAGGAGCTACGCCTGTGAGGGCATGCCTCACCTCCATTGCTGCCATGCCTTCCATGGC
TTGACTGCCTTCCACCTCCTGGATGTGGCAGATAACGCCCTTTCAGACACTAGAGGAAACAGCTTTCCCT
TCTCCAGACAACTGGTACCTTGAGGCTGTCTGGCAACCCCTAACCTGTGACTGCCGCTCCTCTGG
CTGCTCCGGCTCCGCCGCCACCTGGACTTTGGCATGTCCCCCCTGCCTGTGCTGGCCCCCATCATGTC
CAGGGGAAGAGCCTGAAGGAGTTTTAGACATCCTGCCTCCAGGGCACTTCACTGCAAACCAAGCCCTG
ATCCGAAAGTCCGGGCCCTCGATGGGTATTGCAGAGGAGGGCGGGCATGCGGTTTTCTCCTGCTGGA
GATGGAGACCCAGCCCCACTGTCTCCTGGATGAGGCCTCATGGGGCTTGGCTGGCAGGGCTGGGAGA
GTAAGGGTCTAGAGGATGGGACACTGGAGATCCGCTCAGTGCAGCTACGGGACAGAGGGGCTATGTC
TGTGTGGTTAGCAATGTCGCTGGGAATGACTCCCTGAGGACCTGGCTGGAAGTCATCCAGGTGGAACCA
CCAAACGGCACACTTTCTGACCCCAACATACCGTGCCAGGGATCCCAGGGCCTTTTTTCTGGATAGC
AGAGGTGTGGCCATGGTGTGCGAGTGGCTTCTCCCTTCTCACCTCAGTACCCTCTGCTTTGGC
CTGATTGCCCTTTGGAGCAAGGGCAAAGGTCGGGTCAAACATCACATGACCTTTGACTTTGTGGCACCT
CGGCCCTCTGGGATAAAAACCTCTGGGGTAACCGGGTCACTGCCAAGCTCTTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001004432

Insert Size: 1782 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001004432.3](#)

RefSeq Size: 2682 bp

RefSeq ORF: 1782 bp

Locus ID: 339398

UniProt ID: [Q6UY18](#)

Cytogenetics: 1q21.3

Protein Families: Druggable Genome, Transmembrane

MW: 63.8 kDa