

Product datasheet for **SC300685**

ANKRD34B (NM_001004441) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ANKRD34B (NM_001004441) Human Untagged Clone
Tag:	Tag Free
Symbol:	ANKRD34B
Synonyms:	DP58
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC300685 representing NM_001004441.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGAATACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGGATGAAGGTATGGAATTTCAAGTGAAGGAAATTCCTTGATCAAAGCAGTCCATCAGAGCCGGCTT
CGCCTCACAAGACTTTTGTAGAGGCGGTGCCTACATTAATGAGAGCAACGACCGTGGGAAACCCCT
TTAATGATCGCTTGTAAAGACCAACATGTCGATCACCAGAGTGTCAGTAAAGCCAAAATGGTGAATAC
CTGTTAGAGAACAATGCCGATCCCAACATACAGGACAAATCTGGGAAAACGGCTCTGATGCATGCTTGC
TTAGAAAAAGCTGGCCCTGAAGTTGTTTCTTGCTCCTCAAGAGTGGGGCTGACCTCAGCTTGCAAGAC
CATTCTAGTTACTCAGCTCTTGTATGCTATAAATTCAGAAGATACAGAGACCTGAAAGTTCTTCTT
AGTGCTTGCAAGGCAAAAGGAAAGAGGTGTCATCATCACACAGCAAGTTGCCCTGTGGGAAGCAT
ACTACTAAACAATACTTAAATATGCCTCCTGTGGATATAGATGGGTGTCATCCCCAGCTACCTGCACC
ACTCCTTCAGAAATAGACATCAAACTGCCTCATCGCCACTTTCACATTCTTCTGAAACGGAAGTACG
CTTTTTGGCTTTAAAGATCTTGAGCTTGCTGGAAGCAATGATGATACCTGGGACCCAGGTTCCCTGTG
AGGAAACCTGCATTGGCCCTAAGGGGCCCAAGCTCCCCACGCTCCACCTGGGTCAAGAGTCCCCCA
TTATTAATGCACCAGAACAGAGTGGCTTCAATGCAAGAGGAGCTCCAGGATATTACACCAGAGGAAGAA
CTATCCTATAAAACCAATGGGCTGGCACTTTCAGCGGTTCACTAGGCACCAAGCATTGATGTA
AAAGACACTGCACATTTGCTAAGAGCCTTTGATCAGGCCAGCTCAAGGAAGATGTCATATGATGAATA
AATTGTCAATCTTATCTTTCAGAAGGAAATCAGCAATGCATTGAAGTCCCTGTTGACCAGGACCCAGAT
TCTAACAGACAATATTTGCTTCCACCTTAAGAAGTATAGTTCAAAAAGAACTTGGGGGCAATCAC
TACAGCTCTGATCCCAGCTCTCAGCTGGCCTTACCCCTCCAACCTCAGAAGACGGCAAAGCACTTATA
GGAAAGAAAAAGATCCTCTCACCATCTCCTTCCAATTGTCAGAGTCAAAGAATTGTTGGAGAATATC
CCCCAGGTCCCTGAGCAGGAGAAATCATGCAGTTTTAGAAAGGCGAGGTTCCAGGAGCTTCCCTTTA
GATCACAGTGTTACCCAAACCAGACAAGGGTTTCTCCACCCTTAAATGTAAATTCTCACCTCCCATC
TCAGATATCAATGTCAACAACAAGATTTGCAGCCTTCTTTCTTGTTGTTCAAAAAGTGCTTATGCCAACA
GTTCCGATTTTCCCTAAAGAATTCAAAGTAAGAAAATGTTGTTAAGGAGACAATCATTGCAAACTGAA
CAAATTAAGCAATTAGTAACTTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001004441

Insert Size: 1545 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001004441.2](#)

RefSeq Size: 3952 bp

RefSeq ORF: 1545 bp

Locus ID: 340120

UniProt ID: [A5PLL1](#)

Cytogenetics: 5q14.1

MW: 56.4 kDa