

## Product datasheet for SC300796

### OR51S1 (NM\_001004758) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	OR51S1 (NM_001004758) Human Untagged Clone
Tag:	Tag Free
Symbol:	OR51S1
Synonyms:	OR11-24
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC300796 representing NM_001004758. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTCAACATTACCAACTCAGATAGCCCCAATAGCAGCACTTCAATGGCCCCACCTTCTTGTGGTG
GGCATGCCAGGCCTATCAGGTGCACCCTCTGGTGGACATTGCCCTCATTGCTGTCTACCTTCTCTCT
GCACTGGGAAATGGCACCATCCTCTGGATCATTGCCCTGCAGCCCGCCCTGCACCGCCAATGCACTTC
TTCCTTCTTGTAGTGTCTGATATTGGATTGGTCACTGCCCTGATGCCACACTGCTGGGCATC
GCCCTTGTGGTCTCACACTGTCCTGCCTCAGCCTGCCTTCTACAGATGGTTTTATCCATGTCTTT
TCTGTCATGGAGTCTCTGTCTTGTCTGCCATGTCCATTGATCGGGCACTGGCCATCTGCCGACCTCTC
CACTACCCAGCGCTCCTACCAATGGTGTAAATAGCAAAATCAGCCTGGCCATTTCTTTTCGATGCCTG
GGTCTCCATCTGCCCTGCCATTCTGTGGCCTACATGCCCTACTGCCTCCCACAGTCTTAACCCAT
TCTTATTGCTTGCATCCAGATGTGGCTCGTTTGGCCTGCCAGAAGCTTGGGGTGCAGCCTACAGCCTA
TTTGTGGTCTTTTTCAGCCATGGGTTTGGACCCCTGCTTATTTTCTTCTCCTATGGCCTGATTGGCAAG
GTGTTGCAAGGTGTGGAGTCCAGAGAGGATCGTGGAAAGGCTGGTCAAACCTGTGCTGCCACCTCTCT
GCAGTGTCTCTTCTATATCCCTATGATCCTCCTGGCACTGATTAACCATCCTGAGCTGCCAATCACT
CAGCATAACCCATACTCTTCTATCCTATGTCCATTTCTTCTTCTCCATTGATAAACCCATTCTCTAT
AGTGTCAAGATGAAGGAGATTAGAAAGAGAATACTCAACAGGTTGCAGCCAGGAAGGTGGTGGTGTCT
CAGTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001004758
Insert Size:	972 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<a href="#">NM_001004758.1</a>
<b>RefSeq Size:</b>	972 bp
<b>RefSeq ORF:</b>	972 bp
<b>Locus ID:</b>	119692
<b>UniProt ID:</b>	<a href="#">Q8NGJ8</a>
<b>Cytogenetics:</b>	11p15.4
<b>Protein Families:</b>	Transmembrane
<b>Protein Pathways:</b>	Olfactory transduction
<b>MW:</b>	35.3 kDa
<b>Gene Summary:</b>	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]