

RefSeq: [NM_001005354.1](#), [NP_001005354.1](#)

RefSeq Size: 1051 bp

RefSeq ORF: 297 bp

Locus ID: 54458

UniProt ID: [Q9NZ81](#)

Cytogenetics: 12q13.13

Gene Summary: Negatively regulates TSP1 expression at the level of transcription. This down-regulation was shown to reduce taxane-induced apoptosis.[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (3) lacks an in-frame segment in the 5' coding region compared to variant 2. The encoded protein (isoform 2) is shorter than isoform 1.