

Product datasheet for **SC300917**

POTEG (NM_001005356) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	POTEG (NM_001005356) Human Untagged Clone
Tag:	Tag Free
Symbol:	POTEG
Synonyms:	A26C2; ACTBL1; CT104.4; POTE-14; POTE14; POTE14alpha; POTE22
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC300917 representing NM_001005356.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGTGGCTGAGGCTGGTTCAATGCCGGCTGCCTCCTCTGTGAAGAAGCCATTTGGTCTCAGAAGCAAG
ATGGGCAAGTGGTGCCCACTGCTTCCCCTGGTGCAGGGGAGCGGCAAGAGCAACGTGGGCATTCT
GGAGACCACGACGATTCTGCTATGAAGACTCAGGAGCAAGATGGGCAAGTGGTCCGCCACTGCTTC
CCCTGGTGCAGGGGAGCAGCAAGAGCAACGTGGGCATTCTGGAGACCACGACACTCTGCTATGAAG
ACACTCAGGAGCAAGATGGGCAAGTGGTCTGCCACTGCTTCCCCTGCTGCAGGGGAGCGGCAAGAGC
AAAGTGGGCCCTTGGGAGACTACGACGACAGCGCTTTCATGGAGCCGAGGTACCACGTCGCTCGAGAA
GATCTGGACAAGCTCCACAGAGCTGCCTGGTGGGTAAAGTCCCCAGAAAGGATCTCATCGTCATGCTC
AAGGACTGACATGAACAAGAAGGACAAGCAAAAGAGGACTGCTCTACATCTGGCCTCTGCCAATGGA
AATTCAGAAGTAGTAAACTCCTGCTGGACAGACGATGTCAACTTAATATCCTTGACAACAAAAAGAGG
ACAGCTCTGACAAAGGCCGTACAATGCCGGGAAGATGAATGTGCGTTAATGTTGCTGGAACATGGCACT
GATCCGAATATCCAGATGAGTATGAAAATACCGCTCTACACTATGCTATCTACAATGAAGATAAATTA
ATGGCCAAAGCACTGCTCTTATACGGTGTGATATCGAATCAAAAAACAAGCATGGCCTCACACCACTG
TACTTGGTGTACATGAGCAAAAACAGCAAGTGGTGAATTTCTTAATCAAGAAAAAGCAAAATTTAAAT
GCACTGGATAGATATGGAAGAACTGCTCTCATACTTGTGTATGTTGTGGATCGGCAAGTATAGTCAGC
CTTCTACTTGAGCAAAACATTGATGATCTTCTCAAGTCTATCTGGACAGACGGCCAGAGAGTATGCT
GTTTCTAGTCATCATAATGTAATTTGCCAGTACTTTCTGACTACAAAGAAAAACAGATGCTAAAAGTC
TCTTCTGAAAACAGCAATCCAGAACAAGACTTAAAGCTGACATCAGAGGAAGAGTCACAAAGGCTTAAA
GGAAGTAAAAATAGCCAGCCAGAGGAAATGTCTCAAGAACCAGAAATAAATAAGGGTGGTATAGAAG
GTTGAAGAAGAAATGAAGAAGCACGGAAGTACTCATATGGGATTCCAGAAAACCTGCCTAACGGTGCC
ACTGCTGACAATGGTGTGATGATGATTAATTCACCAAGGAAAAGCAGAACACCTGAAAGCCAGCAATTT
CCTGACTGAGAATGAACAGTATCACAGTGTGAACAAAATGATACTCAGAAGCAACTTCTGAAGAA
CAGAACACTGGAATATTACAAGATGAGATTCTGATTTCGAAGAAAAGCAGATAGAAGTGGCTGAAAAT
GAATTCTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites:

SgfI-MluI

Plasmid Map:



ACCN:

NM_001005356

Insert Size:	1527 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001005356.2</u>
RefSeq Size:	1926 bp
RefSeq ORF:	1527 bp
Locus ID:	404785
UniProt ID:	<u>Q6S5H5</u>
Cytogenetics:	14q11.2
MW:	57.1 kDa