

Product datasheet for **SC300930**

LRSAM1 (NM_001005373) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LRSAM1 (NM_001005373) Human Untagged Clone
Tag:	Tag Free
Symbol:	LRSAM1
Synonyms:	CMT2P; RIFLE; TAL
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC300930 representing NM_001005373.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGCCGCTTTCTTCCGGAAGCGGAAACCCAGTGAGGAGGCTCGGAAACGCTGGAGTACCAGATGTGT
TTGGCAAAGAAGCTGGGGCAGATGACATTCTCGACATCTCTAAATGTGAGCTCTCAGAGATTCCATTT
GGAGCTTTTGCAACATGCAAAGTTCTGCAGAAGAAGGTGCTGATCGTCCACACGAATCACCTCACTTCC
CTGCTTCCCAAATCCTGCAGCCTCCTGAGTCTGGCAACCATCAAGTTCTAGATCTCCACGATAATCAG
CTGACAGCCCTTCTGACGATCTGGGGCAGCTGACTGCCCTCCAGGTCTTAAACGTGAAAGGAATCAA
CTGATGCAGCTCCACGTTCCATTGGGAACCTGACCCAGCTCCAGACTCTCAATGTTAAAGACAACAAG
CTGAAGGAGCTTCCAGACACCGTGGGGGAGCTTGAAGCCTGCGTACCCTCAACATCAGTGGAAACGAG
ATCCAGAGATTGCCGAGATGCTGGCTCACGTTCAACCTGGAGATGCTGAGCCTTGACGCTCGGCC
ATGGTCTACCCGCCGCGGAGGTGTGTGGTCCCGGCACTGCGGCCATCTTGCAGTTCTCTGCAAAGAG
TCAGGGCTGGAATACTACCCCTTCTCAGTACTTGTGCAATCTGGAGCAAGATGGAATCGAGAAC
TCTCGGGACAGCCCTGATGGGCCACGGACAGATTCTAAGGGAGGAGTTAGAGTGGCAGAACAGGTTT
TCAGACTATGAGAAGAGGAAGGAACAGAAGATGCTGGAGAACTCGAGTTTGAACGGCGCCTGAACTG
GGCAGCGGGAGCACACCCAGCTCCTTCAGCAGAGCAGCAGCCAGAAGGATGAGATCCTTCAGACGGTC
AAGGAGGAGCAGTCCCGCTGGAGCAGGGCCTGAGTGAGCACCAGCGCCACCTCAACGCAGAGCGGCAG
CGGCTGCAGGAGCAGTGAAGCAGACGGAACAGAACATTTCCAGCCGGATCCAGAAGCTGCTGCAGGAC
AATCAGAGACAAAAGAAAAGCTCCGAGATTTTGAATCGTGGAAAATGAAAGAATAAGAATGGAACAG
TTGATGTCCATAACCCAGGAGGAGACTGAGAGCCTGCGGCGACGTGACGTTGCCTCCGCCATGCAGCAG
ATGCTGACTGAGAGCTGTAAGAACCGGCTCATCCAGATGGCCTACGAATCTCAGAGGCAGAACTGGTC
CAGCAGGCCTGTTCCAGCATGGCCGAAATGGATGAACGATTCCAGCAGATTCTGTGTCGTCAGCAAATG
GATCAGAACAAAGCCATCAGCCAGATCCTGCAGGAGAGCGGATGCAGAAGGCTGCGTTCGAGGCACTC
CAGGTGAAGAAAGACCTGATGCATCGGCAGATCAGGAGCCAGATTAAGTTAATAGAACTGAGTTATTG
CAGCTGACACAGCTGGAGTTAAAGAGGAAGTCCCTGGACACAGAGTCACTCCAGGAGATGATCTCGGAG
CAGCGCTGGGCCCTCAGCTCCCTGCTCCAGCAGCTGCTCAAAGAGAAGCAGCAGGAGGAAGAGCTC
CGGAAATCCTGACGGAGTTAGAAGCCAAAAGTAAACCAGGCAGGAAAATTAAGTGGCTGATTCAGTAT
CAACGGCTTTTGAACCAGAAGCCCTTGTCTTGAAGCTGCAAGAAGAGGGGATGGAGCGCCAGCTGGTG
GCCCTCCTGGAGGAGCTGTCGGCTGAGCACTACCTGCCATCTTTCGCGCACCACCGCCTCTCACTGGAC
CTGCTGAGCCAAATGAGCCAGGGGACCTGGCCAAGGTGGGCGTCTCAGAAGCTGGCCTGCAGCAGCAG
ATCCTCCGGAGAGTCCAGGAACTGCTGGATGCAGCCAGGATCCAGCCAGAGCTGAAACCACCAATGGGT
GAGGTCGTCAACCCTACGGCCCCCAGGAGCCTCCTGAGTCTGTGAGGCCATCCGCTCCCCCTGCAGAG
CTGGAGGTGCAGGCCTCAGAGTGTGTCGTGTGCCTGGAACGGGAGGCCAGATGATCTTCTCAACTGT
GGCCACGTCTGCTGCTGCCAGCAGTGTGCCAGCCACTGCGCACCTGCCCGCTGTGCCCGCAGGACATC
GCCAGCGCCTCCGCATCTACCACAGCAGCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001005373

Insert Size: 2172 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001005373.3](#)

RefSeq Size: 3146 bp

RefSeq ORF: 2172 bp

Locus ID: 90678

UniProt ID: [Q6UWEO](#)

Cytogenetics: 9q33.3-q34.11

Protein Families: Druggable Genome

MW: 83.6 kDa

Gene Summary: This gene encodes a ring finger protein involved in a variety of functions, including regulation of signaling pathways and cell adhesion, mediation of self-ubiquitylation, and involvement in cargo sorting during receptor endocytosis. Mutations in this gene have been associated with Charcot-Marie-Tooth disease. Multiple transcript variants encoding different isoforms have been identified for this gene. [provided by RefSeq, Jan 2012]

Transcript Variant: This variant (2) differs in the 5' UTR compared to variant 1. Variants 1, 2 and 3 encode the same isoform (1).