

## Product datasheet for **SC300952**

### **EBF3 (NM\_001005463) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	EBF3 (NM_001005463) Human Untagged Clone
Tag:	Tag Free
Symbol:	EBF3
Synonyms:	COE3; EBF-3; HADD3; O/E-2; OE-2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC300952 representing NM\_001005463.  
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTTTGGGATT CAGGAGAATATTCGCGCGGGGGACGACCATGAAGGAGGAGCCGCTGGGCAGCGGC
ATGAACCCGGTGGCGCTCGTGGATGCACACGGCGGGCGTGGTGGACGCCAACACGGCCGCCAGAGCGGC
GTGGGGCTGGCGGGGGCGCACTTCGAGAAGCAGCCGCCTTCCAACCTCCGGAAATCCAATTTCTCCAC
TTCGTGCTGGCGCTCTACGATAGGCAGGGGAGCCGGTGGAGATTGAAAGGACCGCTTTTGTGGACTTT
GTGGAGAAAGAGAAAGAGCCAAACAACGAGAAAACCAACAACGGCATCCACTATAAACTCCAGTTATTG
TACAGCAACGGAGTCAGAACAGAGCAAGATCTGTATGTTGCGCTCATAGATTCAATGACCAACAGGCC
ATCGTCTACGAGGGCCAGGACAAGAACCCGGAGATGTGCCGTGTGCTGCTGACCCACGAGATCATGTGC
AGCCGGTGTGTGACAAGAAAAGTTGTGGCAATAGAAACGAAACGCCCTCAGACCCTGTAATCATTGAC
AGATTCTTTCTAAAGTTTTCTCAAGTCAATCAGAAGTGTGGAAGAATGCAGGCAACCCCTCGAGAT
ATCGGGAGATTCCAGTTGTTGTATCGACAACAGTCAACGTGGACGGCCAGTGTGGCCGTGTCAGAC
AACATGTTTGTGACAACAATTCCAACACGGGAGGGCGGCCCGCCCTAGACCCGTGAGAAGCCACT
CCGTGCATCAAGGCCATCAGTCCAGTGAAGGCTGGACCAGGGGGGTGCCACCGTCATCATAATTGGC
GACAACTCTTTGACGGGCTGCAAGTTGTATTCGGAACATGTTGGTGTGGAGCGAGCTGATAACTCCC
CATGCCATCCGAGTCCAGACCCCGCCGAGGCACATTCCTGGCGTCTGCGAAGTGACCCTCTCCTACAAA
TCCAAGCAGTTCTGCAAAGGTGCTCCTGGGCGCTTTGTCTACACCGCCCTTAATGAACCAACCATAGAT
TACGGCTTTCAGAGGTTGCAGAAAGTGATCCCAAGACATCCGGGTGATCCCGAAAGGTTACCCAAGGAG
GTGTTACTGAAGCGGGCGGGACCTGGTGAAGCCTTATACGGAATGCCTCACAACAACCAGGAGATC
ATCTTGAAGCGAGCGGGCGGACATCGCCGAGGCGCTGTACAGCGTTCCTCCGCAATCACAACCAGATCCC
ACCCTGGGCAACAACCCTGCACACACGGGCATGATGGGCGTCAACTCCTTCAGCAGCCAGCTAGCCGTC
AACGTGTCAGAGACGTCAAGCAACGACCAAGTCGGCTACAGTCGCAATACAAGCAGCGTGTCCCGG
CGAGGCTACGTCCCAGCAGTACTCCCCAGCAGTCCAATTACAACACAGTCAGCACTAGCATGAATGGA
TATGGAAGTGGCGCCATGGCCAGTCTAGGGTCCCTGGCTCGCCTGGATTTCTTAATGGCTCCTCCGCT
AACTCTCCCTACGGCATGAAACAGAAGAGCGCCTTCGCGCCCGTGGTCCGGCCCCAAGCCTCTCCTCT
CCTCTCGCACCAGCGCCAACGGGAATGGACTGCAAGCTATGTCTGGGCTGGTAGTCCCGCCAATGTGA
```

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT  
TACAAGGATGACGACGATAAGGTTAAACGGCCGGC

Restriction Sites: Sgfl-Mlul

**Plasmid Map:**


**ACCN:** NM\_001005463

**Insert Size:** 1656 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001005463.2](#)

**RefSeq Size:** 4412 bp

**RefSeq ORF:** 1656 bp

**Locus ID:** 253738

**UniProt ID:** [Q9H4W6](#)

**Cytogenetics:** 10q26.3

**Protein Families:** Transcription Factors

**MW:** 60.4 kDa

**Gene Summary:** This gene encodes a member of the early B-cell factor (EBF) family of DNA binding transcription factors. EBF proteins are involved in B-cell differentiation, bone development and neurogenesis, and may also function as tumor suppressors. The encoded protein inhibits cell survival through the regulation of genes involved in cell cycle arrest and apoptosis, and aberrant methylation or deletion of this gene may play a role in multiple malignancies including glioblastoma multiforme and gastric carcinoma. [provided by RefSeq, Sep 2011]