

Product datasheet for **SC300990**

OR4C45 (NM_001005513) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	OR4C45 (NM_001005513) Human Untagged Clone
Tag:	Tag Free
Symbol:	OR4C45
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC300990 representing NM_001005513. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAACAATGTAACAGAATTCATCCTGCTGGGCCTCACTACAATCCAGAACTGCAGAAATTCCTGTTT
GTTATGTTTTTAATCACCTACTTGATCACACTGGCAGGTAACCTGCTCATCTCAGTCATCATCTTCATC
AGCCCAGCCCTGGGTCCCCATGTACCTTTTTCTGTCTATTTATCCATTATAGATATTTCTACTCT
TCTTCCATAGCCCTAAAATGATCTTTGACTTGATCTCTGAAAACAACACCATATCCTTCAATGGCTGC
ATGACTCAGCTCTTCACAGAACATTTCTTTCGGCAGCTGAGACCATCTTATTAAGTGCATGGCCTAC
GACTGCTATGTGGCCATCTGTAAGCCCTTGCACTATGCAACCATCATGACCCAATCTATGTGTGGATT
CTGATGGTGGTGGCTGGAATTCTGGGATTTGTGCATGGAGGAATCCAGACTTTGTTTCATAGCCCAGTTA
CCATTCTGTGGCCCAATGTCATCGACCACTTTATGTGTGATTTAGTACCTCTTCTGGAGCTGGCCTGC
ACAGACTCACACTTTAGGGCCTCTGATTGCTGCCAACAGTGGGTCAGTGTGTTTCTCATTTTTTCC
ATACTGGATGCTTCTATGTCATCATCTTGTGCTCCCTAAGGTCTCATAGCTCTGAAGGGCATCTCAA
GCTCTGTCTAGTTGTGCCTCTCATATCTTCACTGTCATCTTATTCTTTGTCCTTGTTCATACCTGTAT
CTAAGACCTTAACCTCCTTCCCCTGACAAAGCTGTGACTGTGTTTTGCACCCTATTTACACCTATG
TTGAACCCTTTAATCTACACTGTCAAAAATAAAGCAGTGAAAAATGTCATTAAGAAGCTCTGGAAGCAA
ATAATGACAACTGATGATAAATAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-Mlul



[View online »](#)

Plasmid Map:


ACCN: NM_001005513

Insert Size: 921 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001005513.1](#)

RefSeq Size: 921 bp

RefSeq ORF:	921 bp
Locus ID:	403257
UniProt ID:	A0A5F9ZGQ5
Cytogenetics:	11p11.12
Protein Families:	Transmembrane
Protein Pathways:	Olfactory transduction
MW:	34 kDa
Gene Summary:	<p>Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (1, coding) represents the functional allele encoding the full-length protein. ##RefSeq-Attributes-START## polymorphic pseudogene :: based on alignments, homology RefSeq Select criteria :: based on single protein-coding transcript ##RefSeq-Attributes-END##</p>