

Product datasheet for **SC300993**

OR5K3 (NM_001005516) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	OR5K3 (NM_001005516) Human Untagged Clone
Tag:	Tag Free
Symbol:	OR5K3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC300993 representing NM_001005516. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAATAAGGAGAATCACTCCTTGATAGCTGAGTTCATCCTCACAGGATTTACATATCATCCAAAGCTG
AAGACTGTTCTGTTTGTGGTGTCTTTGCCATCTATCTGATCACCATGGTGGGAAACATTGGTTTGGTG
GCATTGATTTATAGAGCAACGTCTTCACACACCAATGTACATATTTTAGGCAACCTAGTTCTGATG
GATTCCTGCTGTTCTCTGCTATTACTCCCAAGATGTTAGAGAATTTCTTTCTGAGGACAAAAGGATT
ACCCTGTATGAATGTATGGCACAATTTATTTTCTGTCTTGTGAAACTACAGACTGCTTTCTTCTG
GCGGCAATGGCCTATGACTGCTATGTGGCCATATGCAACCCACTGCAGTACCACACCATGATGTCCAAG
ACACTCTGCATTCAGATGACTGCAGGAGCCTACCTAGCTGGCAACCTGCATCCCATGATTGAAGTAGAG
TTTCTGTTGAGGTTAACTTTCTGTGGTCTCATCAAATCAATCATTTTTTCTGTGATGTTCTTCCATTA
TATAGACTTTCTGTATTAACCCTTATATCAATGAATTGGTGTATTTATCTTGGCAGGATCAATTCAA
ATCTTTACTATAGTCTTAGTTTCTATTCTACATCCTTTTCACAATATTTACAATGAAATCCAAGGAG
GGAAGAGGCAAAGCTTTTCTACTTGTGCATCTCACTTTCTCTGTGTCAATATCTGTGATTCCTT
CTTTCATGTATGCTCGACCAGGTGCAGTTAATGAAGGGGATAAAGATATACCTGTTGCTATATTTAT
ACTTTAGTTATTCCTTTATTAATCCTTTTATTTATAGCCTAAGAAATAAGGAAGTAATAAATATTATG
AAAAAATATGAAGAAGAGAAAATTTTGTACATTCTGAAACAAATGTCATCCCCTCTGGCAACTTGA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_001005516

Insert Size: 966 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in *E. coli* are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001005516.1](#)

RefSeq Size: 966 bp

RefSeq ORF: 966 bp

Locus ID: 403277

UniProt ID: [A6NET4](#)

Cytogenetics: 3q11.2

Protein Families: Transmembrane

MW: 36.7 kDa

Gene Summary: Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]