

Product datasheet for **SC301036**

DAZ2 (NM_001005785) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	DAZ2 (NM_001005785) Human Untagged Clone
Tag:	Tag Free
Symbol:	DAZ2
Synonyms:	pDP1678
Vector:	<u>pCMV6 series</u>
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001005785, the custom clone sequence may differ by one or more nucleotides

```

ATGTCTGCTGCAAACTCTGAGACTCCAAACTCAACCATCTCCAGAGAGGCCAGCACCCAG
TCTTCATCAGCTGCAGCTAGCCAAGGCTGGGTGTTACCAGAAGGCAAAATCGTGCCAAAC
ACTGTTTTTGTGGTGAATTGATGCTAGGATGGATGAACTGAGATTGGAAGCTGCTTT
GGTAGATACGGTTCAGTGAAAGAAGTGAAGATAATCACGAATCGAACTGGTGTGTCCAAA
GGCTATGGATTTGTTTCGTTTGTAATGACGTGGATGTCCAGAAGATAGTAGGATCACAG
ATACATTTCCATGGTAAAAAGCTGAAGCTGGGCCCTGCAATCAGGAAACAAAAGTTATGT
GCTCGTCATGTGCAGCCACGTCCTTTGGTAGTTAATCCTCCTCCTCCACCACAGTTTCAG
AACGTCTGGCGGAATCCAAACACTGAAACCTACCTGCAGCCCCAAATCACGCCGAATCCT
GTAACCTCAGCACGTTTCAGGCTTATTCTGCTTATCCACATTACACAGGTCAGGTCATCACT
GGATGTGCACTTGCTTGTATATAATTATCAGGAATATCCTACTTATCCCGATTTCAGCATTT
CAGGTCACCACTGGATATCAGTTGCCTGTATATAATTATCAGCCATTTCTCTGCTTATCCA
AGATCACCATTTTCAGGTCAGTCTGGATATCAGTTGCCTGTATATAATTATCAGGCATTT
CCTGCTTATCCAAATTCACCATTTCAAGTCGCCACTGGATATCAGTTCCCTGTATACAAT
TATCAGCCATTTCTCTGCTTATCCAAGTTCAACATTTTCAGGTCAGTCTGGATATCAGTTG
CCTGTATATAATTATCAGGCATTTCTCTGCTTATCCAAATTCACCATTTCAAGTCGCCACT
GGATATCAGTTCCCTGTATACAATTATCAGGCATTTCTCTGCTTATCCAAATTCACCACTT
CAGGTCACCACTGGATATCAGTTGCCTGTATACAATTATCAGGCATTTCTCTGCTTATCCA
AGTTCAACATTTTCAGGTCACCACTGGATATCAGTTGCCTGTATATAATTATCAGGCATTT
CCTGCTTATCCAAGTTCAACATTTTCAGGTCACCACTGGATATCAGTTGCCTGTATATAAT
TATCAGGCATTTCTCTGCTTATCCAAGTTCAACATTTTCAGGTCACCACTGGATATCAGTTG
CCTGTATATAATTATCAGGCATTTCTCTGCTTATCCAAGTTCAACATTTTCAGGTCACCACT
GGATATCAGTTGCCTGTATATAATTATCAGGCATTTCTCTGCTTATCCAAGTTCAACATTT
CAGGTCACCACTGGATATCAGTTGCCTGTATATAATTATCAGGCATTTCTCTGCTTATCCA
AGTTCAACATTTTCAGGTCACCACTGGATATCAGTTGCCTGTATATAATTATCAGGCATTT
CCTGCTTATCCAATTCAGCAGTTTCAGGTCACCACTGGATATCAGTTCCATGTATACAAT
TACCAGATGCCACCGCAGTGCCTGTTGGGGAGCAAAGGAGAAATCTGTGGACCGAAGCA
TACAAATGGTGGTATCTTGTCTGTTTAATCCAGAGAAGAGACTGA

```

Restriction Sites: Please inquire



[View online »](#)

ACCN:	NM_001005785
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001005785.1, NP_001005785.1</u>
RefSeq Size:	3756 bp
RefSeq ORF:	1605 bp
Locus ID:	57055
UniProt ID:	<u>Q13117</u>
Cytogenetics:	Yq11.223
Gene Summary:	This gene is a member of the DAZ gene family and is a candidate for the human Y-chromosomal azoospermia factor (AZF). Its expression is restricted to premeiotic germ cells, particularly in spermatogonia. It encodes an RNA-binding protein that is important for spermatogenesis. Four copies of this gene are found on chromosome Y within palindromic duplications; one pair of genes is part of the P2 palindrome and the second pair is part of the P1 palindrome. Each gene contains a 2.4 kb repeat including a 72-bp exon, called the DAZ repeat; the number of DAZ repeats is variable and there are several variations in the sequence of the DAZ repeat. Each copy of the gene also contains a 10.8 kb region that may be amplified; this region includes five exons that encode an RNA recognition motif (RRM) domain. This gene contains one copy of the 10.8 kb repeat. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jul 2008] Transcript Variant: This variant (2) has multiple differences in the coding region but maintains the reading frame, compared to variant 1. Variant 2 encodes isoform 2 which lacks one DAZ repeat, compared to isoform 1.