

## Product datasheet for **SC301080**

### **SIN1 (MAPKAP1) (NM\_001006619) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	SIN1 (MAPKAP1) (NM_001006619) Human Untagged Clone
Tag:	Tag Free
Symbol:	MAPKAP1
Synonyms:	JC310; MIP1; SIN1; SIN1b; SIN1g
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


**Fully Sequenced ORF:** >SC301080 representing NM\_001006619.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCCTTCTTGGACAATCCAACATCATTCTAGCTCATATTCGACAGTCACATGTGACCAGTGATGAC
ACGGGAATGTGTGAGATGGTTCTCATTGATCATGATGTTGACCTAGAGAAGATTCATCCTCCTTCAATG
CCTGGAGACAGTGGGTGAGAAATTCAGGGAAGCAATGGTGAAGCTCAGGGCTATGTATATGCCAGTCA
GTGATATTACCTCAAGTTGGGACTTTGGTATTAGAAGACGCTCAAACACAGCTCAAAGATTAGAACGA
CTCCGAAAAGAGAGACAAAACCAGATCAAATGCAAAAATATTCAGTGGAAAGAAAGAAATTTAAGCAA
TCAGCCCAGGAGTTAAAGTCACTGTTTGAATAAATCTCTCAAAGAGAAGCCTCCAATTTCTGGGAAG
CAGTCGATATTATCTGTACGCCTAGAACAGTGCCTCTGCAGCTGAATAACCTTTTAAACGAGTATTCC
AAATTTGATGGCAAGGTCATGTAGGTACAACAGCAACCAAGAAGATCGATGTCTACCTCCCTCTGCAC
TCGAGCCAGGACAGACTGCTGCCAATGACCGTGGTGACAATGGCCAGCGCCAGGGTGCAGGACCTGATC
GGGCTCATCTGCTGGCAGTATAACAAGCGAAGGACGGGAGCCGAAGCTCAATGACAATGTCACTGCCTAC
TGCTGCATATTGCTGAGGATGATGGGGAGGTGGACACCGATTTCCCGGCTGGATTCCAATGAGCCC
ATTCATAAGTTTGGCTTCAGTACTTTGGCCCTGGTTGAAAAGTACTCATCTCTGGTCTGCATCCAAA
GAGTCACTCTTTGTTTGAATAAATGCTGCTCATGGATTCTCCCTTATTCAGGTGGACAACACAAAGGTT
ACCATGAAGGAAATCTTACTGAAGGCAGTGAAGCGAAGAAAAGGATCCCAGAAAGTTTCAGGCCCTCAG
TACCGCCTGGAGAAGCAGAGCGAGCCCAATGTCGCCGTTGACCTGGACAGCACTTTGGAGAGCCAGAGC
GCATGGGAGTTCTGCCTGGTCCGCGAGAACAGTATCTCTGGAGACAAAGTAGAGATAGACCCTGTTACG
AATCAGAAAGCCAGCACTAAGTTTGGATTAAGCAGAAACCATCTCAATCGATTCCGACCTGCTCTGT
GCCTGTGACCTTGTGAAGAGAAAAGCCCAAGTACGCAATATTTAAACTCACGTATCTAAGCAATCAC
GACTATAAACACCTCTACTTTGAATCGGACGCTGCTACCGTCAATGAAATTTGTGCTCAAGGTTAACTAC
ATCCTGGAATCGCGAGCTAGCACTGCCCGGCTGACTACTTTGCTCAAAAACAAAGAAAAGTGAACAGA
CGTACGAGCTTCAGTTCCAGAAGGAGAAGAAATCCGGGCAGCAGTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**Plasmid Map:**



**ACCN:** NM\_001006619

**Insert Size:** 1428 bp

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_001006619.1</a></u>
<b>RefSeq Size:</b>	3254 bp
<b>RefSeq ORF:</b>	1428 bp
<b>Locus ID:</b>	79109
<b>UniProt ID:</b>	<u><a href="#">Q9BPZ7</a></u>
<b>Cytogenetics:</b>	9q33.3
<b>Protein Families:</b>	Druggable Genome
<b>MW:</b>	53.7 kDa
<b>Gene Summary:</b>	<p>This gene encodes a protein that is highly similar to the yeast SIN1 protein, a stress-activated protein kinase. Alternatively spliced transcript variants encoding distinct isoforms have been described. Alternate polyadenylation sites as well as alternate 3' UTRs have been identified for transcripts of this gene. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (3), also known as Sin1g, lacks an alternate in-frame exon in the 3' coding region, compared to variant 1. The resulting protein (isoform 3), also known as the gamma isoform, is shorter, compared to isoform 1.</p>