

Product datasheet for **SC301083**

WDR33 (NM_001006622) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	WDR33 (NM_001006622) Human Untagged Clone
Tag:	Tag Free
Symbol:	WDR33
Synonyms:	NET14; WDC146
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC301083 representing NM_001006622. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCTACAGAAATTGGTTCTCCTCCTCGTTTTTCCATATGCCAAGTTCAGCACCAGGCACCTCGA
CAGCTGTTTTATAAGCGACCTGATTTTGCACAACAGCAAGCAATGCAACAGCTTACTTTTATGGAAAA
CGAATGAGAAAAGCTGTGAACCGAAAAACCATAGACTACAATCCATCTGTAATTAAGTATTTGGAGAAC
AGAATATGGCAAAGAGACCAGAGAGATATGCGGGCAATTCAGCCTGATGCAGGTTATTACAATGATCTG
GTCCCACCTATAGGAATGTTGAATAATCCTATGAATGCAGTAACAACAAAATTTGTTCCGACATCAACA
AATAAAGTAAAGTGTCTGTATTTGTTGTTAGGTGGACTCCAGAAGGAAGACGCTTGGTCACTGGAGCT
TCTAGTGGGGAGTTTACCCTGTGGAATGGACTCACTTTCAATTTTGAAACAATATTACAGGCTCAGGAC
AGCCCAGTGAGGGCCATGACGTGGTCACATAATGACATGTGGATGTTGACAGCAGACCACGGAGGATAT
GTGAAATATTGGCAGTCGAACATGAACAACGTCAAGATGTTCCAGGCACATAAGGAGGCGATTAGAGAG
GCCAGGTTTATACACAATATACCATTTTCTGTAGTCCCTATTGTCATGGTTAAATATTCTCTAAGTGT
ATTCTGGGTGCAGAGATGCATGGGCTCTGTGAGTTTCTGGGAACTTTCTGCACCCTATAAAACAATA
TTTTCTTTGTTTTACACATTCACCATTTTGTGGCACCTTTCTGAAGTAGTGTGTCCCGGATACAG
CCTTTGCAATATGTTAGAGATGTACTGTCTGCCCATTTTGCAGTGGTTTTCTCTTTTCATTTATGATT
AATAATGTGTATACGTTATTCCTTTTTATTATCTACTGTGAAGACAAGAATATTTTCATTCCAAATAAA
GAATTCAGTCTTTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



Plasmid Map:



ACCN: NM_001006622

Insert Size: 981 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001006622.2](https://www.ncbi.nlm.nih.gov/RefSeq/record/NM_001006622.2)

RefSeq Size: 3183 bp

RefSeq ORF: 981 bp

Locus ID: 55339

UniProt ID: [Q9C0J8](#)

Cytogenetics: 2q14.3

Protein Families: Stem cell - Pluripotency

MW: 38.3 kDa

Gene Summary: This gene encodes a member of the WD repeat protein family. WD repeats are minimally conserved regions of approximately 40 amino acids typically bracketed by gly-his and trp-asp (GH-WD), which may facilitate formation of heterotrimeric or multiprotein complexes. Members of this family are involved in a variety of cellular processes, including cell cycle progression, signal transduction, apoptosis, and gene regulation. This gene is highly expressed in testis and the protein is localized to the nucleus. This gene may play important roles in the mechanisms of cytodifferentiation and/or DNA recombination. Multiple alternatively spliced transcript variants encoding distinct isoforms have been found for this gene. [provided by RefSeq, Jul 2008]

Transcript Variant: This variant (2) represents the shortest transcript. It lacks multiple 3' exons but has an alternate 3' segment, as compared to variant 1. The encoded isoform 2 has a shorter and distinct C-terminus, has only two WD repeats, and lacks the collagen-like and GPR domains, compared to isoform 1.