

Product datasheet for **SC301136**

ALG8 (NM_001007027) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ALG8 (NM_001007027) Human Untagged Clone
Tag:	Tag Free
Symbol:	ALG8
Synonyms:	CDG1H; PCLD3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC301136 representing NM_001007027.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCGGCGCTCACAATTGCCACGGGTACTGGCAATTGGTTTTCGGCTTTGGCGCTCGGGGTGACTCTT
CTCAAATGCCTTCTCATCCCCACATACCATCCACAGATTTTGAAGTACACCGAAACTGGCTTGCTATC
ACTCACAGTTTTGCCAATATCACAGTGGTATTATGAGGCAACTTCAGAGTGACGTTGGATTACCCCTT
TTCTTTGCATGGTTTGAGTATATCCTGTACATGTTGCCAAATATTTTGTCAAGAAATGCTGAATGTC
CATAATTTGAATTACTCCAGCTCAAGGACCTTACTTTTCCAGAGATTTTCCGTCATCTTTATGGATGTA
CTCTTTGTGTATGCTGTCGGTGTGCTGTAATGCATTGATGGAAAAAAGTGGGTAAAGAACTTACA
GAAAAGCCAAAATTTATTCTGTCGGTATTACTTCTGTGGAACCTCGGTTATTAATTGTGGACCATATT
CATTTTCAGTACAATGGCTTTTTATTTGGATTAATGCTACTCTCCATTGCAGGATATTTTCAGAAAAGG
CATATGGAAGGAGCATTTCTCTTTGCTGTTCTCCTACATTTCAAGCATATCTACCTCTATGTAGCACCA
GCTTATGGTGTATATCTGCTGCGATCCTACTGTTTCACTGCAAAATAACCAGATGGTCTATTGATGG
AAGAGTTTCAGCTTTGTTTCGTGTTATTTCCCTGGGACTGGTTGTTTTCTAGTTTCTGCTTTTCATTG
GGTCTTTTCTGGCCTTGAATCAGCTGCCTCAAGTCTTTTCCCGACTCTTTCTTTCAAGAGGGGCTC
TGTGATGCATATTTGGGCTCCAACTTCTGGGCTTTGTACAATGCTTTGGACAAAGTGTGCTGTGTCATC
GGTTTGAATTTGAAATTTCTTGATCCCAACAATATCCCAAGGCCTCAATGACAAGTGGTTTGGTTTCAG
CAGTTCCAACACACAGTCCTTCCCTCAGTGACTCCCTGGCAACCTCATCTGCACACTGATTGCCATA
TTGCCCTCTATTTCTGTCTTTGGTTTAAACCCCAAGGGCCAGAGGCTTTCTCCGATGTCTAACTCTT
TGTGCCCTTGAGCTCCTTTATGTTTGGGTGGCATGTTTCATGAAAAAGCCATACTTCTAGCAATTTCCCA
ATGAGCCTTTTGTCTGTGGAAAAAGCAGGAGACGCTTCGATTTTTCTGATTCTGACCACAACAGGACAT
TATCCCTCTTTCTCTGCTCTTCACTGCACCAGAACTTCCCATTAATAATCTTACTCATGTTACTATTC
ACCATATATAGTATTTTCGTCACCTGAAGACTTTATTCAGACGGAGTTTACCCCTGTTGCCAGGCTGGA
GTGCAATGGCACGATCTCAGTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites:

SgfI-MluI

Plasmid Map:



ACCN: NM_001007027

Insert Size: 1404 bp

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001007027.2</u>
RefSeq Size:	1743 bp
RefSeq ORF:	1404 bp
Locus ID:	79053
UniProt ID:	<u>Q9BVK2</u>
Cytogenetics:	11q14.1
Protein Families:	Transmembrane
Protein Pathways:	Metabolic pathways, N-Glycan biosynthesis
MW:	53.2 kDa
Gene Summary:	<p>This gene encodes a member of the ALG6/ALG8 glucosyltransferase family. The encoded protein catalyzes the addition of the second glucose residue to the lipid-linked oligosaccharide precursor for N-linked glycosylation of proteins. Mutations in this gene have been associated with congenital disorder of glycosylation type 1h (CDG-1h). Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (2) contains an additional segment in the coding region compared to variant 1. The resulting isoform (b) contains a shorter and distinct C-terminus compared to isoform a.</p>