

Product datasheet for **SC301158**

FSD2 (NM_001007122) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FSD2 (NM_001007122) Human Untagged Clone
Tag:	Tag Free
Symbol:	FSD2
Synonyms:	SPRYD1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC301158 representing NM_001007122.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGAGGAGGAATCGGGGAGGAAGTGGGGCTGGACAGGTCTACTCCAAGGATTTCCACTTTTACCAC
ATGGACCTGTATGACTCTGAAGACAGACTGCACCTCTTCCAGAAGAGAACTAGGATGAGGAAAGTA
GTCCAAGCTGAAATGGCCAATGAGTCAAGAGGAGCAGGGGATGGTAAGGCTCAAAGAGACCTTCAAGAG
GAAGTGGATGAACCTGTCCACTTATATGGACTTGAAGATGATCATGAATTAGGGGATGAGTTTGTGAT
GAAACATACCCAGAACAGGGGTTTCAAGATATCCTCCTTATATGATGAAGAGGAGAGACCCAGCCAGG
GAGCAGAGAGACTGGAGACTTAGTGGAGAGGCAGCGGAGGCCGAGGACCTGGGCTTCGGAGGGTGGGGC
TCAGCAGGCCAGTGCCAGGACTTGGGGAAGCCTATAGGTACACACACGGCCGTGCCAGCGAGGAGTAT
GAATGCTATGTCATCCCGAGGAAGAGGATGAAGAAGAAGCTGCTGATGTCTTGTGTCACTTGTAAAG
ACTCCAATAAGAGCTTTTCAAGAGTGGTGGTGAACATAAGGAGCATGAAGTATCCCACTCAATGAA
GCACTGGAAGTGGCAAGGATGAAATTCACAAAAACATGTACAAATTGGAAGAGCAGATAATTGAGATG
GAAACCTTTGCAATCATTTGGAGGAGGTTTTCATCACTGTGGAGGAGAAATTTGGAAGAACAGAACAA
AACTTTGAGTCACATTACAACGAGATCTTGGAACACTTGCTCAAAAATATGAGGAAAAATACAAGCT
CTAGGGGAGAAAAAGAAAGAGAGCTGGAAGCCTTATATGGACAAGTGGTCAGCTGTGGAGAAAACCTA
GATACTTGCAAGAACTGATGGAACAATAGAGGAGATGTGTACGAGGAGAAAGTGGATTTTATAAAG
GATGCAGTGGCTATGGCTGACAGACTCGGGAATTCCTGAAAACAAAAACAGACGTGGAATCTCTGCA
CAGCCTGAGTTTGAAGACCAGACCTTGGATTTCTGTGATGTGGAGCAGCTGATGGGCTCCATTAACACC
ATTCCAGCTCCTTCTGCTCCAGTCATAAATCCACAGGTCCTAACTCAGCCACAGGTTCTCAGTCCGA
GTGTGCTGGAGCTTATACTCCGACGACACTGTGGAGAGCTATCAGCTGTCTACCGGCCAGTGCAGGAC
AGCTCACCTGGGACGGACCAAGCAGAGTTTACAGTGACTGTCAAAGAAACATATTGCTCAGTGACAAAC
CTTGTCGCAATACCCAGTATGAATTTTGGGTCACAGCTCACAACAGGGCTGGCCCCAGCCCTCTAGC
GAGCGTGAGTGTACATGACAGCACCTTCTCCCCCATTATTAACCAAGAGATAAGGAGCTGTGAA
GAGGCTGTGCTGATTTGCTGGGAGTCTGGGAACCTGAATCCTGTGGACTCGTACACTGTGGAGCTGACC
CAGGCTGAAAGTCCAGAAGCCTCGGGTGAAGTGAAGTCTGTTGTGGGCATCCCGACCTGCGAGTCCGTG
GTGCAGCTGCAGCCGGGCGGAGCTACATTATCTATGTGCGAGCCCTCAATATGGGGGGCCCCAGCGTG
AGGAGCGAGCCAGCTACAGTCCACACCATAGGAAGCTACTTTCGCTAAACAAGGACACTTGCCATCCC
TGGCTGACCATTTCTGAAGACGAGCTTACGGCTGTACGAAGTGAAGGAGAACCCCGCCAGGGAGCTG
TACCCAGCGACACTCACTTACCAGGTGTGTTGCTGTATGGGAAATCTAATCCAGTCCGAGGGCAC
CATTACTGGGAGGTGGAGGTGGATGAGCATTGGACTACAGAGTTGGTGTGGCCTTTGCAGATGTCCGT
AAACAGGAGGATCTGGGAGCAAATTCCTCTCCTGGTGCATGAGGCACACATTTGCATCATCAAGGCAT
AAGTATGAATTTCTGCACAACAGAACTCCAGATATAAGAATAACAGTTCCACCAAAGAAGATTGGC
ATTCTATTAGACTATGAACATTCAGAGTTGTATTTTCAATGTGGACCTTTCTCAGCATCTATATACA
TTAGTTGTGACGTTACGAATTTGTGCATCCCTGTTTTCTTTGGAAAGCCTGGGTGTCTAAAGGTA
CATAATGGCATTCAATGCCAAACATGTCACTTTCTATAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: SgfI-MluI

ACCN: NM_001007122

Insert Size: 2250 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001007122.3</u>
RefSeq Size:	6436 bp
RefSeq ORF:	2250 bp
Locus ID:	123722
UniProt ID:	<u>A1L4K1</u>
Cytogenetics:	15q25.2
Protein Families:	Druggable Genome
MW:	85.4 kDa
Gene Summary:	This gene encodes a protein that belongs to the FN3/SPRY family of proteins. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Aug 2013] Transcript Variant: This variant (1) represents the longest transcript and encodes the longer isoform (a). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.