

Product datasheet for **SC301223**

LMBRD2 (NM_001007527) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LMBRD2 (NM_001007527) Human Untagged Clone
Tag:	Tag Free
Symbol:	LMBRD2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC301223 representing NM_001007527.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGAGTGGTGCAGCTTTGGGACTTGAGATTGTTTTGTCTTTTTCTGGCATTATTTCTGCTTCATCGA
TATGGAGACTTTAAGAAACAGCATAGACTTGTGATTATTGGAACACTGCTTGCTTGGTATCTCTGCTTT
CTTATTGTCTTCATACTGCCTCTGGATGTTAGTACGACAATATACAACCGGTGCAAGCATGCTGCTGCA
AATTCAAGCCCTCCTGAGAATAGCAACATTACAGGATTGTACGCAACTGTAACCCTGTTCCAAGCCAG
CATCCTTGTTCAGCCATGGAGTTACATTCTGATGGAATCATGCCAATTTCTGGAGGGTAGTGTAT
TGGACGTCACAATTTTAAACATGGATTCTTACCTTTTATGCAGTCATATGCAAGATCAGGAGGGTTT
TCCATCACTGGAAAGATCAAACTGCACATAATTGAGAATGCAATCTACTATGGCACCTATTTGCTGATT
TTTGGAGCATTTTAAATTTATGTAGCTGTAAACCCACATTTACATTTAGAATGGAACCAGCTTCAGACA
ATTGGGATAGCTGCTGCAATACATGGGTCTGTTTCTTCTGTGTTGTTGTTGGGTATGGCTTGGTG
GAAATTCCTCGATCACTGGAATGGAGCAAAAAGGGTTATCTACTTATGAAAACGTATTTAAGGCA
GCCAACTGATGACAGAGAAAGCAGATGCAGAAGAGAATTTGGAAGATGCCATGGAGGAGTTTCGTAAA
GTGAATGAAAGCATCAAGTATAATCACCCATTGAGAAAATGTGTTGATACAATACTTAAAAAGTGCCTT
ACAGAGTATCAGGAAAAAATGGGTAGGAACATGGATGATTATGAAGATTTTATGAAAAGCATAGTATC
TATCCAAGTAAAAAAGTCTGGTAAACTGCATAAACAGGTGATTTATTCAGTTCAGAGACACCGTCGA
ACTCAAGTACAATGGCAGATTCTTTGGAACAAGCATTTTATCTAGAAGATGTAGCAAAAAATGAACT
AGTGCTACTCATCAGTTTGTTCACACCTTTCAATCGCCAGAGCCAGAAAAATCGATTTATCCAATATTTT
TATAATCCTACATTTGAATGGTACTGGGAATGTCTTTGCGACCATGGTTTTACAAGATACTTGTGTG
GTTCTGTCCATCTTCTGTGATTGTTGTGGTCAGAGTGCACATTCTTAGCACTACTCCTGTCTTA
TCCCTCTTTGCGGTCTTCATACAGCTGGCAGAAAAACATATAATTATATTTATATCGAGATTGCCTGC
TTCCTTTCCATCTTCTTCTAAGTATCTGTGTTTATTCTACTGTGTTTCAGGATTCGTGTATTTAACTAT
TATTATTTGGCCTCACATCACCAGACTGATGCTTATAGCCTTCTTTTCACTGGCATGCTATTTTGGCGT
TTGACACCTCCTTTATGTCTTAATTTCTGGGTTTGACCCATATGGATTCATCTATCTCACAAGAAT
ACTCAACCACTGCTTATACATCTATTATGGGTTCCATGAAAGTTTTATCCTTTATTGCAGATGGATTC
TATATATATTATCCTATGTTGGTGGTAATTCTCTGCATTGCTACTATTTTAGTTTGGGAACCCGTTGT
TTGAATCTGCTCGGTTTCCAGCAGTTTATGGGAGATGATGATATGACATCAGACTTAGTTAATGAAGGA
AAAGAATTAATCAGAAAAGAGAAGAAAAAGGCAAGGAAGAAGGTGAAAATCGAAGAAGAGAA
TGGAAAGAACGTTATGGACACAATAGAGAAGATTCCACTAGGAACAGAAATATTCATACTGACCCAAA
GAGTCAAACCTCTCAGATGTTAATACCAACCGTTCTGCATTCAAAATATACCAGGGCTAATAACAGGACT
GAAAGGGACCGGATAGAACTTCTCAAGATGCAGAACCTTTGGATTTAATGCAGAAACATTCACTGAT
GATCCTCTTGAATCTGAATCAGGAAGGTATCAGCCTGGTGGACGATATCTCTCGATGTCTCGCAGTGAC
ATATTTAATGATGTTTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001007527

Insert Size: 2088 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001007527.1</u>
RefSeq Size:	3237 bp
RefSeq ORF:	2088 bp
Locus ID:	92255
UniProt ID:	<u>Q68DH5</u>
Cytogenetics:	5p13.2
Protein Families:	Transmembrane
MW:	81.2 kDa
Gene Summary:	Recruited to ligand-activated beta-2 adrenergic receptor/ADRB2, it negatively regulates the adrenergic receptor signaling pathway (PubMed:28388415). May also regulate other G-protein coupled receptors including type-1 angiotensin II receptor/AGTR1 (Probable). [UniProtKB/Swiss-Prot Function]