

Product datasheet for SC301228

Dysferlin interacting protein 1 (PPP1R27) (NM_001007533) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Dysferlin interacting protein 1 (PPP1R27) (NM_001007533) Human Untagged Clone
Tag:	Tag Free
Symbol:	Dysferlin interacting protein 1
Synonyms:	DYSFIP1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC301228 representing NM_001007533. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
 GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
 ATGCCTAGCAGAACTGCCGCTATGCCGCTACAGCCACGGCAGCGCGCGCGGATGCTGGCTGAT
 CGCAGCGTGCGTTTCCCTAATGATGTCCTGTTCTTGACCACATCCGGCAGGGTGACCTGGAGCAGGTG
 GGGCGCTTCATCCGACTCGGAAAGTCTCCCTGGCCACCATCCACCCCTCAGGCCTGGCCGCTTGCAT
 GAAGCCGTGCTCTTGAAACCTGGAATGCGTGAAGCTGCTGGTCAAATACGGGGCTGACATTACCAG
 CGAGATGAGGCGGGCTGGACACCCCTGCACATTGCCTGCAGCGATGGGTACCCTGACATAGCCAGGTAC
 CTTATCTCCCTGGGAGCGGACAGGGATGCAACCAACGACGATGGCGACCTGCCCTCCGACCTCATCGAC
 CCGGACTACAAGGAGCTGGTGGAGCTTTCAAAGGGACCACGATGGACTGA
 ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTAAACGGCCGGC

Restriction Sites:	SgfI-MluI
ACCN:	NM_001007533
Insert Size:	465 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).


[View online »](#)

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001007533.3</u>
RefSeq Size:	805 bp
RefSeq ORF:	465 bp
Locus ID:	116729
UniProt ID:	<u>Q86WC6</u>
Cytogenetics:	17q25.3
MW:	17.4 kDa
Gene Summary:	Inhibits phosphatase activity of protein phosphatase 1 (PP1) complexes.[UniProtKB/Swiss-Prot Function]