

Product datasheet for **SC301246**

IRGQ (NM_001007561) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	IRGQ (NM_001007561) Human Untagged Clone
Tag:	Tag Free
Symbol:	IRGQ
Synonyms:	FKSG27; IRGQ1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC301246 representing NM_001007561.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCTCCGCCCGCAGGGTGACGTGACCGCTTGTCTCTGGGGCTCCGGGCTTGGGGAAGTCCGCGCTG
ATCGCAGCGCTGTGCGACAAGGATGTGGAGACGCTCGAGGCCCGGAGGGACGGCCGACTCCGGGGTT
CCCAGCCTAAGAGCTGCGGGCCAGGCCCTTTTCTGGGCGAGCTGAGTGCCACCCGACGCGCCGGGG
CCCTGGGCGGCGAAGCCAACGTGCTGGTACTGGTGTGCTGCCGGACCCGAGGGGAACGGGAACCGTTG
GCTCCAGCCCTGGGAGAGGCAGCGCTGGCCGCCCTGGCCGAGGGACCCGCTACTGGTGTGCGGAAC
CTCCGCTCTGGGATTACAGAGTGCAGCCAGGCCCGTGATCAGACAGCAGCTCTGCTGAACAGCGCG
GGGTTAGGAGCTGCGGATCTGTTTGTGCTACCGGCGAAGTGGCGAGCAGCGACGGCTGCGAGGAGCTA
GAGCGCTCCGGGCGGCGCTGCAGAGCCAGGCAGAAGCGCTGCGGAGGCTCTGCCACCGGCGCAGGAT
GGCTTCGAGGTGTTGGGTGCAGCAGAGCTAGAGGCTGTGCTGAGGCCTTTGAGACCGGCGCGCTTGG
GCTGCGCTGTGCTGGGTGCGCTCAGGCCTGGAGCGCTGGGAGCGCACGCTAGACCTGGCCGTGGCT
GGCAAGGCTGACGTGGGCCCTTGTGGTGGACATGCTGCTTGGATTGGATCTGGCGACCCAGGCGCTGCG
CCTGCTTCGGTCCCCACAGCACCCTCCCTTCCAGCCCCAGAGCGCCGAATGTGGTGTCTTGACC
GTGCTCTGGGCCACACGGGCACTGCCACCACCGCGGCCGCCCTCTCACCAACGCACTACGACGCC
CTCATCTCGTACCCCTGGGGCCCCACTGAGAAGGACTGGGCCAGGTCCAGGCTTGTGCTACCA
GATGCGCTCTTGTCTGCGTGCGCACAGACGGCGAGGGCGAGGATCCGGAGTGTCTGGGAGAAGCAAG
ATGGAGAATCCCAAGGGCGAGAGCTTAAGAAGCAGGTGGAGGGGATTGGAGAATGCACTCAGTAAG
GGAAGGGAGAAATGTAGCGCTGGATCGCAGAAAGCAGGCAGCGGGGAAGTCTGGGAAAGCTGGCAGC
GAGGGTTTGCAGCAGTTGTGCGCATGAAGAAATCAGGTGGTGGCGACTCAGAGCGGGCCGCTGCGTTA
AGCCCGGAGGACGAGACGTGGGAGGTGCTGGAGGAGCGCCGCCAGTGTTCCTTACGGCCCTGGC
GGACTCCAGGGCTATGCGAATGGCTGCGGCGAGCGCTCCCCCAGCCAGGCGAGGGCACTGCTGCTG
GCGTTGCCACCAGCATCTCCAGCGCTGCCGAACCAAGGCTGCGGCGTTGCGAGCCGGGCGTGGAGG
CCAGCTCTGCTGGCTAGTCTGGCGGCGGCGGCGCACCCTCCAGGGCTGGGCTGGGATGCGACGTG
GCACCTTCTGCGGGGTGAGCTGGCGGAGTGGCGACGGGGCTGGGCTGGAACCCACGGCACTGGCTCGA
CGTGAGCGTGCCCTGGGCTGGCTTCTGGAGAGCTGGCAGCGCGCTCATTTCCAGGCCCGGTGACG
CGCGCCGAGGTGGAAGCAAGACTGGGCGCTGGGCGGGCGAGGGCACTGCTGGGGGCGCAGCACTGGG
GCTCTCTCTCTGTGGCTGCGGGTGGTGCAGCGGCGACAGGTGGCTGGGCTACCGAGCGGCTCAC
GGCGTCTGCTGAGGCTCTCGATGAGATGCGGGCTGATGCTGAGGCTGTGCTGGCACCCCTGAGCCT
GCCAGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001007561

Insert Size: 1872 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001007561.2](#)

RefSeq Size: 9703 bp

RefSeq ORF: 1872 bp

Locus ID: 126298

UniProt ID: [Q8WZA9](#)

Cytogenetics: 19q13.31

MW: 62.7 kDa