

Product datasheet for SC301248

TRK fused gene (TFG) (NM_001007565) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	TRK fused gene (TFG) (NM_001007565) Human Untagged Clone
Tag:	Tag Free
Symbol:	TFG
Synonyms:	HMSNP; SPG57; TF6; TRKT3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC301248 representing NM_001007565. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTT TAGTGAACCGTCAGAATTTTGT AATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCC GCGATCGCC
ATGAACGGACAGTTGGATCTAAGTGGGAAGCTAATCATCAAAGCTCAACTGGGGAGGATATTCGGCGA
ATTCCTATT CATAATGAAGATA TACTTATGATGAATTAGTGCTAATGATGCAACGAGTTTT CAGAGGA
AAACTTCTGAGTAATGATGAAGTAA CAATAAAGTATAAAGATGAAGATGGAGATCTTATAACAATTTT
GATAGTTCTGACCTTTCTTTGCAATTCAGTGCAGTAGGATACTGAAACTGACATTATTTGTTAATGGC
CAGCCAAGACCCCTTGAATCAAGTCAGGTGAAATATCTCCGTCGAGA ACTGATAGAACTTCGAAATAAA
GTGAATCGTTTATTGGATAGCTTGAACCACCTGGAGAACCAGGACCTTCCACCAATATTCCTGAAAT
GATACTGTGGATGGTAGGGAAGAAAAGTCTGCTTCTGATCTTCTGGAAAACAGTCTACTCAGGTTATG
GCAGCAAGTATGTCTGCTTTTGATCCTTTAAAAAACAAGATGAAATCAATAAAAAATGTTATGTCAGCG
TTTGGCTTAACAGATGATCAGGTTTCAGGGCCACCCAGTGCTCCTGCAGAAGATCGTT CAGGAACACCC
GACAGCATTGCTTCTCCTCCTCAGCAGCTCACCACAGGCGTT CAGCCACAGCAGCCACCATATA
GGAGCTCAGACTCAAGCAGGT CAGATTGAAGGT CAGATGTACCAACAGTACCAGCAACAGGCCGGCTAT
GGTGCACAGCAGCCGAGGCTCCACCTCAGCAGCTCAACAGTATGGTATTCAGTATTCAGCAAGCTAT
AGTCAGCAGACTGGACCTCAACAACCTCAGCAGTTCAGGGATATGGCCAGCAACCAACTCCCAGGCA
CCAGCTCCTGCC TTTTCTGGTCAGCCTCAACA ACTGCCTGCTCAGCCGCCACAGCAGTACCAGCGAGC
AATTATCCTGCACAAACTTACACTGCCAAACTTCTCAGCCTACTAATTATACTGTGGCTCTGCCTCT
CAACCTGGAATGCTCCAAGCCAACCTGGGGCCTATCAACCAAGACCAGGTTTTACTTCACTTCTGGA
AGTACCATGACCCTCCTCCAAGTGGGCCTAATCCTTATGCGCGTAACCGTCTCCTCCTTTGGTCAGGGC
TATACCAACCTGGACCTGGTTATCGA TAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



[View online >](#)

Plasmid Map:


ACCN: NM_001007565

Insert Size: 1203 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001007565.2](#)

RefSeq Size: 1841 bp

RefSeq ORF: 1203 bp

Locus ID: 10342

UniProt ID: [Q92734](#)

Cytogenetics: 3q12.2

Protein Pathways: Pathways in cancer, Thyroid cancer

MW: 43.4 kDa

Gene Summary: There are several documented fusion oncoproteins encoded partially by this gene. This gene also participates in several oncogenic rearrangements resulting in anaplastic lymphoma and mixoid chondrosarcoma, and may play a role in the NF-kappaB pathway. Multiple transcript variants have been found for this gene. [provided by RefSeq, Sep 2010]
Transcript Variant: This variant (2) differs in the 5' UTR, compared to variant 1. Variants 1, 2, and 3 encode the same protein (isoform 1).