

Product datasheet for SC301292

BAP29 (BCAP29) (NM_001008405) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	BAP29 (BCAP29) (NM_001008405) Human Untagged Clone
Tag:	Tag Free
Symbol:	BCAP29
Synonyms:	BAP29
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC301292 representing NM_001008405. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGACACTCCAATGGGCTGCAGTGGCAACCTTTCTTTATGCCGAAATAGGACTCATTTTAACTCTTCTGC
CTACCTTTTATTCCTCCTCAGAGATGGCAGAAGATTTTTTCATTTAATGTCTGGGGTAAAATTGCAACT
TTTTGGAACAAGGCTTTCCTTACCATTATCATCCTATTGATTGTTCTATTCTAGATGCTGTGAGAGAA
GTAAGGAAATATTCCTCAGTTCATACCATTGAGAAGAGCTCCACCAGCAGACCTGATGCCTATGAACAC
ACACAGATGAAACTTTTTAGGTCTCAAAGAAATCTTTACATTTCTGGATTTTCCCTATTTTTTTGGCTA
GTTTTGAGACGTCTGGTTACGCTTATTACTCAACTGGCAAAGAAGTGTCAAACAAGGTGTACTTAAA
ACTCAAGCAGAAAATACTAACAAGGCTGCCAAAAATTTATGGAAGAAAACGAAAACTAAAAGGATT
TTGAAAAGCCATGGTAAAGATGAAGAATGTGTTTTGGAAGCAGAAAATAAAAACTAGTAGAAGACCAG
GAGAACTGAAAAGTGAATTAAGGAAGACTTCAGATGCCCTTTCTAAGGCACAAAATGATGTGATGGAA
ATGAAGATGCAGTCAGAGAGACTTTGAAAGAATATGATCAACTCCTGAAAGAAGACTCTGAACTTCAG
CATTCCAGTTTTGGTGAATTTTTAAGCAAAAAGGCCACAAAATGGCAGCATTGGAAAACAGACTGGT
TCTAGAAAAGGATCGTTTAGAAAAGGGCAACAAGAAAAGACTGTGAACTTTATAAAAAGACTTGAAT
ATACTGTGTCAAAATGATAATTTTTGTTATGTTAGCCTCTAGAAAATTTAAGTTTCAGAAAATGCACTAT
GACCGTTTCGTAATTTTTTAATGCCACACATAGGTTGTATTGTAAATGGCATTATCAAAAATTTTGATG
ATGTTTCAGATATATTGCAAAGTCTGTATTCCAGCTCTTAAGAAAATATAAGCATGTTAAATACCATA
TTTACATATTGA
ACGCGTACGCGGCCGCTCGAGCAGAAAACATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_001008405

Insert Size: 1047 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001008405.2](#)

RefSeq Size: 5788 bp

RefSeq ORF: 1047 bp

Locus ID: 55973

UniProt ID: [Q9UHQ4](#)

Cytogenetics: 7q22.3

Protein Families: Druggable Genome, Transmembrane

MW: 40.8 kDa

Gene Summary: May play a role in anterograde transport of membrane proteins from the endoplasmic reticulum to the Golgi. May be involved in CASP8-mediated apoptosis (By similarity).
[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (1) encodes the longer isoform (a).