

Product datasheet for **SC301337**

CCBL2 (KYAT3) (NM_001008661) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CCBL2 (KYAT3) (NM_001008661) Human Untagged Clone
Tag:	Tag Free
Symbol:	KYAT3
Synonyms:	CCBL2; KAT3; KATIII
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC301337 representing NM_001008661. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTGTAGTGAACCGTCAGAATTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTTTTGGCCAGAGGAGCCTCTGCTCTTAGCGGTAGAGCAAAATCCTGAAGACAATTTCTCT
TCCAAAATCCTCGGATTCTACTTCTGCTAAAATGCTACTGAAATTCACAAATGCAAAACGGATTGAA
GGACTTGATAGTAATGTGTGGATTGAATTTACCAAATGGCTGCAGACCCCTTCTGTTGTGAATCTGGC
CAAGGCTTCCAGATATATCCCCTCTACATATGTAAGAAGAATTATCAAAGATTGCAGCAATCGAT
AGCCTGAATCAGTATACACGAGGCTTTGGCCATCCATCACTTGTGAAAGCTCTGCCTATCTGTATGAA
AAGCTTTATCAAAGCAAATTGATCAAATAAAGAAATCCTTGTGACAGTAGGAGCATATGGATCTCTT
TTTAACACCATTCAAGCATTAAATTGATGAGGGAGATGAAGTCATACTAATAGTGCCTTTCTATGACTGC
TATGAGCCCATGGTGAGAATGGCTGGAGCAACACCTGTTTTTATCCCCTGAGATCTAAACCTGTTTAT
GGAAAAAGATGGTCTAGTTCTGACTGGACATTAGATCCTCAAGAACTGGAAAGTAAATTTAATCCAAA
ACCAAAGCTATTACTAAATACTCCACATAACCCACTTGGCAAGGTGTATAACAGAGAGGAACTGCAA
GTAAATGCTGACCTTTGCATCAAATATGACACACTCTGCATCAGCGATGAGGTTTATGAATGGCTTGTA
TATTCTGGAAATAAGCACTTAAAAATAGCTACTTTCCAGGTATGTGGGAGAGAACAATAACAATAGGA
AGTGCTGGAAAGACTTTTCAGTGAACCTGGCTGGAAGCTTGGCTGGTCCATTGGTCCAAATCATTGATA
AAACATTTACAGACAGTTCAACAAAACAGATTTATACTTGTGCAACTCCTTTACAGGAAGCCTTGGCT
CAAGCTTTCTGGATTGACATCAAGCGCATGGATGACCCAGAATGTTACTTTAATTTGCCAAAAGAG
TTAGAAGTAAAAAGAGATCGGATGGTACGTTTACTTGAAGTGTGGCCCTAAAACCCATAGTTTCTGAT
GGAGGATACTTATCATCGCTGATGTGTCTTGTAGTCCAGACCTCTGATATGAAGAATAATGAG
CCTTATGACTATAAGTTTGTGAAATGGATGACTAAACATAAGAACTATCAGCCATCCCCGTTTCAGCA
TTCTGTAACCTCAGAGACTAAATCACAGTTTGTGAAAGTTTGTGCGTTTTTGTCTTCAATAAAAAAGACAG
CACTGGATGCTGCTGAAGAAATCATCAAGGCATGGAGTGTACAGAAGTCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```



[View online »](#)

Restriction Sites: Sgfl-MluI

Plasmid Map:



ACCN: NM_001008661

Insert Size: 1365 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001008661.2](#)

RefSeq Size: 2165 bp

RefSeq ORF: 1365 bp

Locus ID: 56267

UniProt ID: [Q6YP21](#)

Cytogenetics: 1p22.2

MW: 51.4 kDa

Gene Summary: This gene encodes an aminotransferase that transaminates kynurenine to form kynurenic acid, which is a metabolite of tryptophan. Multiple alternatively spliced transcript variants that encode different proteins have been described for this gene. This gene shares 5' exon structure with the RNA binding motif protein, X-linked-like 1 locus on chromosome 1, but the coding sequences are non-overlapping. [provided by RefSeq, Mar 2017]
Transcript Variant: This variant (1) represents the longest transcript and encodes the longest protein (isoform 1).