

## Product datasheet for **SC301377**

### **CCT3 (NM\_001008800) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	CCT3 (NM_001008800) Human Untagged Clone
Tag:	Tag Free
Symbol:	CCT3
Synonyms:	CCT-gamma; CCTG; PIG48; TCP-1-gamma; TRIC5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC301377 representing NM\_001008800.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGATGGGCCATCGTCCAGTGCCTGCTCAGCCAGAACACAAAGCGTGAATCCGGAAGAAAAGTTCAA
TCTGGAAACATCAATGCTGCCAAGATTCAGTCCAGCATCCAGCGGCCAAGTCCATGATCGAAATTAGC
CGGACCCAGGATGAAGAGGTTGGAGATGGGACCACATCAGTAATTATTCTTGACGGGAAAATGCTGTCT
GTAGCTGAGCACTTCTTGAGCAGCAGATGCACCCAACAGTGGTATCAGTGCTTACCGCAAGGCATTG
GATGATATGATCAGCACCCATAAGAAAATAAGTATCCCAGTCGACATCAGTGACAGTGATATGATGCTG
AACATCATCAACAGCTCTATTACTACCAAAGCCATCAGTCGGTGGTCATCTTTGGCTTGCAACATTGCC
CTGGATGCTGTCAAGATGGTACAGTTTGAGGAGAATGGTCGAAAAGAGATTGACATAAAAAAATATGCA
AGAGTGGAAAAGATACCTGGAGGCATCATTGAAGACTCCTGTGTCTTGCCTGGAGTCATGATTAACAAG
GATGTGACCCATCCAGTATGCGGGCTATATCAAGAACCCTCGATTGTGCTGCTGGATTCTTCTCTG
GAATACAAGAAAGGAGAAAGCCAGACTGACATTGAGATTACACGAGAGGAGGACTTCACCCGAATTCTC
CAGATGGAGGAAGAGTACATCCAGCAGCTCTGTGAGGACATTATCCAACCTGAAGCCCGATGTGGTCATC
ACTGAAAAGGGCATCTCAGATTTAGCTCAGCACTACCTTATGCGGGCCAATATCACAGCCATCCGAGA
GTCCGGAAGACAGACAATAATCGCATTGCTAGAGCCTGTGGGGCCGGATAGTCAGCCGACCAGAGGAA
CTGAGAGAAGATGATGTTGGAACAGGAGCAGGCCTGTTGGAAATCAAGAAAATTGGAGATGAATACTTT
ACTTTCATCACTGACTGCAAAGACCCCAAGGCCTGCACCATTCTCTCCGGGGGGCTAGCAAAGAGATT
CTCTCGGAAGTAGAACGCAACCTCCAGGATGCCATGCAAGTGTGTCGCAATGTTCTCTGGACCCCTCAG
CTGGTGCCAGGGGGTGGGGCTCCGAGATGGCTGTGGCCATGCCTTGACAGAAAAATCCAAGCCATG
ACTGGTGTGGAACAATGGCCATACAGGGCTGTTGCCAGGCCCTAGAGGTATTCTCTGACCTGATC
CAGAACTGTGGGGCCAGCACCATCCGTCTACTTACCTCCCTTCGGGCCAAGCACACCCAGGAGAAGTGT
GAGACCTGGGGTGTAAATGGTGAGACGGTACTTTGGTGACATGAAGGAAGTGGGCATATGGGAGCCA
TTGGCTGTGAAGCTGCAGACTTATAAGACAGCAGTGGAGACGGCAGTTCTGCTACTGCGAATTGATGAC
ATCGTTTCAGGCCACAAAAGAAAGGGATGACCAGAGCCGGCAAGGCGGGGCTCTGATGCTGGCCAG
GAGTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:**

Sgfl-MluI

**Plasmid Map:**



**ACCN:**

NM\_001008800

<b>Insert Size:</b>	1524 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001008800.2</a>
<b>RefSeq Size:</b>	1993 bp
<b>RefSeq ORF:</b>	1524 bp
<b>Locus ID:</b>	7203
<b>UniProt ID:</b>	<a href="#">P49368</a>
<b>Cytogenetics:</b>	1q22
<b>MW:</b>	56.4 kDa
<b>Gene Summary:</b>	<p>The protein encoded by this gene is a molecular chaperone that is a member of the chaperonin containing TCP1 complex (CCT), also known as the TCP1 ring complex (TRiC). This complex consists of two identical stacked rings, each containing eight different proteins. Unfolded polypeptides enter the central cavity of the complex and are folded in an ATP-dependent manner. The complex folds various proteins, including actin and tubulin. Alternate transcriptional splice variants have been characterized for this gene. In addition, a pseudogene of this gene has been found on chromosome 8. [provided by RefSeq, Aug 2010]</p> <p>Transcript Variant: This variant (3) lacks two alternate in-frame exons, compared to variant 1, resulting in a shorter protein (isoform c), compared to isoform a.</p>